

Aug 19, 2025

🌐 SMART-9N amplification enabling rapid nanopore sequencing of RNA viruses

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ld6yr8g5b/v1>

Vanessa Smilansky¹

¹SecureBio

SecureBio Detection



Vanessa Smilansky

SecureBio

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ld6yr8g5b/v1>

External link: <https://naobservatory.org>

Protocol Citation: Vanessa Smilansky 2025. SMART-9N amplification enabling rapid nanopore sequencing of RNA viruses . protocols.io <https://dx.doi.org/10.17504/protocols.io.n92ld6yr8g5b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 15, 2025

Last Modified: August 19, 2025

Protocol Integer ID: 222535

Keywords: wastewater, wastewater sample, viromics, metagenomics, nanopore sequencing, untargeted sequencing, RNA, viral RNA, untargeted RNA, total rna extraction of virus, diverse range of rna virus, nucleic acid extraction of virus, rna virus, total rna extraction, enabling rapid nanopore, rapid oxford nanopore technology, viral rna extract, rapid nanopore, nucleic acid extraction, rna, nucleic acid observatory, nasal swab, high genome coverage, untargeted sequencing, virus, workflow for untargeted sequencing, wastewater

Funders Acknowledgements:

Open Philanthropy

Grant ID: NA

Abstract

This protocol describes a rapid Oxford Nanopore Technologies (ONT) workflow for untargeted sequencing of RNA viruses. It has been optimized for wastewater and nasal swabs processed using the Nucleic Acid Observatory (NAO) protocols [Concentration and nucleic acid extraction of viruses from wastewater influent V.3](#) and [Concentration and total RNA extraction of viruses from nasal swabs V.1](#), respectively. It assumes that viral RNA extracts have undergone background depletion and does not detail those steps. In brief, RNA is reverse transcribed using SMART-9N amplification and prepared for sequencing with ONT's Rapid PCR Barcoding Kit 24 V14 (SQK-RPB114.24). In our experience, this approach consistently yields sufficient material for ONT's recommended inputs and recovers a diverse range of RNA viruses from both wastewater and nasal swabs, with long reads, high genome coverage, and reduced background noise.

Guidelines

RNA processing and handling: Please review [Protocol Note: Working with RNA Samples](#) before handling RNA samples.

Materials

Consumables:

- 100 μ M RLB RT 9N primer (IDT - Sequence: 5'-TTTTTCGTGCGCCGCTTCAACNNNNNNNNN-3')
- 100 μ M RLB TSOmG template-switching oligo (IDT - Sequence: 5'-GCTAATCATTGCTTTTTTCGTGCGCCGCTTCAACATmGmGmG-3')
- LunaScript RT Master Mix Kit (NEB, E3025)
- Sequenase Version 2.0 DNA Polymerase (ABI, 70775Y200UN)
- Rapid PCR Barcoding Kit 24 V14 (ONT, SQK-RPB114.24)
- 2X LongAmpTaq Master Mix (NEB, M0287)
- Agencourt AMPure XP beads (Beckman CoulterTM, A63881)
- Nuclease-free water (e.g. ThermoFisher, AM9937)
- 100% ethanol
- 1.5 mL Eppendorf DNA LoBind tubes
- 0.2 mL thin-walled PCR tubes
- QubitTM dsDNA HS Assay Kit (ThermoFisher, Q32851)
- QubitTM Assay Tubes (Invitrogen, Q32856)

Equipment:









- Ice bucket with ice
- Microfuge
- Micropipettes (1000 μ L, 200 μ L, 20 μ L, 2 μ L) and holder
- Timer
- Vortex mixer
- Thermal cycler
- Magnetic separation racks: compatible with 0.2 mL tubes and 1.5 mL tubes
- Hula mixer (gentle rotator mixer)
- QubitTM fluorometer

Before start

Prepare a biosafety cabinet or dead air box for RNA processing and handling. Gather materials and reagents. Ensure proper PPE.

Stage 1. cDNA synthesis and clean-up

1h 16m

- 1 Prepare  12 μL of each viral RNA extract in a separate clean 0.2 mL tube. 1m
- 2 Add  2 μL of RLB RT 9N (100 μM) to each 0.2 mL tube containing viral RNA extract. 1m
- 3 Incubate the reaction at  65 $^{\circ}\text{C}$ for  00:05:00 , then transfer the samples to ice immediately. Keep the samples on ice for  00:02:00 . 7m
- 4 Add  2 μL of TSOmG (100 μM) followed by  4 μL of LunaScript RT Master Mix Kit (Primer-free) to each 0.2 mL tube for a total volume of  20 μL . 2m
- 5 Incubate in thermal cycler using the following protocol: 13m














Step	Temperature	Time
Primer Annealing	25 $^{\circ}\text{C}$	2 minutes
cDNA Synthesis	55 $^{\circ}\text{C}$	10 minutes
Heat Inactivation	95 $^{\circ}\text{C}$	1 minute

- 6 Prepare the following Sequenase master mix (adjusted to +10% volume) in a 1.5 mL tube, according to the total number of samples: 5m

A	B	C
		Enter the number of samples:
Reagents	Volume (μL) per sample	2
5X Sequenase Buffer	2	4.4

	A	B	C
Sequenase Dilution Buffer		0.9	1.98
Sequenase		0.6	1.32
NFW		7.7	16.94
TOTAL		11.2	24.64

Thoroughly mix the reaction by gently pipetting and briefly spinning down.

- 7 Add  11.2 μL of Sequenase master mix to each sample. 1m
- 8 Incubate at  37 °C for  00:08:00 . 8m
- 9 Resuspend AMPure XP Beads (AXP) by vortexing. Add  25 μL (0.8X) of AXP beads to each sample and mix by flicking the tube. 2m
- 10 Incubate on a Hula mixer at  Room temperature for  00:05:00 . 5m
- 11 Spin down and pellet on a magnet for  00:05:00 ; remove supernatant. 6m
- 12 Add  200 μL freshly prepared 80% ethanol (do not remove from magnet); incubate  00:00:30 , remove ethanol. Repeat. 2m
- 13 Dry spin and remove residual ethanol. Allow pellet to air dry for about  00:00:30 on the magnet, but do not dry the pellet to the point of cracking. 2m
- 14 Remove the tube from the magnet and resuspend pellet in  11 μL NFW. 1m
- 15 Incubate on a Hula mixer at  Room temperature for  00:10:00 . 10m



- 16 Spin down, pellet on magnet, and transfer $10\ \mu\text{L}$ of eluate to a clean 0.2 mL PCR tube. 5m
- 17 Use $2\ \mu\text{L}$ of each sample to quantify with Qubit dsDNA HS assay. Use remaining $8\ \mu\text{L}$ sample for Rapid Barcode PCR. 5m

Note

Take your cDNA sample(s) forward to Stage 2 of the protocol.
Optional: if you are not using your sample(s) immediately, store at -20°C .

Stage 2. Rapid Barcode PCR

2h 10m

- 18 Prepare $5\ \mu\text{L}$ each of viral cDNA (diluted to a concentration of $1\ \text{ng}/\mu\text{L}$, for a total of 5 ng) in a separate clean 0.2 ml tube. 5m

Note

If the concentration of viral cDNA is $<1\ \text{ng}/\mu\text{L}$, then the volume can be increased to reach a total of 5 ng. Additional volume beyond $5\ \mu\text{L}$ must be subtracted from the $20\ \mu\text{L}$ of NFW added in step 19.

- 19 To each tube, add $20\ \mu\text{L}$ of NFW, $1\ \mu\text{L}$ of Rapid Barcode Primer (RLB01-24, at $10\ \mu\text{M}$), and $25\ \mu\text{L}$ 2X LongAmp *Taq* Master Mix. Mix by gently flicking the tube and spin down. 5m
- 20 Amplify using the following cycling conditions: 95°C for $00:00:45$; 25 cycles of 95°C for $00:00:15$, 56°C for $00:00:15$ and 65°C for $00:04:00$ followed by a final extension step of 65°C for $00:06:00$. Hold at 10°C . 2h

Note

Take your PCR product(s) forward to the clean-up, quantification, and adapter attachment steps.
Optional: PCR amplification may be performed and held at the final hold temperature overnight.

Stage 3. Clean-up, quantification, and adapter attachment

47m

21 Quantify the sample tubes (PCR products from the previous step) using the Qubit dsDNA HS Assay Kit. 5m


22 In a new 1.5 mL Eppendorf DNA LoBind tube, pool all barcoded samples in equal ratios to a combined final concentration 800 ng. 5m



Note


For example: if 10 barcodes were used, add 80 ng from each sample.


23 Resuspend the AXP beads by vortexing. Add 0.6X volume of AXP beads to the pooled samples. 2m


24 Incubate on a Hula mixer at  Room temperature for  00:05:00 . 5m

25 Spin down and pellet on a magnet for  00:05:00 ; remove supernatant. 6m



26 Add  500 μL freshly prepared 80% ethanol (do not remove from magnet); incubate  00:00:30 , remove ethanol. Repeat. 2m

27 Dry spin and remove residual ethanol. Allow pellet to air dry for about  00:00:30 on the magnet, but do not dry the pellet to the point of cracking. 2m

28 Remove the tube from the magnet and resuspend pellet in  14 μL Elution Buffer (EB). 1m

29 Incubate at  Room temperature for  00:02:00 . 2m

30 Spin down, pellet on magnet, and transfer  13 μL of eluate to a clean 1.5 mL tube. 5m


31 Use  2 μL of eluate to quantify with Qubit dsDNA HS assay. Use remaining  11 μL for adapter attachment. 5m

Expected result


The eluate should be approximately 10–40 ng/ μL .

32 In a new 1.5 ml Eppendorf DNA LoBind tube, prepare the Rapid Adapter (RA) dilution as described below and mix by pipetting: 1m

Reagent	Volume
Rapid Adapter (RA)	1.5 μL
Adapter Buffer (ADB)	3.5 μL
TOTAL	5 μL

33 Add  1 μL of diluted RA to the library pool. Mix gently by flicking the tube and spin down. 1m

34 Incubate at  Room temperature RT for  00:05:00 . 5m

35 The prepared library pool is used for loading into a MinION/GridION/PromethION Flow Cell. Keep the pool on ice or at  4 $^{\circ}\text{C}$ until ready to load.



Protocol references

[Rapid metagenomic sequencing for surveillance of bacterial, fungal and viral pathogens using SQK-RPB114.24](#)