

Oct 11, 2022

## Snap-Frozen Tissue Preparation

DOI

[dx.doi.org/10.17504/protocols.io.n92ldzoyxv5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ldzoyxv5b/v1)

Stephen Fisher<sup>1</sup>, Marielena Grijalva<sup>1</sup>, Rong Guo<sup>1</sup>, Sarah A Johnston<sup>1</sup>, Hieu Nguyen<sup>1</sup>, John Renz<sup>2</sup>, Jean G Rosario<sup>1</sup>, Steven Rudich<sup>2</sup>, Brian Gregory<sup>1</sup>, Junhyong Kim<sup>1</sup>, Kate O'Neill<sup>1</sup>

<sup>1</sup>University of Pennsylvania; <sup>2</sup>Gift of Life Donor Program



Stephen Fisher

University of Pennsylvania

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldzoyxv5b/v1>

**Protocol Citation:** Stephen Fisher, Marielena Grijalva, Rong Guo, Sarah A Johnston, Hieu Nguyen, John Renz, Jean G Rosario, Steven Rudich, Brian Gregory, Junhyong Kim, Kate O'Neill 2022. Snap-Frozen Tissue Preparation. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n92ldzoyxv5b/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 13, 2022

**Last Modified:** October 11, 2022

**Protocol Integer ID:** 68611

**Keywords:** biobanking of fresh frozen tissue, frozen tissue preparation this protocol, frozen tissue preparation, fresh frozen tissue, frozen tissue, preservation of tissue, snap freezing, stable in nonfixed surgical specimen, nci biospecimen evidence, intact tissue morphology, freezing protocol, nonfixed surgical specimen, freezing, national cancer institute biospecimen evidence, biobanking, placing tissue, dry ice, tissue, ideal for downstream assay, biopreservation, downstream assay, biospecimen, preservation, liquid nitrogen, rna, tissue preparation

**Funders Acknowledgements:**

NIH

Grant ID: U54HD104392



## Abstract

This protocol describes preservation of tissue by snap-freezing. Biospecimens preserved with this protocol will be suitable for downstream analysis of DNA, RNA, protein, and morphology endpoints. Multiple workflows can be used for snap freezing including using placing tissue in a vessel either with or without media and exposing it to liquid nitrogen vapor, immersing it in liquid nitrogen, or placing it on dry ice with a cooling device (Micke et al., 2006). The workflow pursued depends on the desired balance between feasibility, time, and cost as well as the possible downstream analyses performed. The snap-freezing protocol detailed below may not be ideal for downstream assays that require intact tissue morphology (Engel, Vaught, & Moore, 2014) (see [NCI Biospecimen Evidence-Based Practices](#)).

Micke, P., Ohshima, M., Tahmasebpour, S., Ren, Z.-P., Ostman, A., Pontén, F., & Botling, J. (2006). Biobanking of fresh frozen tissue: RNA is stable in nonfixed surgical specimens. *Laboratory Investigation; a Journal of Technical Methods and Pathology*, 86(2), 202–211. <https://doi.org/10.1038/labinvest.3700372>

Engel, K. B., Vaught, J., & Moore, H. M. (2014). National Cancer Institute Biospecimen Evidence-Based Practices: A novel approach to pre-analytical standardization. *Biopreservation and Biobanking*, 12(2), 148–150. <https://doi.org/10.1089/bio.2013.0091>

## Materials

- Freshly dissected tissue block
-  Dry Ice
-  1.5 mL LoBind tubes **Eppendorf Catalog #022431021**
- Access to -80°C ultra-low freezer

## Troubleshooting



- 1 Place pre-weighed piece of tissue in either a 1.5mL Eppendorf tube or 1.5mL cryo tube.
- 2 Quick freeze by immediately placing tube on dry ice.
- 3 Transfer tube to a -80°C ultra-low freezer for Biobanking or until ready for downstream processing.