



Aug 21, 2023

Soil DNA Extraction Modified Protocol for Dryland Agroecosystems

DOI

dx.doi.org/10.17504/protocols.io.eq2ly7domlx9/v1

McKenzie L. Stock¹, Paul Gabriel¹, Jennifer J. Randall¹, Nicole Pietrasiak²

¹New Mexico State University; ²University of Nevada Las Vegas



McKenzie L. Stock

New Mexico State University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly7domlx9/v1>

Protocol Citation: McKenzie L. Stock, Paul Gabriel, Jennifer J. Randall, Nicole Pietrasiak 2023. Soil DNA Extraction Modified Protocol for Dryland Agroecosystems. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2ly7domlx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 20, 2022



Last Modified: August 21, 2023

Protocol Integer ID: 70247

Keywords: soil dna extraction modified protocol for dryland agroecosystem, soil dna extraction modified protocol, soil dna extraction, dryland soil, agricultural soil, extracted dna, quality dna, dna extraction, high soil salinity, qiagen dneasy powerlyzer powersoil kit, dryland agroecosystem, dna, qiagen manufacturer protocol, quality checks for next generation amplicon

Funders Acknowledgements:

New Mexico Department of Agriculture Specialty Crop Block Grant

Grant ID: AM#21SCBPNM1027

Abstract

The Qiagen DNEasy PowerLyzer PowerSoil Kit is widely used in DNA based sample processing and extraction procedures. However, dryland soils are typically high in salts and secondary metabolites, which can cause interferences with A 260/230 quality of extracted DNA. We modified the Qiagen manufacturer protocol to account for high soil salinity by including additional washing steps and extra centrifugation time. Quality DNA was extracted from agricultural soils using this protocol which passed the quality checks for next generation amplicon sequencing.



Troubleshooting

Sample Prep

- 1 Soil samples should be stored cold, at least -20 C (ideally at -80 C).

Homogenization

3m

- 2 Fill a **12 ml homogenization vial** approximately 2/3 full with soil (**Spex SamplePrep 6133PC-T**). Include three 6.35 mm diameter chrome steel beads (**BioSpec Products Cat. No. 11079635c**).
- 3 Homogenize at 4000 rpm for  00:00:10 using the **SPEX SamplePrep 1200 Genolyte** homogenizer with the **single-vial attachment for 12mL vials**. Re-homogenize as needed at 4000 rpm for  00:00:10 to avoid the soil from thawing and sticking together.

20s

DNA Extraction

17m 10s

- 4 This protocol is modified from the Qiagen manufacturer protocol for the **DNeasy PowerLyzer PowerSoil Kit**.

Equipment

DNeasy PowerLyzer PowerSoil Kit

NAME

DNA Extraction Kit

TYPE

Qiagen

BRAND













12855-50

SKU














<https://www.qiagen.com/us/products/discovery-and-translational-research/dna-rna-purification/dna-purification/microbial-dna/dneasy-powerlyzer-powersoil-kit/>

LINK



- 5 Add approximately  80 μL of soil sample by volume to the PowerBead Tube provided.
- 6 Add  750 μL of PowerBead Solution to the PowerBead Tube.
- 7 Add  60 μL of Solution C1 and invert several times or vortex briefly.
- 8 Bead beating: Spex 1200 Genolyte homogenizer: Place the PowerLyzer Glass Bead Tubes into the **tube holder attachment for the homogenizer (holds 6)**. The PowerBead Tubes must be balanced on the tube holder. Run the samples at a time and RPM suitable for your soil type. For our samples, 3000 rpm for  00:00:35 . 35s
- 9 Centrifuge Bead Tubes at 10,000 x rcf for  00:03:00 . Note: We increased the centrifugation time to ensure that no sediment/soil is left in suspension, improving DNA quality downstream. 3m
- 10 Transfer the supernatant to a clean 2 ml collection tube. Note: Expect  400-500 μL . Supernatant may still contain some soil particles.
- 11 Add  250 μL of Solution C2 and vortex for 5 sec. Incubate at 2-8°C for 5 min.
- 12 Centrifuge the tubes at 10,000 x rcf for 1 min. Avoiding the pellet, transfer up to  600 μL of supernatant to a clean 2 ml collection tube.
- 13 Add  200 μL of Solution C3 and vortex briefly. Incubate at 2-8°C for  00:05:00 . 5m
- 14 Centrifuge the tubes at 10,000 x rcf for  00:03:00 . Avoiding the pellet, transfer up to  750 μL of supernatant to a clean 2 ml collection tube. Note: We increased the centrifugation time to ensure that no sediment/soil is left in suspension, improving DNA quality downstream. 3m



- 15 Add  1200 μL of Solution C4 to the supernatant and vortex for  00:00:05 . 5s
- 16 Load  675 μL of the supernatant onto a MB Spin Column and centrifuge at 10,000 x rcf for  00:01:00 . Discard the flow through and add an additional  675 μL of supernatant. 1m
- 17 Centrifuge at 10,000 x rcf for  00:01:00 . Load the remaining supernatant onto the MB Spin Column and centrifuge at 10,000 rcf for  00:01:00 . Note: A total of three loads for each sample processed are required. 2m
- 18 Add 500 μL of Solution C5 and centrifuge at 10,000 x rcf for  00:00:30 . Discard the flow through. Repeat this step once more: add  500 μL of Solution C5 and centrifuge at 10,000 x rcf for  00:00:30 . Discard the flow through. Note: We repeated this washing step to minimize salt contamination, which helps improve DNA A260/230 values. 1m
- 19 Centrifuge again at 10,000 x rcf for  00:01:00 . 1m
- 20 Carefully place the spin filter in a clean 2 ml collection tube. Avoid splashing any Solution C5 onto the MB Spin Column.
- 21 Add  50 μL of Solution C6 to the center of the white filter membrane. Note: We used the minimum recommended amount of Solution C6 to yield greatest DNA concentrations.
- 22 Centrifuge at 10,000 x rcf for  00:00:30 . Discard the MB Spin Column. 30s
- 23 The DNA is now ready for downstream applications.