



Feb 01, 2024

Solid phase binding assay - Clusterin binding to Very Low-Density Lipoprotein Receptor (VLDLR)

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm36kol3p/v1>

Patricia Yuste-Checa¹, Andreas Bracher¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm36kol3p/v1>

Protocol Citation: Patricia Yuste-Checa, Andreas Bracher, F Ulrich Hartl 2024. Solid phase binding assay - Clusterin binding to Very Low-Density Lipoprotein Receptor (VLDLR). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm36kol3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 24, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 94536

Keywords: ASAPCRN, density lipoprotein receptor, lipoprotein receptor, clusterin, vldlr, binding assay, linked immunosorbent assay

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol details how to monitor Clusterin binding to the Very Low-Density Lipoprotein Receptor (VLDLR) by Enzyme-linked immunosorbent assays (ELISA) adapted from Leeb et al. (2014).

Attachments



Solid phase binding ...

21KB

Materials

Buffers

- TBS-C: Tris-Buffered Saline $\text{pH } 7.4$, $[M] 2 \text{ millimolar (mM)}$ CaCl_2
- Blocking solution: 2% BSA, 0.05% Tween in TBS-C buffer.
- Quenching solution: $[M] 2 \text{ Mass Percent}$ sulfuric acid

 Clusterin- α Antibody (B-5) **Santa Cruz Biotechnology Catalog #sc-5289**

 RAP Antibody (E-7) **Santa Cruz Biotechnology Catalog #sc-515625**

 1-Step[®] Ultra TMB-ELISA Substrate Solution **Thermo Fisher Catalog #34028**



Solid phase binding assay - Clusterin binding to Very Low-Density Lipoprotein Receptor (VLDLR)

13h

- 1 Coat the corresponding wells of a 96-well plate (Nunc-Immuno MicroWell 96 well solid plate, MERCK) with  100 μL of TBS-C containing  10 μL VLDLR ectodomain
 Overnight at  4 $^{\circ}\text{C}$.



Note

The same number of wells should be incubated with TBS-C without VLDLR. These wells will be used as ligand background binding (addition of ligand to wells without immobilized receptor).

- 2 Wash the plate once with TBS-C.



Note

The washing step should be quick to avoid dilution and detachment of the receptor.

- 3 Add Blocking solution and incubate the plate for  02:00:00 at  Room temperature (25 $^{\circ}\text{C}$). The wells without receptor are now coated with BSA.

2h



- 4 Remove Blocking solution and apply a series of increasing concentrations of ligand diluted in Blocking solution, each in a final volume of  100 μL . Each ligand concentration should be added to one well with immobilized VLDLR and one well coated with BSA (Blocking solution) for ligand background binding. One well with VLDLR and one well coated with BSA should be incubated without ligand to determine the general plate background signal.

Note

- For Clusterin, a concentration range from 50 nM to 10000 nM is recommended (approximate K_D = 80-140 nM).
- Low Density Lipoprotein-Related Protein-Associated Protein 1 (LRPAP1 or RAP) is a molecular chaperone for LDL receptor-related proteins and therefore it can be used as positive control and as a competitor binder. For RAP binding, a concentration range from 1 nM to 60 nM is recommended (approximate K_D = 1-2 nM). For competition assays, mix a fixed concentration of the ligand with increasing concentrations of the competitor (RAP).



5 Incubate  01:00:00 at  Room temperature ( 25 °C).

1h



6 Wash the plate three times with Blocking solution.



Note

If testing the effect of pH on ligand binding, wash the wells once with TBS-C Blocking solution ( 7.4) or a low pH buffer like SA-C Blocking solution (

 10 millimolar (mM) Na-acetate  5.2 ,  150 millimolar (mM) NaCl,

 3 millimolar (mM) CaCl_2 , 2% BSA, 0.05% Tween) and incubate with the

corresponding buffers for  01:00:00 at  Room temperature ( 25 °C). After the incubation time, wash the plate once with the same buffers.

7 Add the corresponding primary antibodies diluted 1/100 in Blocking solution and incubate  01:00:00 at  Room temperature ( 25 °C).

1h



Note

anti-Clusterin (sc-5289 Santa Cruz Biotechnologies) and anti-RAP (sc-515625 Santa Cruz Biotechnologies) can be used for Clusterin and RAP detection, respectively.

8 Wash the plate three times with TBS-C Blocking solution.



9 Add the corresponding secondary antibody (horseradish peroxidase (HRP) conjugated) diluted 1/10,000 in Blocking solution and incubate  01:00:00 at  Room temperature ( 25 °C).

1h



10 Wash the plate three times with TBS-C Blocking solution.





- 11 Add  100 μ L per well of the HRP substrate 1-Step Ultra TMB ELISA Substrate Solutions (Thermo Fisher Scientific, 34028) to develop the plate and incubate until the desired color develops.



Note

3 and 10 minutes incubation time with the developing solution are normally enough under these conditions to develop RAP and Clusterin signal, respectively.

- 12 Add  100 μ L per well of quenching solution to stop the reaction.



- 13 Measure absorbance at 450 nm.

Note

First, subtract the background signal of each sample (VLDLR coated well – BSA coated well). Next, subtract plate background from each sample (wells incubated without ligand).

Protocol references

Leeb C, Eresheim C, Nimpf J. Clusterin is a ligand for apolipoprotein E receptor 2 (ApoER2) and very low density lipoprotein receptor (VLDLR) and signals via the Reelin-signaling pathway. *J Biol Chem*. 2014 Feb 14;289(7):4161-72. doi: 10.1074/jbc.M113.529271. Epub 2013 Dec 31. PMID: 24381170; PMCID: PMC3924281.