



Nov 21, 2020

## SOLUTION- 02 - Phosphate Buffered Saline (PBS)

DOI

[dx.doi.org/10.17504/protocols.io.bpximpke](https://dx.doi.org/10.17504/protocols.io.bpximpke)

Marco Cosentino<sup>1</sup>, Elisa Storelli<sup>1</sup>, Alessandra Luini<sup>1</sup>, Massimiliano LM Legnaro<sup>1</sup>, Emanuela Rasini<sup>1</sup>, Marco Ferrari<sup>1</sup>, Franca Marino<sup>1</sup>

<sup>1</sup>Center for Research in Medical Pharmacology, University of Insubria (Varese)



Farmacologia Medica

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bpximpke](https://dx.doi.org/10.17504/protocols.io.bpximpke)

**Document Citation:** Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino 2020. SOLUTION- 02 - Phosphate Buffered Saline (PBS). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bpximpke>

**License:** This is an open access document distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Created:** November 21, 2020

**Last Modified:** November 21, 2020

**Document Integer ID:** 44746



## Abstract

This recipe is used in the following protocols:

- Separation and purification of human PBMC from FRESH BLOOD
- Separation and purification of human PBMC from BUFFY COAT
- Magnetic bead-based CD4+ T cell isolation from PBMCs with Dynabeads: CD4 Positive Isolation Kit
- Magnetic bead-based TREG-TEFF cell isolation from PBMC with Miltenyi CD4+CD25+ Regulatory T cell Isolation Kit
- Staining of human PBMC or ISOLATED SUBSETS with Cell Proliferation Dye-eFluor™ 670 (CPD-eFluor670) for cell proliferation evaluation by Flow Cytometry



### **Phosphate-buffered saline (PSB 1X)**

**NaCl**- 8 g (137 mM)

**KCl**- 1.5 g (2 mM)

**Na<sub>2</sub>HPO<sub>4</sub>**- 1.44 g (10 mM)

**KH<sub>2</sub>PO<sub>4</sub>**- 0.24 g (1.8 mM)

To prepare 1 L of PBS 1X, dissolve the reagents listed above in 800 mL of ultrapure H<sub>2</sub>O.

Adjust the pH to 7.4 with HCl or NaOH, and then add H<sub>2</sub>O to 1 L.

**NaCl**code: S9625, Sigma

**KCl**code: P9541, Sigma

**Na<sub>2</sub>HPO<sub>4</sub>**code: 1.06585, Sigma

**KH<sub>2</sub>PO<sub>4</sub>**code: P0662, Sigma