



Nov 21, 2020

SOLUTION- 02 - Phosphate Buffered Saline (PBS)

DOI

dx.doi.org/10.17504/protocols.io.bpximpke

Marco Cosentino¹, Elisa Storelli¹, Alessandra Luini¹, Massimiliano LM Legnaro¹, Emanuela Rasini¹, Marco Ferrari¹, Franca Marino¹

¹Center for Research in Medical Pharmacology, University of Insubria (Varese)



Farmacologia Medica

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bpximpke>

Document Citation: Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino 2020. SOLUTION- 02 - Phosphate Buffered Saline (PBS). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bpximpke>

License: This is an open access document distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: November 21, 2020



Last Modified: November 21, 2020

Document Integer ID: 44746

Keywords: purification of human pbmc, teff cell isolation from pbmc, human pbmc, teff cell isolation, flow cytometry, pbmc, cell proliferation evaluation, phosphate buffered saline, cd4 positive isolation kit, isolated subsets with cell proliferation dye, dynabeads cd4 positive isolation kit, cell proliferation dye, fresh blood

Abstract

This recipe is used in the following protocols:

- Separation and purification of human PBMC from FRESH BLOOD
- Separation and purification of human PBMC from BUFFY COAT
- Magnetic bead-based CD4+ T cell isolation from PBMCs with Dynabeads: CD4 Positive Isolation Kit
- Magnetic bead-based TREG-TEFF cell isolation from PBMC with Miltenyi CD4+CD25+ Regulatory T cell Isolation Kit
- Staining of human PBMC or ISOLATED SUBSETS with Cell Proliferation Dye-eFluor™ 670 (CPD-eFluor670) for cell proliferation evaluation by Flow Cytometry

Troubleshooting



Phosphate-buffered saline (PSB 1X)

NaCl- 8 g (137 mM)

KCl- 1.5 g (2 mM)

Na₂HPO₄- 1.44 g (10 mM)

KH₂PO₄- 0.24 g (1.8 mM)

To prepare 1 L of PBS 1X, dissolve the reagents listed above in 800 mL of ultrapure H₂O.
Adjust the pH to 7.4 with HCl or NaOH, and then add H₂O to 1 L.

NaClcode: S9625, Sigma

KClcode: P9541, Sigma

Na₂HPO₄code: 1.06585, Sigma

KH₂PO₄code: P0662, Sigma