



Jul 31, 2020

SOLUTION- 11 - CD4+T cells isolation buffer

DOI

dx.doi.org/10.17504/protocols.io.bi79khr6

Marco Cosentino¹, Elisa Storelli¹, Alessandra Luini¹, Massimiliano LM Legnaro¹, Emanuela Rasini¹, Marco Ferrari¹, Franca Marino¹

¹Center for Research in Medical Pharmacology, University of Insubria (Varese, Italy)



Farmacologia Medica

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bi79khr6

Document Citation: Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino 2020. SOLUTION- 11 - CD4+T cells isolation buffer. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bi79khr6>

License: This is an open access document distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: July 31, 2020

Last Modified: July 31, 2020

Document Integer ID: 39905



Composition g/ml:

BSA 0.250 g

EDTA 0.2235 g

Add PBS(Ca²⁺ and Mg²⁺ free) to 240 ml, adjust to pH7.4 with NaOH or HCl, bring to 250 ml with PBS.

The substances are in the closet 1 and in refrigerator 1 (Room TSO8), weigh the substances with Gibertini balance (Room TSO8).

Filter this buffer by using a sterile syringe and filter.

Sterile filters and syringes are in closet 1 (Room PSO3).

Storage:  4 °C refrigerator 2 - (Room PSO3).

BSA: catalog number A2153 , Sigma, Italy

EDTA: catalog number ED2SS, Sigma, Italy

Prymo siringe filter: catalog number EPSPE2230 Euroclone, Italy.