



Feb 02, 2024

SOP for Immune cells isolation mouse adult brain with CD45 beads

DOI

dx.doi.org/10.17504/protocols.io.n92ldmxxnl5b/v1

Malu G Tansey¹

¹College of Medicine | University of Florida



Senthilkumar Karuppagounder

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldmxxnl5b/v1>

Protocol Citation: Malu G Tansey 2024. SOP for Immune cells isolation mouse adult brain with CD45 beads. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ldmxxnl5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 30, 2024



Last Modified: May 31, 2024

Protocol Integer ID: 94359

Keywords: ASAPCRN, sop for immune cells isolation mouse adult brain, immune cells isolation mouse adult brain, cd45 bead, immune cell, cd45, sop, cell, brain

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) Collaborative Research Network

Abstract

SOP for Immune cells isolation mouse adult brain with CD45 beads

Materials

Material:

Miltenyi GMHO machine (ref 130-096-427).

Adult brain dissociation kit mouse and rat (ref 130-107-677; 50 digestions).

GentleMACS C tubes (ref 130-093-271 for 25 tubes or 130-096-334 for 100 tubes).

D-PBS with calcium, magnesium, glucose, and pyruvate. (ref 14-287-080)

PB buffer (always prepare fresh): Prepare a solution containing D-PBS, pH 7.2, and

0.5% bovine serum albumin (BSA) by diluting MACS BSA Stock Solution (# 130-091-376)

1:20 with D-PBS (+/+).

70 um filters.

Small petri dishes and scalpels to cut mouse brains.

Double-distilled water (ddH₂O) or MilliQ.

MACS MultiStand and magnetic field.

LS columns (ref 130-042-401 for 25 columns).

Columns' size depends on the amount of cells (MS 2×10^8 total cells and 10^7 expected

CD11b+ cells; LS 2×10^9 total

cells and 10^8 expected CD11b+ cells)

Miltenyi mouse CD45 beads (#ref 130-052-301).

Beads Buffer (always prepare fresh): MACS BSA Stock Solution (#130-091-376) 1:20 with AutoMacs Rinsing Solution (#130-091-222).

Ice Buckets

Troubleshooting



Before:

- 1 EVERYTHING DONE ON ICE!!!
- 2 Label all the tubes and make sure all the material is available.
- 3 Make sure all the buffers are prepared and cold!! Can be done the day before, or day of.
- 4 Set up centrifuge at 4°C (use centrifuge on the right; close to the external entrance of the lab).
- 5 Prepare the hood and check if there is enough material in the tissue culture room.
- 6 Prepare enzyme mixes and store in the fridge:
*Tip: buffer Z can be added first in C-tubes-then start cutting brains-add them in C-tube-add enzyme P.
Enzyme mix 1: 50 ul of enzyme P and 1900 ul of buffer Z per sample/C tube.*
Enzyme mix 2: 10 ul of enzyme A and 20 ul of buffer Y per sample/C tube.

A. Brain dissociation

- 7 Remove mice brains and put them in pre-prepared 15mL falcon tubes filled with 3mL cold D-PBS.
- 8 Add 1900 ul of buffer Z in each C-tube.
- 9 Cut brain in ~8 equal pieces and add them in C-tubes (do not cut too much or viability will be affected).
- 10 Add 50 ul of enzyme P in each C-tube.
- 11 Add 30 ul of enzyme mix 2 in each C-tube.



- 12 Close lids tightly until lid clicks into place
- 13 Place C-tubes in GMHO machine-make sure all the brain pieces are at the bottom of the tubes; otherwise, they will not be properly dissociated.
- 14 Start program: 37_ABDK_01 for 1 brain >100 mg; 37_ABDK_02 for smaller pieces 20-100 mg

Tip: if any piece is stuck in the upper part, stop the program, mix brain pieces and start the program again.
- 15 Centrifuge briefly or mix up and down manually to ensure all the sample is at the bottom of the tube.
- 16 Prepare 50 ml tube with 70 um filters. Tip: add 2ml of PBS to the filter before.
- 17 Add the cell suspension into the C-tube.
- 18 Add 10 ml of cold D-PBS to the C-tubes, close the C-tubes, shake gently, and apply the cell suspension into a 50 ml falcon with a 70 um filter.
- 19 Add an extra 10 ml of D-PBS in the C-tube to make sure everything is taken and transfer again to the 50 ml falcon tube with the 70um filter. Optional: Use an extra 10 ml of D-PBS to dissociate small pieces that remain in the filter as much as possible.
- 20 Centrifuge 300xg for 10 min at 4°C. B. Debris removal

B. Debris removal

- 21 Remove supernatant, resuspend the cell pellet with 3100 ul of D-PBS and add 900 ul of debris removal solution.
- 22 Transfer cell solution to 15 ml tubes. If there are a lot of samples to one at a time. Avoid having the samples with the debris removal solution sitting on ice for long periods of time.
- 23 Overlay with 4 ml of D-PBS (debris removal solution is critical to have the overlay). The overlay volumes vary when using 1 or 2 brains:



	Debris Removal Solution	D-PBS	Overlay (D-PBS)	Reagent tube
20–100 mg	450 μ L	1550 μ L	2 mL	5 mL
400–500 mg (~ 1 brain)	900 μ L	3100 μ L	4 mL	15 mL
800–1000 mg (~ 2 brains)	1800 μ L	6200 μ L	4 mL	15 mL

- 24 Centrifuge 3000xg for 10 min at 4°C with full acceleration and brake of 7. Tip1: use the centrifuge on the right, close to the external entrance of the lab. Place the tubes on the inside of the centrifuge, the positions close to the rotor. Tip2: if after the centrifugation the debris layer is in the bottom of the tube with the cells, remove the upper layers, add again 900 μ L of debris removal solution, overlay with respective D-PBS, and centrifuge again. Optional: If centrifuges give suboptimal centrifugation, the acceleration and brake can be reduced.
- 25 Three phases should be formed; aspirate the two top phases and discard them.
- 26 Add D-PBS up to 15 ml and gently invert the tube three times- do not vortex!
- 27 Centrifuge 1000xg for 10 min at 4°C with full acceleration and full brake. C. Red blood cell removal

C. Red blood cell removal

- 28 Prepare red blood cell removal (RBCR; prepare fresh every time). RBCR: 1:10 RBCR solution (10x) with ddH₂O (i.e., 1 ml of RBCR with 9ml of cold ddH₂O).
- 29 Resuspend the cell pellet with the appropriate RBCR volume:
- 30 ml for 100–1000 mg (one mouse brain).
- 31 5 ml for 20–100 mg.



- 32 Incubate the cells for 10 min at 2°C -8°C (fridge at 4°C).
- 33 Add PB buffer:
- 34 ml for 100-1000 mg (one mouse brain)
- 35 ml for 20-100 mg
- 36 Centrifuge 300xg for 10 min at 4°C. If need to count cells: resuspend in 10 ml of Beads Buffer◊ take an aliquote to count cells◊ centrifuge 300g 10 min at 4°C. D. CD45 beads isolation

D. CD45 beads isolation

- 37 Prepare Eppendorf tubes to incubate cells with CD45 microbeads.
- 38 Resuspend the cell pellet in 90 ul of Beads Buffer and 10 ul of CD45 microbeads for 107 cells Tip: cells should be counted to resuspend the cell pellet with appropriate vol of beads buffer and beads; however, we never expect more than 107 total cells per brain, thus directly resuspend the pellet with min vol to speed up the protocol.
- 39 Incubate the cells in the fridge (4°C) for 15 min. Tip: mix up and down after 7 min.
- 40 Add 2 ml of Beads Buffer (for 107 total cells) to each sample and centrifuge 300xg for 10 min. Tip: no temperature recommended, but I would use at 4°C.
- 41 Prepare the MACS MultiStand and magnetic field under the hood.
- 42 Place the columns (make sure they are correctly placed) and 15 ml tubes for trash.
- 43 Pre-wash columns with 500 ul Beads Buffer (MS columns).



- 44 Resuspend cell pellet with 500 ul of Beads Buffer (up to 108 cells). Tip: make sure the columns don't dry from now on.
- 45 Transfer the cell suspension to the columns (CD45+ cells are attached in the magnetic part of the column).
- 46 Wash the columns x3 times with 500 ul of Beads Buffer, add buffer only when column is empty each time.
- 47 Place columns in new 15 ml tubes.
- 48 Add 1 ml (MS) or 5 ml (LS) of Beads Buffer to each column and flush out the columns by pushing the plunger into the column.
- 49 Take an aliquot to count cells.
- 50 Centrifuge 300xg for 10 min at 4°C.
- 51 Resuspend according to the follow-up experiment.

Protocol references

SOP adapted from Miltenyi's instructions