



Oct 02, 2024

SOP16v1_TGD_Qubit_dsDNA_BR_Assays_for_DNAQuantification

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lyppqelx9/v1>

Varsha Rajesh¹

¹Stanford University



Varsha Rajesh

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyppqelx9/v1>

Protocol Citation: Varsha Rajesh 2024. SOP16v1_TGD_Qubit_dsDNA_BR_Assays_for_DNAQuantification. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2lyppqelx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 23, 2020

Last Modified: October 02, 2024

Protocol Integer ID: 43571

Keywords: dna bench, blue dna, pcr, dna, lab attire, buffer mix

Abstract

Guidelines

Bench work should be completed on the **pre-PCR/DNA bench**.

Blue DNA-only labeled materials should be used for this protocol to avoid cross-contamination.

Lab attire: gloves, lab coat, and safety glasses/goggles.

Before Start

Grab ice.

Use this spreadsheet template to calculate reagent + buffer mix.



Qubit Prep Template - dsDNA BR A...

Guidelines

Bench work should be completed on the **pre-PCR/DNA bench**.

Blue DNA-only labeled materials should be used for this protocol to avoid cross-contamination.

Lab attire: gloves, lab coat, and safety glasses/goggles.

Materials

	A	B	C
	Material	Catalog	Supplier
	Qubit Assay Tubes	Q32856	Life Technologies
	Qubit dsDNA BR Assay Reagent Kit	Q32853	Life Technologies
	Qubit Fluorometer	(Snyder Lab, Rm 2265)	Life Technologies

Safety warnings


 Waste Disposal: all components can be thrown in the normal trash.

Before start

Grab ice.

Use this spreadsheet template to calculate reagent + buffer mix.


 Qubit Prep Template - dsDNA BR A...

- 1 Prepare the correct amount of mastermix of Qubit reagent in buffer for your samples and standards.

Note: A calibration is required for every usage of the Qubit, so always include standards in the calculations.

Reagent and buffer mastermix is light sensitive. Wrap the tube in foil to protect from light.

- 1.1 Reagent must be diluted 1:200 in buffer.
Each standard requires 190 uL of mastermix.
Each sample requires 190-199 uL of mastermix. (You can put anywhere from 1 to 10 uL of sample in each qubit assay).

- 2 Select the reagent calculator on the Qubit Fluorometer main menu.
Ensure that your calculations match up with the numbers calculated by the Qubit for mastermix preparation.
- 3 Prepare samples: Aliquot the appropriate amount of mastermix (190-199 uL) into each 0.5 mL qubit assay tube. Add the appropriate amount of sample (1-10 uL) into each tube, and let sit for 2 minutes to incubate.

After incubation, the sample fluorescence is stable for 3 hours.

- 4 Prepare Standards 1 and 2. Each requires 10 uL of standard and 190 uL of mastermix. Keep the standards on ice.
- 5 Go back to home screen and select the dsDNA assay: Broad Range. Select "Read Standards".
- 6 Input Standard 1, then read. Then input Standard 2. Generate calibration curve.
- 7 Enter what sample volume you added to each tube.
- 8 After reading one tube, just take out the sample and put in the next tube. This way all samples will be saved under the same assay data set.
- 9 When finished, upload the data onto a USB.