

Oct 09, 2024

SOP52v1_TGD_IIDP-HIGISampleProcessing_Public

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg32joev25/v1>

Varsha Rajesh¹, Swaraj Thaman¹

¹Stanford University



Swaraj Thaman

LSRP1, Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg32joev25/v1>

Protocol Citation: Varsha Rajesh, Swaraj Thaman 2024. SOP52v1_TGD_IIDP-HIGISampleProcessing_Public. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg32joev25/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 08, 2024



Last Modified: October 09, 2024

Protocol Integer ID: 109320

Keywords: translational genomics of diabetes lab, integrated islet distribution program, translational genomic, diabetes lab, pancreatic tissue, diabetes, higisampleprocessing-public this protocol detail, genotyping

Abstract

This protocol details the workflow to process pancreatic tissues (received from the Integrated Islet Distribution Program's distribution centers) for genotyping at the Translational Genomics of Diabetes Lab at Stanford University.

Receiving Samples

- 1 When acinar samples are received from the various isolation centers, match what is physically sent with sample info list sent with package and make sure all samples are accounted for. Take note of the condition of the samples.
- 2 Let isolation center know samples were received and in what condition.
- 3 Print QR barcodes for the acinar samples in this format and tape onto the tubes.

HIGI#ACIN (# being the next sequential number)
- 4 Store tubes in Yalow -80°C freezer in "IIDP - HIGI Acinar Samples" box.

Extracting DNA

- 5 Thaw acinar tubes on ice.
- 6 Use DNA extraction protocol for DNeasy Blood & Tissue Kit (<https://protocols.io/view/dna-extraction-qiagen-dneasy-blood-amp-tissue-kit-bjngkmbw>).

QC

- 7 After extraction, measure the concentration by nanodrop, taking account of A260/280 and A260/230 values as well. A260/280 value should be close to 1.8 and A260/230 should be around 1.8-2.0.
- 8 Take a portion of the sample to measure concentration by Qubit. Protocol can be found here: <https://protocols.io/view/qubit-dsdna-assays-bnstmeen>. For genotyping, Qubit concentration should be above 50 ng/uL and total DNA should be above 750 ng. In case the DNA concentration is below 50 ng/uL repeat the DNA extraction with a fresh piece of acinar tissue from the same donor and QC the extracted DNA. Keep repeating this till you are able to extract at least 750ng of DNA for that donor (it is okay if the concentration of the extracted DNA is still less than 50ng/uL) - for these samples, when you submit them for genotyping (step 10), make sure that you pool together and submit all the DNA you have extracted for that donor and request for the sample to be 'speed-vac'ed at the SFGF core - this will help bring the concentration up to 50ng/uL, while still making sure that there is 750ng of DNA, thereby meeting the QC requirements for genotyping.



- 9 Create QR codes that include sample name in this format: HIGI#DNA. Label should also include nanodrop concentration and QC info.

Genotyping

- 10 Record sample info, nanodrop QC and Qubit information in genotyping QC spreadsheet, which will be sent to the SFGF core whenever the entire genotyping batch is ready (follow this protocol: <https://www.protocols.io/edit/sop54v1-tgd-genotypingsubmissionforcollaborators-b8jaruie>).
- 11 As of 2023-12-01, samples were sent to be genotyped on Illumina's Infinium Omni Exome 2.5 arrays.
As of 2024-10-31, samples were sent to be genotyped on Illumina's Infinium Global Diversity Array-8+ v1.0 Omni2.5-Exome-8v1.6 arrays.
- 12 Following the QC of the genotyping data, if the data for a donor fails QC, resend the DNA from that same donor on another genotyping batch.