


Apr 20, 2024

Spatial Localization of 3D Scanned EEG Electrodes with MeshLab

DOI

<https://dx.doi.org/10.17504/protocols.io.261geo8zjl47/v1>



Joel P Diaz-Fong^{1,2}, Jordana Glouberman¹, Agatha Lenartowicz¹

¹Semel Institute for Neuroscience and Human Behavior, University of California Los Angeles;

²Institute of Medical Science, University of Toronto



Joel P Diaz-Fong

Brain Health Imaging Centre at Centre for Addiction and Ment...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261geo8zjl47/v1>

Protocol Citation: Joel P Diaz-Fong, Jordana Glouberman, Agatha Lenartowicz 2024. Spatial Localization of 3D Scanned EEG Electrodes with MeshLab. [protocols.io https://dx.doi.org/10.17504/protocols.io.261geo8zjl47/v1](https://dx.doi.org/10.17504/protocols.io.261geo8zjl47/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 09, 2020

Last Modified: April 20, 2024

Protocol Integer ID: 36806

Keywords: digitization, electrode localization, cognitive neuroscience, spatial localization of 3d scanned eeg electrode, 3d scanned eeg electrode, spatial localization of electroencephalogram, electrodes from 3d, utilizing meshlab, meshlab, electroencephalogram, pickpoints tool, meshlab this protocol, mesh model, electrode, lidar scanner, head model, eeg, spatial localization, 3d

Funders Acknowledgements:

National Institute of Mental Health

Grant ID: R01MH116268

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol outlines step-by-step instructions for utilizing Meshlab's *PickPoints* tool to facilitate spatial localization of electroencephalogram (EEG) electrodes from 3D scanned head models. These include mesh models created from structured-light/infrared, photogrammetry, and LiDAR scanners, to name a few.

Image Attribution

All images were generated using the MeshLab opensource software. The example 3D mesh model used in this protocol is a scanned mannequin head wearing an EGI HydroCel Geodesic Sensor Net (Electrical Geodesics Inc., Eugene, Oregon) that was custom ordered without the four facial electrodes (total of 125 electrodes). The 3D model was scanned using the Structure Sensor (Occipital Inc., Boulder, CO), a structured-light depth camera, mounted on an iPad 6 (Apple, Cupertino, CA). Data acquisition was performed with the itSeez3D software (Itseez3D, Inc., Santa Clara, CA).



Guidelines

Follow the Manufacturer's Instructions. Please refer to the user manual provided by the manufacturer of your 3D scanning device for instructions on how to obtain good quality scans with their device. The accuracy of the coordinates obtained using the method described in this protocol relies heavily on the quality of the 3D mesh model you provide. For details about the accuracy and precision of your device, please consult your device's technical specifications.

Materials

EEG, 3D mesh model, computer, meshlab software

Before start

Download the latest version of the MeshLab software using the installer found on their website (www.meshlab.net) or download their source code directly from GitHub (<https://github.com/cnr-isti-vclab/meshlab>). For 3D scanning sample data, please visit: <https://osf.io/87av2/>



3D Scanning Setup

10m

1 Participant Preparation

There are a few necessary steps to ensure that the participant is prepared in a way that guarantees a good quality scan. This includes:

1.1 **Fiducial Visibility**

- *Identify Key Fiducials.* Determine the anatomical landmarks (fiducials) to be used for coordinate registration (e.g., bridge of the nose, ears).
- *Ensure Unobstructed View.* Instruct participants to remove any accessories or items that may obstruct the visibility of fiducials during the scan. Also, verify that the bridge of the nose and ears are clearly visible in the final 3D mesh model.

1.2 **Electrode Visibility**

- *Check Electrode Placement.* Verify the correct placement of electrodes on the participant's scalp according to the EEG system used. Ensure that no electrodes are obstructed by hair or any other objects.
- *Tuck Away Hair.* Advise participants to ensure all strands of hair are neatly tucked away from the scanner's field of view of the electrodes.

1.3 **Reduce Reflective Surfaces**

- *Identify Reflective Surfaces.* Identify any reflective surfaces in the scanning environment that may impact electrodes visibility.
- *Apply Non-Reflective Coverings.* Use non-reflective stickers or a wax marker to cover reflective surfaces, especially those in close proximity to the electrodes.

Import Mesh Model

2m

2 Import the Mesh Model

After acquiring a 3D mesh model, open the *MeshLab* software and import the 3D mesh file.

Software

MeshLab

NAME

Visual Computing Lab - ISTI - CNR

DEVELOPER

<https://github.com/cnr-isti-vclab/meshlab>

REPOSITORY

<https://www.meshlab.net/>

SOURCE LINK

**Note**

In case you do not have a 3D model, please use the dataset below for practice.

Dataset

3D Scanning Sample Data for EEG Electrode Localization^{NAME}

<https://osf.io/87av2/>

LINK

2.1 Import mesh

- To do so, select *File > Import Mesh...*


Note

The *MeshLab* software can import several types of 3D triangular mesh files

- .obj, .ply, .stl, .3ds, etc.

2.2 Optional: Once the file is imported into the software, remove the trackball to get a better view of the mesh model.

- Select *View > Show Trackball*
- Unselect *Show Trackball* to remove the trackball



Imported mesh model with overlaid trackball

Spatial Localization

35m

3 Open the *PickPoints* Tool

2m

- Once the model is imported into the *MeshLab* software, you can begin to localize the fiducial and electrode coordinates with the *PickPoints* tool.
- To start, find the *PickPoints* icon in the toolbar or select *Edit > PickPoints*



PickPoints Icon

Note

PickPoints Window

When first opening the *PickPoints* tool, the window will initially be empty. The *Point Name* should correspond to the electrodes you are trying to localize.

The screenshot shows the 'Form' window of the PickPoints tool. It has a title bar with a close button. The interface is divided into several sections:

- Mode:** Three radio buttons: 'Pick Point' (selected), 'Move Point', and 'Select Point'.
- Buttons:** 'Load Points From File' and 'Save' are located to the right of the mode buttons.
- Table:** A table with five columns: 'Point Name', 'X', 'Y', 'Z', and 'active'. The table is currently empty.
- Actions:** A row of buttons below the table: 'Rename Point', 'Remove Point', 'Clear Point', 'Remove All Points', and 'Undo last move'.
- Template Controls:** A section with a checkbox 'Save this as your default template', and buttons 'Save', 'Load', 'Clear', and 'Add Point'.
- Template Name:** A text field showing 'No Template Loaded'.
- Normal Options:** A section with a checked checkbox 'Show Normal?', and a 'Draw as a:' section with 'Pin' (selected) and 'Line' radio buttons.

3.1 **Enter Point Names or load an existing template**

- Before you start selecting points on the mesh, manually enter each electrode label as a *Point Name* (see 5.2 about saving templates).
- If you have an existing template, go to *Template Controls*, select *Load* and choose the .pptpl template file.

3m

Note

PickPoints Window with Template


Point Name	X	Y	Z	active
NAS				☒
LHJ				☐
RHJ				☐
E1				☐
E2				☐
E3				☐
E4				☐
E5				☐
E6				☐
E7				☐
E8				☐
E9				☐
E10				☐
E11				☐
E12				☐
E13				☐

An example of a loaded template. The Point Name now includes the three fiducial labels and the electrode labels for EGI's HydroCell EEG net.

4 Begin selecting points on the mesh

- Locate the fiducial or electrode on the mesh model and **right-click** the location with your cursor. A green pin should appear at that location.
- To rotate the mesh model, unselect *PickPoints* and rotate/drag the model with your cursor. Once you've adjusted the model to the desired position, reselect *PickPoints* to continue.

30m



Green pin placed on Nasion (NAS).

Mode:					Load Points From File				
☒ Pick Point ☐ Move Point ☐ Select Point					Save				
Point Name	X	Y	Z	active					
NAS	0.0165417	0.123304	0.07327	☒					
LHJ	0	0	0	☐					
RHJ	0	0	0	☐					
E1	0	0	0	☐					
E2	0	0	0	☐					
E3	0	0	0	☐					
E4	0	0	0	☐					
E5	0	0	0	☐					

After right-clicking the model, XYZ coordinate points should appear next to the Point Name, and the checkbox in the active column will be selected. The next row will automatically be highlight.

4.1 **Selecting the fiducials**

- As mentioned in 1.1, the fiducials (anatomical landmarks) should be determined beforehand. The fiducials should also be the same across all your scans for a given

2m



study.

- For the example in this protocol, we select the nasion (NAS), the left helix-tragus junction (LHJ), and the right helix-tragus junction (RHJ).

4.2 **Selecting the electrodes**

28m

- Select each of the electrodes, point by point, following the order of the *Point Names* column.
- For a good approximation of the electrode location, try to select the point at the center of the electrode. The green pin should appear perpendicular to the scalp if the surface of the electrode is flat.

Save PickedPoints

2m

5 **Save electrode coordinates**


- Select *Save* at the upper right-hand corner of the *PickPoints* tool.

Note

Coordinate points are saved as an XML formatted picked point (.pp) file.

5.1 **Optional:** save a snapshot


- Snapshots of the mesh model with the green coordinate pins can be saved as an image file, which is useful for documenting the quality of mesh model and how the fiducials were defined for the selected coordinate points.
- Go to *File > Save snapshot*



Example snapshot of the mesh model with all the selected points. The model texture (or surface detail) can be removed for anonymity.

5.2 **Optional:** save the template of *Point Names* so the electrode labels can be loaded for processing multiple participants

- Under *Template Controls*, select *Save*. Templates are saved as a *.pptpl* filetype.



Extract Coordinates

5m

6 Finally, the coordinates should be extracted from the *PickedPoints* file and converted to a preferred format.

6.1 **Convert *PickedPoints* XML file to text**

Extract coordinates from XML file and convert them to text.



6.2 **Convert coordinates to a standard space**

The coordinates that are extracted from the *PickedPoints* file are typically the same unit of measurement and coordinate space as the mesh model. Use the fiducials as a reference point to convert the coordinate space to a standard format.

Protocol references

P. Cignoni, M. Callieri, M. Corsini, M. Dellepiane, F. Ganovelli, G. Ranzuglia. Cignoni, P., Callieri, M., Corsini, M., Dellepiane, M., Ganovelli, F., & Ranzuglia, G. (2008, July). Meshlab: an open-source mesh processing tool. In *Eurographics Italian chapter conference* (Vol. 2008, pp. 129-136). Sixth Eurographics Italian Chapter Conference, page 129-136, 2008