

Oct 22, 2024

spatial-Mux-seq

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g71n79gwz/v1>


Pengfei Guo¹, Yanxiang Deng¹

¹University of Pennsylvania



Pengfei Guo

University of Pennsylvania



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g71n79gwz/v1>

Protocol Citation: Pengfei Guo, Yanxiang Deng 2024. spatial-Mux-seq. protocols.io

<https://dx.doi.org/10.17504/protocols.io.q26g71n79gwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 25, 2024

Last Modified: October 22, 2024

Protocol Integer ID: 106407

Keywords: Spatial multi-omics, nano-CUT&Tag, co-profiling of five modalities, microfluidic in situ, open chromatin at the same time, including multiple epigenomic modality, open chromatin, multiple epigenomic modality, whole transcriptome, histone modification, cellular level, generation sequencing, mux, tissue scale, sequencing

Abstract

Spatial-Mux-seq is a multi-modal spatial technology that allows simultaneous profiling of up to five different modalities, including multiple epigenomic modalities (any combination of two histone modifications and open chromatin at the same time), whole transcriptome, and a panel of proteins at tissue scale and cellular level in a spatially resolved manner. This was achieved by integrating microfluidic *in situ* barcoding and the nanobody-tethered transposition chemistry directly in tissue followed by high-throughput Next-Generation Sequencing.

Materials

Preparation of buffers:

DNA oligos used for PCR and preparation of sequencing library.

	RT-primer	/5Phos/CATCGGCGTACGA CTNNNNNNNNNN/iBiodT/ TTTTTTTTTTTTTTTTVN
	template switch primer	AAGCAGTGGTATCAACGCA GAGTGAATrGrG+G
	Ligation linker 1	AGTCGTACGCCGATGCGA AACATCGGCCAC
	Ligation linker 2	CGAATGCTCTGGCCTCTC AAGCACGTGGAT
	PCR primer 1	CAAGCGTTGGCTTCTCGC ATCT
	PCR primer 2	AAGCAGTGGTATCAACGCA GAGT
	primer 3 (cite-seq)	CCTTGGCACCCGAGAATT *C*C
	P5 oligo (cite-seq)	AATGATACGGCGACCACC GAGATCTACACTAGATCGC TCGTCGGCAGCGTCAGAT



		GTGTATAAGAGACAGCCTT GGCACCCGAGAATTCC
	N501	AATGATACGGCGACCACC GAGATCTACACTAGATCGC TCGTCCGCAGCGTCAGAT GTGTATAAGAGACAG
	N701	CAAGCAGAAGACGGCATA CGAGATTTCGCCTTAGTCTC GTGGGCTCGGAGATGTGTA TAAGAGACAGCAAGCGTTG GCTTCTCGCATCT
	N702	CAAGCAGAAGACGGCATA CGAGATCTAGTACGGTCTC GTGGGCTCGGAGATGTGTA TAAGAGACAGCAAGCGTTG GCTTCTCGCATCT
	N703	CAAGCAGAAGACGGCATA CGAGATTTCTGCCTGTCTC GTGGGCTCGGAGATGTGTA TAAGAGACAGCAAGCGTTG GCTTCTCGCATCT
	N704	CAAGCAGAAGACGGCATA CGAGATGCTCAGGAGTCTC GTGGGCTCGGAGATGTGTA TAAGAGACAGCAAGCGTTG GCTTCTCGCATCT
	N705	CAAGCAGAAGACGGCATA CGAGATAGGAGTCCGTCTC GTGGGCTCGGAGATGTGTA TAAGAGACAGCAAGCGTTG GCTTCTCGCATCT
	N706	CAAGCAGAAGACGGCATA CGAGATCATGCCTAGTCTC GTGGGCTCGGAGATGTGTA TAAGAGACAGCAAGCGTTG GCTTCTCGCATCT
	N707	CAAGCAGAAGACGGCATA



		CGAGATGTAGAGAGGTCTC GTGGGCTCGGAGATGTGTA TAAGAGACAGCAAGCGTTG GCTTCTCGCATCT
	N709	CAAGCAGAAGACGGCATA CGAGATAGCGTAGCGTCTC GTGGGCTCGGAGATGTGTA TAAGAGACAGCAAGCGTTG GCTTCTCGCATCT
	N710	CAAGCAGAAGACGGCATA CGAGATCAGCCTCGGTCT CGTGGGCTCGGAGATGTG TATAAGAGACAGCAAGCGT TGGCTTCTCGCATCT
	N711	CAAGCAGAAGACGGCATA CGAGATTGCCTCTTGTCTC GTGGGCTCGGAGATGTGTA TAAGAGACAGCAAGCGTTG GCTTCTCGCATCT
	N712	CAAGCAGAAGACGGCATA CGAGATTCTCTACGTCTC GTGGGCTCGGAGATGTGTA TAAGAGACAGCAAGCGTTG GCTTCTCGCATCT

Supplementary Table 5.**Barcode A Sequence**

	Barcode A-1	/5Phos/AGGCCAGAGCATT CGAACGTGATGTGGCCGA TGTTTCG
	Barcode A-2	/5Phos/AGGCCAGAGCATT CGAAACATCGGTGGCCGA TGTTTCG
	Barcode A-3	/5Phos/AGGCCAGAGCATT CGATGCCTAAGTGGCCGAT



		GTTTCG
	Barcode A-4	/5Phos/AGGCCAGAGCATT CGAGTGGTCAGTGGCCGA TGTTTCG
	Barcode A-5	/5Phos/AGGCCAGAGCATT CGACCACTGTGTGGCCGA TGTTTCG
	Barcode A-6	/5Phos/AGGCCAGAGCATT CGACATTGGCGTGGCCGA TGTTTCG
	Barcode A-7	/5Phos/AGGCCAGAGCATT CGCAGATCTGGTGGCCGA TGTTTCG
	Barcode A-8	/5Phos/AGGCCAGAGCATT CGCATCAAGTGTGGCCGAT GTTTCG
	Barcode A-9	/5Phos/AGGCCAGAGCATT CGCGCTGATCGTGGCCGA TGTTTCG
	Barcode A-10	/5Phos/AGGCCAGAGCATT CGACAAGCTAGTGGCCGA TGTTTCG
	Barcode A-11	/5Phos/AGGCCAGAGCATT CGCTGTAGCCGTGGCCGA TGTTTCG
	Barcode A-12	/5Phos/AGGCCAGAGCATT CGAGTACAAGGTGGCCGA TGTTTCG
	Barcode A-13	/5Phos/AGGCCAGAGCATT CGAACAACCAAGTGGCCGA



		TGTTTCG
	Barcode A-14	/5Phos/AGGCCAGAGCATT CGAACCGAGAGTGGCCGA TGTTTCG
	Barcode A-15	/5Phos/AGGCCAGAGCATT CGAACGCTTAGTGGCCGA TGTTTCG
	Barcode A-16	/5Phos/AGGCCAGAGCATT CGAAGACGGAGTGGCCGA TGTTTCG
	Barcode A-17	/5Phos/AGGCCAGAGCATT CGAAGGTACAGTGGCCGA TGTTTCG
	Barcode A-18	/5Phos/AGGCCAGAGCATT CGACACAGAAGTGGCCGA TGTTTCG
	Barcode A-19	/5Phos/AGGCCAGAGCATT CGACAGCAGAGTGGCCGA TGTTTCG
	Barcode A-20	/5Phos/AGGCCAGAGCATT CGACCTCCAAGTGGCCGA TGTTTCG
	Barcode A-21	/5Phos/AGGCCAGAGCATT CGACGCTCGAGTGGCCGA TGTTTCG
	Barcode A-22	/5Phos/AGGCCAGAGCATT CGACGTATCAGTGGCCGAT GTTTCG
	Barcode A-23	/5Phos/AGGCCAGAGCATT CGACTATGCAGTGGCCGAT



		GTTTCG
	Barcode A-24	/5Phos/AGGCCAGAGCATT CGAGAGTCAAGTGGCCGA TGTTTCG
	Barcode A-25	/5Phos/AGGCCAGAGCATT CGAGATCGCAGTGGCCGA TGTTTCG
	Barcode A-26	/5Phos/AGGCCAGAGCATT CGAGCAGGAAGTGGCCGA TGTTTCG
	Barcode A-27	/5Phos/AGGCCAGAGCATT CGAGTCACTAGTGGCCGAT GTTTCG
	Barcode A-28	/5Phos/AGGCCAGAGCATT CGATCCTGTAGTGGCCGAT GTTTCG
	Barcode A-29	/5Phos/AGGCCAGAGCATT CGATTGAGGAGTGGCCGA TGTTTCG
	Barcode A-30	/5Phos/AGGCCAGAGCATT CGCAACCACAGTGGCCGA TGTTTCG
	Barcode A-31	/5Phos/AGGCCAGAGCATT CGGACTAGTAGTGGCCGAT GTTTCG
	Barcode A-32	/5Phos/AGGCCAGAGCATT CGCAATGGAAGTGGCCGA TGTTTCG
	Barcode A-33	/5Phos/AGGCCAGAGCATT CGCACTTCGAGTGGCCGA



		TGTTTCG
	Barcode A-34	/5Phos/AGGCCAGAGCATT CGCAGCGTTAGTGGCCGA TGTTTCG
	Barcode A-35	/5Phos/AGGCCAGAGCATT CGCATACCAAGTGGCCGAT GTTTCG
	Barcode A-36	/5Phos/AGGCCAGAGCATT CGCCAGTTTCACTGGCCGA TGTTTCG
	Barcode A-37	/5Phos/AGGCCAGAGCATT CGCCGAAGTAGTGGCCGA TGTTTCG
	Barcode A-38	/5Phos/AGGCCAGAGCATT CGCCGTGAGAGTGGCCGA TGTTTCG
	Barcode A-39	/5Phos/AGGCCAGAGCATT CGCCTCCTGAGTGGCCGA TGTTTCG
	Barcode A-40	/5Phos/AGGCCAGAGCATT CGCGAACTTAGTGGCCGA TGTTTCG
	Barcode A-41	/5Phos/AGGCCAGAGCATT CGCGACTGGAGTGGCCGA TGTTTCG
	Barcode A-42	/5Phos/AGGCCAGAGCATT CGCGCATACAGTGGCCGA TGTTTCG
	Barcode A-43	/5Phos/AGGCCAGAGCATT CGCTCAATGAGTGGCCGA



		TGTTTCG
	Barcode A-44	/5Phos/AGGCCAGAGCATT CGCTGAGCCAGTGGCCGA TGTTTCG
	Barcode A-45	/5Phos/AGGCCAGAGCATT CGCTGGCATA GTGGCCGA TGTTTCG
	Barcode A-46	/5Phos/AGGCCAGAGCATT CGGAATCTGAGTGGCCGA TGTTTCG
	Barcode A-47	/5Phos/AGGCCAGAGCATT CGCAAGACTAGTGGCCGA TGTTTCG
	Barcode A-48	/5Phos/AGGCCAGAGCATT CGGAGCTGAAGTGGCCGA TGTTTCG
	Barcode A-49	/5Phos/AGGCCAGAGCATT CGGATAGACAGTGGCCGAT GTTTCG
	Barcode A-50	/5Phos/AGGCCAGAGCATT CGGCCACATAGTGGCCGA TGTTTCG
	Barcode A-51	/5Phos/AGGCCAGAGCATT CGGCGAGTAAGTGGCCGA TGTTTCG
	Barcode A-52	/5Phos/AGGCCAGAGCATT CGGCTAACGAGTGGCCGA TGTTTCG
	Barcode A-53	/5Phos/AGGCCAGAGCATT CGGCTCGGTAGTGGCCGA



		TGTTTCG
	Barcode A-54	/5Phos/AGGCCAGAGCATT CGGGAGAACAGTGGCCGA TGTTTCG
	Barcode A-55	/5Phos/AGGCCAGAGCATT CGGGTGCGAAGTGGCCGA TGTTTCG
	Barcode A-56	/5Phos/AGGCCAGAGCATT CGGTACGCAAGTGGCCGA TGTTTCG
	Barcode A-57	/5Phos/AGGCCAGAGCATT CGGTCGTAGAGTGGCCGA TGTTTCG
	Barcode A-58	/5Phos/AGGCCAGAGCATT CGGTCTGTCAAGTGGCCGA TGTTTCG
	Barcode A-59	/5Phos/AGGCCAGAGCATT CGGTGTTCTAGTGGCCGAT GTTTCG
	Barcode A-60	/5Phos/AGGCCAGAGCATT CGTAGGATGAGTGGCCGAT GTTTCG
	Barcode A-61	/5Phos/AGGCCAGAGCATT CGTATCAGCAGTGGCCGAT GTTTCG
	Barcode A-62	/5Phos/AGGCCAGAGCATT CGTCCGTCTAGTGGCCGAT GTTTCG
	Barcode A-63	/5Phos/AGGCCAGAGCATT CGTCTTCACAGTGGCCGAT



		GTTTCG
	Barcode A-64	/5Phos/AGGCCAGAGCATT CGTGAAGAGAGTGGCCGA TGTTTCG
	Barcode A-65	/5Phos/AGGCCAGAGCATT CGTGGAACAAGTGGCCGA TGTTTCG
	Barcode A-66	/5Phos/AGGCCAGAGCATT CGTGGCTTCAGTGGCCGA TGTTTCG
	Barcode A-67	/5Phos/AGGCCAGAGCATT CGTGGTGGTAGTGGCCGA TGTTTCG
	Barcode A-68	/5Phos/AGGCCAGAGCATT CGTTCACGCAGTGGCCGA TGTTTCG
	Barcode A-69	/5Phos/AGGCCAGAGCATT CGAACTCACCGTGGCCGA TGTTTCG
	Barcode A-70	/5Phos/AGGCCAGAGCATT CGAAGAGATCGTGGCCGA TGTTTCG
	Barcode A-71	/5Phos/AGGCCAGAGCATT CGAAGGACACGTGGCCGA TGTTTCG
	Barcode A-72	/5Phos/AGGCCAGAGCATT CGAATCCGTCGTGGCCGA TGTTTCG
	Barcode A-73	/5Phos/AGGCCAGAGCATT CGAATGTTGCGTGGCCGA



		TGTTTCG
	Barcode A-74	/5Phos/AGGCCAGAGCATT CGACACGACCGTGGCCGA TGTTTCG
	Barcode A-75	/5Phos/AGGCCAGAGCATT CGACAGATTTCGTGGCCGAT GTTTCG
	Barcode A-76	/5Phos/AGGCCAGAGCATT CGAGATGTACGTGGCCGAT GTTTCG
	Barcode A-77	/5Phos/AGGCCAGAGCATT CGAGCACCTCGTGGCCGA TGTTTCG
	Barcode A-78	/5Phos/AGGCCAGAGCATT CGAGCCATGCGTGGCCGA TGTTTCG
	Barcode A-79	/5Phos/AGGCCAGAGCATT CGAGGCTAACGTGGCCGA TGTTTCG
	Barcode A-80	/5Phos/AGGCCAGAGCATT CGATAGCGACGTGGCCGA TGTTTCG
	Barcode A-81	/5Phos/AGGCCAGAGCATT CGATCATTCCGTGGCCGAT GTTTCG
	Barcode A-82	/5Phos/AGGCCAGAGCATT CGATTGGCTCGTGGCCGA TGTTTCG
	Barcode A-83	/5Phos/AGGCCAGAGCATT CGCAAGGAGCGTGGCCGA



		TGTTTCG
	Barcode A-84	/5Phos/AGGCCAGAGCATT CGCACCTTACGTGGCCGA TGTTTCG
	Barcode A-85	/5Phos/AGGCCAGAGCATT CGCCATCCTCGTGGCCGA TGTTTCG
	Barcode A-86	/5Phos/AGGCCAGAGCATT CGCCGACAACGTGGCCGA TGTTTCG
	Barcode A-87	/5Phos/AGGCCAGAGCATT CGCCTAATCCGTGGCCGA TGTTTCG
	Barcode A-88	/5Phos/AGGCCAGAGCATT CGCCTCTATCGTGGCCGAT GTTTCG
	Barcode A-89	/5Phos/AGGCCAGAGCATT CGCGACACACGTGGCCGA TGTTTCG
	Barcode A-90	/5Phos/AGGCCAGAGCATT CGCGGATTGCGTGGCCGA TGTTTCG
	Barcode A-91	/5Phos/AGGCCAGAGCATT CGCTAAGGTCGTGGCCGA TGTTTCG
	Barcode A-92	/5Phos/AGGCCAGAGCATT CGGAACAGGCGTGGCCGA TGTTTCG
	Barcode A-93	/5Phos/AGGCCAGAGCATT CGGACAGTGCGTGGCCGA



		TGTTTCG
	Barcode A-94	/5Phos/AGGCCAGAGCATT CGGAGTTAGCGTGGCCGA TGTTTCG
	Barcode A-95	/5Phos/AGGCCAGAGCATT CGGATGAATCGTGGCCGAT GTTTCG
	Barcode A-96	/5Phos/AGGCCAGAGCATT CGGCCAAGACGTGGCCGA TGTTTCG
	Barcode A-97	/5Phos/AGGCCAGAGCATT CGCGGAAGAAGTGGCCGA TGTTTCG
	Barcode A-98	/5Phos/AGGCCAGAGCATT CGGTGACAAGGTGGCCGA TGTTTCG
	Barcode A-99	/5Phos/AGGCCAGAGCATT CGGAACCAGAGTGGCCGA TGTTTCG
	Barcode A-100	/5Phos/AGGCCAGAGCATT CGTTGCTGGAGTGGCCGA TGTTTCG

Supplementary Table 6.**Barcode B Sequence**

	Barcode B-1	AAGCGTTGGCTTCTCGCA TCTAACGTGATATCCACGT GCTTGAG
	Barcode B-2	CAAGCGTTGGCTTCTCGC



		ATCTAAACATCGATCCACG TGCTTGAG
	Barcode B-3	CAAGCGTTGGCTTCTCGC ATCTATGCCTAAATCCACG TGCTTGAG
	Barcode B-4	CAAGCGTTGGCTTCTCGC ATCTAGTGGTCAATCCACG TGCTTGAG
	Barcode B-5	CAAGCGTTGGCTTCTCGC ATCTACCACTGTATCCACG TGCTTGAG
	Barcode B-6	CAAGCGTTGGCTTCTCGC ATCTACATTGGCATCCACG TGCTTGAG
	Barcode B-7	CAAGCGTTGGCTTCTCGC ATCTCAGATCTGATCCACG TGCTTGAG
	Barcode B-8	CAAGCGTTGGCTTCTCGC ATCTCATCAAGTATCCACG TGCTTGAG
	Barcode B-9	CAAGCGTTGGCTTCTCGC ATCTCGCTGATCATCCACG TGCTTGAG
	Barcode B-10	CAAGCGTTGGCTTCTCGC ATCTACAAGCTAATCCACG TGCTTGAG
	Barcode B-11	CAAGCGTTGGCTTCTCGC ATCTCTGTAGCCATCCACG TGCTTGAG



	Barcode B-12	CAAGCGTTGGCTTCTCGC ATCTAGTACAAGATCCACG TGCTTGAG
	Barcode B-13	CAAGCGTTGGCTTCTCGC ATCTAACAACCAATCCACG TGCTTGAG
	Barcode B-14	CAAGCGTTGGCTTCTCGC ATCTAACCGAGAATCCACG TGCTTGAG
	Barcode B-15	CAAGCGTTGGCTTCTCGC ATCTAACGCTTAATCCACG TGCTTGAG
	Barcode B-16	CAAGCGTTGGCTTCTCGC ATCTAAGACGGAATCCACG TGCTTGAG
	Barcode B-17	CAAGCGTTGGCTTCTCGC ATCTAAGGTACAATCCACG TGCTTGAG
	Barcode B-18	CAAGCGTTGGCTTCTCGC ATCTACACAGAAATCCACG TGCTTGAG
	Barcode B-19	CAAGCGTTGGCTTCTCGC ATCTACAGCAGAATCCACG TGCTTGAG
	Barcode B-20	CAAGCGTTGGCTTCTCGC ATCTACCTCCAAATCCACG TGCTTGAG
	Barcode B-21	CAAGCGTTGGCTTCTCGC ATCTACGCTCGAATCCACG TGCTTGAG



	Barcode B-22	CAAGCGTTGGCTTCTCGC ATCTACGTATCAATCCACG TGCTTGAG
	Barcode B-23	CAAGCGTTGGCTTCTCGC ATCTACTATGCAATCCACG TGCTTGAG
	Barcode B-24	CAAGCGTTGGCTTCTCGC ATCTAGAGTCAAATCCACG TGCTTGAG
	Barcode B-25	CAAGCGTTGGCTTCTCGC ATCTAGATCGCAATCCACG TGCTTGAG
	Barcode B-26	CAAGCGTTGGCTTCTCGC ATCTAGCAGGAAATCCACG TGCTTGAG
	Barcode B-27	CAAGCGTTGGCTTCTCGC ATCTAGTCACTAATCCACG TGCTTGAG
	Barcode B-28	CAAGCGTTGGCTTCTCGC ATCTATCCTGTAATCCACG TGCTTGAG
	Barcode B-29	CAAGCGTTGGCTTCTCGC ATCTATTGAGGAATCCACG TGCTTGAG
	Barcode B-30	CAAGCGTTGGCTTCTCGC ATCTCAACCACAATCCACG TGCTTGAG
	Barcode B-31	CAAGCGTTGGCTTCTCGC ATCTGACTAGTAATCCACG TGCTTGAG



	Barcode B-32	CAAGCGTTGGCTTCTCGC ATCTCAATGGAAATCCACG TGCTTGAG
	Barcode B-33	CAAGCGTTGGCTTCTCGC ATCTCACTTCGAATCCACG TGCTTGAG
	Barcode B-34	CAAGCGTTGGCTTCTCGC ATCTCAGCGTTAATCCACG TGCTTGAG
	Barcode B-35	CAAGCGTTGGCTTCTCGC ATCTCATACCAAATCCACG TGCTTGAG
	Barcode B-36	CAAGCGTTGGCTTCTCGC ATCTCCAGTTCAATCCACG TGCTTGAG
	Barcode B-37	CAAGCGTTGGCTTCTCGC ATCTCCGAAGTAATCCACG TGCTTGAG
	Barcode B-38	CAAGCGTTGGCTTCTCGC ATCTCCGTGAGAATCCACG TGCTTGAG
	Barcode B-39	CAAGCGTTGGCTTCTCGC ATCTCCTCCTGAATCCACG TGCTTGAG
	Barcode B-40	CAAGCGTTGGCTTCTCGC ATCTCGAACTTAATCCACG TGCTTGAG
	Barcode B-41	CAAGCGTTGGCTTCTCGC ATCTCGACTGGAATCCACG TGCTTGAG



	Barcode B-42	CAAGCGTTGGCTTCTCGC ATCTCGCATACAATCCACG TGCTTGAG
	Barcode B-43	CAAGCGTTGGCTTCTCGC ATCTCTCAATGAATCCACG TGCTTGAG
	Barcode B-44	CAAGCGTTGGCTTCTCGC ATCTCTGAGCCAATCCACG TGCTTGAG
	Barcode B-45	CAAGCGTTGGCTTCTCGC ATCTCTGGCATAATCCACG TGCTTGAG
	Barcode B-46	CAAGCGTTGGCTTCTCGC ATCTGAATCTGAATCCACG TGCTTGAG
	Barcode B-47	CAAGCGTTGGCTTCTCGC ATCTCAAGACTAATCCACG TGCTTGAG
	Barcode B-48	CAAGCGTTGGCTTCTCGC ATCTGAGCTGAAATCCACG TGCTTGAG
	Barcode B-49	CAAGCGTTGGCTTCTCGC ATCTGATAGACAATCCACG TGCTTGAG
	Barcode B-50	CAAGCGTTGGCTTCTCGC ATCTGCCACATAATCCACG TGCTTGAG
	Barcode B-51	CAAGCGTTGGCTTCTCGC ATCTGCGAGTAAATCCACG TGCTTGAG



	Barcode B-52	CAAGCGTTGGCTTCTCGC ATCTGCTAACGAATCCACG TGCTTGAG
	Barcode B-53	CAAGCGTTGGCTTCTCGC ATCTGCTCGGTAATCCACG TGCTTGAG
	Barcode B-54	CAAGCGTTGGCTTCTCGC ATCTGGAGAACAATCCACG TGCTTGAG
	Barcode B-55	CAAGCGTTGGCTTCTCGC ATCTGGTGCGAAATCCACG TGCTTGAG
	Barcode B-56	CAAGCGTTGGCTTCTCGC ATCTGTACGCAAATCCACG TGCTTGAG
	Barcode B-57	CAAGCGTTGGCTTCTCGC ATCTGTCTGTAGAATCCACG TGCTTGAG
	Barcode B-58	CAAGCGTTGGCTTCTCGC ATCTGTCTGTCAATCCACG TGCTTGAG
	Barcode B-59	CAAGCGTTGGCTTCTCGC ATCTGTGTTCTAATCCACG TGCTTGAG
	Barcode B-60	CAAGCGTTGGCTTCTCGC ATCTTAGGATGAATCCACG TGCTTGAG
	Barcode B-61	CAAGCGTTGGCTTCTCGC ATCTTATCAGCAATCCACG TGCTTGAG



	Barcode B-62	CAAGCGTTGGCTTCTCGC ATCTTCCGTCTAATCCACG TGCTTGAG
	Barcode B-63	CAAGCGTTGGCTTCTCGC ATCTTCTTCACAATCCACG TGCTTGAG
	Barcode B-64	CAAGCGTTGGCTTCTCGC ATCTTGAAGAGAATCCACG TGCTTGAG
	Barcode B-65	CAAGCGTTGGCTTCTCGC ATCTTGGAACAAATCCACG TGCTTGAG
	Barcode B-66	CAAGCGTTGGCTTCTCGC ATCTTGGCTTCAATCCACG TGCTTGAG
	Barcode B-67	CAAGCGTTGGCTTCTCGC ATCTTGGTGGTAATCCACG TGCTTGAG
	Barcode B-68	CAAGCGTTGGCTTCTCGC ATCTTTCACGCAATCCACG TGCTTGAG
	Barcode B-69	CAAGCGTTGGCTTCTCGC ATCTAACTCACCATCCACG TGCTTGAG
	Barcode B-70	CAAGCGTTGGCTTCTCGC ATCTAAGAGATCATCCACG TGCTTGAG
	Barcode B-71	CAAGCGTTGGCTTCTCGC ATCTAAGGACACATCCACG TGCTTGAG



	Barcode B-72	CAAGCGTTGGCTTCTCGC ATCTAATCCGTCATCCACG TGCTTGAG
	Barcode B-73	CAAGCGTTGGCTTCTCGC ATCTAATGTTGCATCCACG TGCTTGAG
	Barcode B-74	CAAGCGTTGGCTTCTCGC ATCTACACGACCATCCACG TGCTTGAG
	Barcode B-75	CAAGCGTTGGCTTCTCGC ATCTACAGATTCATCCACG TGCTTGAG
	Barcode B-76	CAAGCGTTGGCTTCTCGC ATCTAGATGTACATCCACG TGCTTGAG
	Barcode B-77	CAAGCGTTGGCTTCTCGC ATCTAGCACCTCATCCACG TGCTTGAG
	Barcode B-78	CAAGCGTTGGCTTCTCGC ATCTAGCCATGCATCCACG TGCTTGAG
	Barcode B-79	CAAGCGTTGGCTTCTCGC ATCTAGGCTAACATCCACG TGCTTGAG
	Barcode B-80	CAAGCGTTGGCTTCTCGC ATCTATAGCGACATCCACG TGCTTGAG
	Barcode B-81	CAAGCGTTGGCTTCTCGC ATCTATCATTCCATCCACG TGCTTGAG



	Barcode B-82	CAAGCGTTGGCTTCTCGC ATCTATTGGCTCATCCACG TGCTTGAG
	Barcode B-83	CAAGCGTTGGCTTCTCGC ATCTCAAGGAGCATCCACG TGCTTGAG
	Barcode B-84	CAAGCGTTGGCTTCTCGC ATCTCACCTTACATCCACG TGCTTGAG
	Barcode B-85	CAAGCGTTGGCTTCTCGC ATCTCCATCCTCATCCACG TGCTTGAG
	Barcode B-86	CAAGCGTTGGCTTCTCGC ATCTCCGACAACATCCACG TGCTTGAG
	Barcode B-87	CAAGCGTTGGCTTCTCGC ATCTCCTAATCCATCCACG TGCTTGAG
	Barcode B-88	CAAGCGTTGGCTTCTCGC ATCTCCTCTATCATCCACG TGCTTGAG
	Barcode B-89	CAAGCGTTGGCTTCTCGC ATCTCGACACACATCCACG TGCTTGAG
	Barcode B-90	CAAGCGTTGGCTTCTCGC ATCTCGGATTGCATCCACG TGCTTGAG
	Barcode B-91	CAAGCGTTGGCTTCTCGC ATCTCTAAGGTCATCCACG TGCTTGAG



	Barcode B-92	CAAGCGTTGGCTTCTCGC ATCTGAACAGGCATCCACG TGCTTGAG
	Barcode B-93	CAAGCGTTGGCTTCTCGC ATCTGACAGTGCATCCACG TGCTTGAG
	Barcode B-94	CAAGCGTTGGCTTCTCGC ATCTGAGTTAGCATCCACG TGCTTGAG
	Barcode B-95	CAAGCGTTGGCTTCTCGC ATCTGATGAATCATCCACG TGCTTGAG
	Barcode B-96	CAAGCGTTGGCTTCTCGC ATCTGCCAAGACATCCACG TGCTTGAG
	Barcode B-97	CAAGCGTTGGCTTCTCGC ATCTCGGAAGAAATCCACG TGCTTGAG
	Barcode B-98	CAAGCGTTGGCTTCTCGC ATCTGTGACAAGATCCACG TGCTTGAG
	Barcode B-99	CAAGCGTTGGCTTCTCGC ATCTGAACCAGAATCCACG TGCTTGAG
	Barcode B-100	CAAGCGTTGGCTTCTCGC ATCTTTGCTGGAATCCACG TGCTTGAG

Supplementary Table 7.
Chemicals and reagents



	A	B	C
	Name	Catalog number	Vender
	Formaldehyde solution	PI28906	Thermo Fisher Scientific
	HEPES pH 7.5	BBH-75-250	Boston BioProducts
	Glycine	50046	Sigma-Aldrich
	NaCl	AM9760G	Thermo Fisher Scientific
	Digitonin	G9441	Promega
	MgCl ₂	AM9530G	Thermo Fisher Scientific
	Spermidine	S0266	Sigma-Aldrich
	EDTA-free Protease Inhibitor Cocktail	11873580001	Millipore Sigma
	NP40	11332473001	Sigma-Aldrich
	EDTA Solution pH 8.0	AB00502	AmericanBio
	Bovine Serum Albumin (BSA)	A8806	Sigma-Aldrich
	Anti-H3K27me3 antibody	Ab6002	Abcam
	Anti-H3K27ac antibody	8173	Cell Signaling Technology
	Anti-H3K4me3 antibody	9751	Cell Signaling Technology



	A	B	C
	TotalSeq% ₀₀₀ -A Mouse Universal Cocktail	199901	Biolegend
	Anti mouse CD4	A0001	Biolegend
	Anti mouse CD3	A0182	Biolegend
	Anti mouse CD34	A0857	Biolegend
	Anti mouse CD140a	A0573	Biolegend
	Anti mouse CD133	A1037	Biolegend
	Anti mouse CD90.1	A0380	Biolegend
	Anti mouse CD90.2	A0075	Biolegend
	Anti mouse B220	A0103	Biolegend
	TruStain FcX% ₀₀₀ PLUS (anti-mouse CD16/32) antibody	156604	Biolegend
	Cell Staining Buffer	420201	Biolegend
	Fab Fragment Goat Anti-Mouse IgG	115-007-003	Jackson ImmunoResearch
	Triton X-100	T8787	Sigma-Aldrich
	T4 DNA Ligase	M0202L	New England Biolabs
	T4 DNA Ligase Reaction Buffer	B0202S	New England Biolabs
	NEBuffer 3.1	B7203S	New England Biolabs
	DPBS	14190144	Thermo Fisher Scientific



	A	B	C
	Proteinase K	E00491	Thermo Fisher Scientific
	SPRI beads	A63880	Beckman Coulter
	NEBNext High-Fidelity 2X PCR Master Mix	M0541L	New England Biolabs
	SYBR Green I Nucleic Acid Gel Stain	S7563	Thermo Fisher Scientific
	DNA Clean & Concentrator-5	D4014	Zymo Research
	Tn5 Transposase - unloaded	C01070010	Diagenode
	Tagmentation Buffer (2x)	C01019043	Diagenode
	Sodium dodecyl sulfate	71736	Sigma-Aldrich
	Maxima H Minus Reverse Transcriptase (200 U/L)	EP0751	Thermo Fisher Scientific
	dNTP mix	R0192	Thermo Fisher Scientific
	Dynabeads [®] MyOne [®] Streptavidin C1	65001	Thermo Fisher Scientific
	RNase Inhibitor	Y9240L	Enzymatics
	Kapa Hotstart HiFi ReadyMix	KK2601	Kapa Biosystems



	A	B	C
	Nextera XT DNA Library Preparation Kit	FC-131-1024	Illumina

Safety warnings

- ⚠ Digitonin is toxic and care should be taken especially when weighing out the powder. Use full PPE including a mask, lab coat and gloves while handling any amount of digitonin.



Cryosectioning Procedure

- 1 Specimens should be handled with RNase Away.

20m

Note

This protocol is to be used with tissue preserved in OCT and stored in -80 °C .

- Cool Cryostat to -20 °C . Clean the work surfaces with RNase away and change a new cutting blade.
- Place a clean slide box inside the cryostat chamber to store slides.
- Clean the tissue holder with RNase away and put it in the cryostat chamber . Adhere the OCT specimen to a tissue holder and allow it to equilibrate for 00:20:00 to reach the chamber temperature and strengthen the adhesion between the OCT block and the holder.
- Place the specimen in the cryostat with the tissue holder and cut at Thickness 8-10 µm.
- Store completed slides at -80 °C .

Tn5 loading

10m

- 2 **Annealing adaptor sequences:**

1h 35m

Tn5MErev:

5'-/Phos/CTGTCTCTTATACACATCT-3'

Tn5ME-A:

5'-/Phos/CATCGGCGTACGACTAGATGTGTATAAGAGACAG-3'

Tn5ME-B1:

5'-/Phos/CATCGGCGTACGACT**TAGCCT**AGATGTGTATAAGAGACAG-3'

Tn5ME-B2:

5'-/Phos/CATCGGCGTACGACT**ATAGAG**AGATGTGTATAAGAGACAG-3'

Tn5ME-B3:

5'-/Phos/CATCGGCGTACGACT**CCTATC**AGATGTGTATAAGAGACAG-3'

- Resuspend the oligonucleotides (Tn5MErev, Tn5ME-A, Tn5ME-B) in water to a final concentration of 100 µM each.

- Mix equimolar amounts of Tn5MErev/Tn5ME-A and Tn5MErev/Tn5ME-B in separate 200 µl PCR tubes.

e.g:

tube1: 10 µL Tn5ME-A + 10 µL Tn5MErev
 tube2: 10 µL Tn5ME-B + 10 µL Tn5MErev

- Denature in the thermocycler for 00:05:00 at 95 °C , and cool down slowly on the thermocycler by ramping down to 12 °C by 0.1 °C per second

Pause point: Store the annealed oligos at -20 °C

- Mix annealing mix.

The final Tn5 concentration is 2 µM of Tn5 dimer.

Use unique barcodes for specific nano-Tn5 or WT-Tn5.

2x Tn5 loading buffer

2x Tn5 loading buffer		Can be stored at 4C without DTT *	
final conc	stock	stock	for 10ml
100 mM	HEPES pH 7.2	250 mM	4 ml
200 mM	NaCl	5M	400 µl
0.2 mM	EDTA	0.5M	4µl
0.2 %	Triton-X	10%	200 µl
20%	Glycerol	100%	2ml
2 mM	DTT	200 mM	*
	water		3.4 ml

Mouse Nano-Tn5 (MeA/B1/R loaded):

4 µL Annealed, barcoded MeA/Me-Rev oligos Concentration 50 micromolar (µM)
 4 µL Annealed, barcoded MeB1/Me-Rev oligos Concentration 50 micromolar (µM)
 42 µL Glycerol
 6 µL Nano-Tn5 (mouse, 5mg/ml, MW = 73941 g/mol; Concentration 67.6 micromolar (µM))
 44 µL 2x Nano-Tn5 loading buffer

**Rabbit Nano-Tn5 (MeA/B2/R loaded):**

-  4 μL Annealed, barcoded MeA/Me-Rev oligos Concentration50 micromolar (μM)
-  4 μL Annealed, barcoded MeB2/Me-Rev oligos Concentration50 micromolar (μM)
-  42 μL Glycerol
-  4.4 μL Nano-Tn5 (rabbit, 6.8mg/ml, MW = 73013 g/mol Concentration93.1 micromolar (μM))
-  45.6 μL 2x Nano-Tn5 loading buffer

For nano-Tn5, incubate for  01:00:00 at  23 °C .

WT Tn5-MeA with barcode (MeA/B3/R for ATAC-seq)

-  5 μL Annealed, barcoded MeA/Me-Rev oligos Concentration50 micromolar (μM)
-  5 μL Annealed, barcoded MeB3/Me-Rev oligos Concentration50 micromolar (μM)
-  10 μL Tn5 (Order from Diagenode)

For WT-Tn5, incubate for  00:30:00 at  23 °C .

Then mix with  10 μL Glycerol.

Note

Store the loaded nano-Tn5 or WT-Tn5 at Temperature-20 °C.

Fixation

- 3 Frozen tissue slides were first thawed for  00:01:00 at  37 °C . 1m
- 4 Tissue was fixed with  1 mL [M] 0.2 % volume formaldehyde in PBS (with 0.05 U μL^{-1} RNase Inhibitor) at  Room temperature for  00:05:00 . 5m

- 5

Discard reagent and quench with  1 mL [M] 1.25 Molarity (M) glycine (with 0.05 U μl^{-1} RNase Inhibitor) for  00:05:00 .

5m
- 6

Discard reagent and wash slide with 1x PBS (with 0.05 U μl^{-1} RNase Inhibitor) for  00:01:00 .

1m
- 7

Discard reagent and rinse slide sequentially in RNase free water. Wipe away excess liquid.
- 8

Allow slide to air dry and image the slide on a Keyence BZ-X810 scanning microscope at 10x or 20x resolution.

ATAC-seq

10m

9 Buffer preparation for ATAC-seq

Wash Buffer I				
Reagent	1 ml	2 ml	Stock Concentration	Final Concentration
Tris-HCl (pH 7.4)	10 μl	20 μl	1 M	10 mM
NaCl	2 μl	4 μl	5 M	10 mM
MgCl ₂	3 μl	6 μl	1 M	3 mM
BSA	200 μl	400 μl	5%	1%
Tween-20	10 μl	20 μl	10%	0.1%
Nuclease-free Water	775 μl	2 x 775 μl	-	-



1X Lysis Buffer				
Reagent	1 ml	2 ml	Stock Concentration	Final Concentration
Tris-HCl (pH 7.4)	10 μl	20 μl	1 M	10 mM
NaCl	2 μl	4 μl	5 M	10 mM
MgCl ₂	3 μl	6 μl	1 M	3 mM
Tween-20	10 μl	20 μl	10%	0.1%
NP-40	10 μl	20 μl	10%	0.1%
Digitonin	10 μl	20 μl	1%	0.01%
BSA	200 μl	400 μl	5%	1%
Nuclease-free Water	755 μl	2 x 755 μl	-	-

Lysis Dilution Buffer				
Reagent	1ml	2 ml	Stock Concentration	Final Concentration
Tris-HCl (pH 7.4)	10 μl	20 μl	1 M	10 mM
NaCl	2 μl	4 μl	5 M	10 mM
MgCl ₂	3 μl	6 μl	1 M	3 mM
BSA	200 μl	400 μl	5%	1%
Nuclease-free Water	785 μl	2 x 785 μl	-	-

0.1X Lysis Buffer			
Reagent	2 ml	Stock Concentration	Final Concentration
1X Lysis Buffer	200 μl	1X	0.1X
Lysis Dilution Buffer	2 ml – 200 μl	-	-

- 10

Attach the PDMS reservoir and 12mm diameter clamp to the tissue.



Note

The PDMS reservoir should cover the entire tissue region of interest on the slide, but small enough to save reagents.

- 11 Permeabilize the tissue with  200 μL 0.1X lysis buffer and incubate at

15m

 Room temperature for  00:15:00 .

- 12 Discard reagent and Wash the tissue with  200 μL wash buffer I at

5m

 Room temperature for  00:05:00 .

- 13 Discard reagent and quickly wash the tissue with  200 μL wash buffer.

- 14 Prepare transposition mix:

 50 μL 2X TD Buffer (see materials section)

 33 μL 1X PBS

 1 μL 10% Tween-20 (final 0.1% v/v)

-  1 μL

1% Digitonin (final 0.01% v/v)
-  5 μL

Tn5 Transposase (WT Tn5)
-  10 μL

nuclease-free H_2O

15

Discard reagent and add

 100 μL

 transposition mix.

16

Incubate for

 00:30:00

 at

 37 $^{\circ}\text{C}$

 .

30m

17

Discard reagent and stop the tagmentation by adding

 200 μL

 of

[M] 40 millimolar (mM) EDTA at

 Room temperature

 for

 00:05:00

 .

5m

Nano-CUT&Tag

18

Buffer preparation for Nano-CUT&Tag

2 M spermidine:

Reagent

Final concentration

Final 2 M volume (50 μL)

ddH₂O

N/A

34.3 μL

spermidine

15.7 μL

Wash bufferII:

mix 1 ml 1 M HEPES pH 7.5, 1.5 ml 5 M NaCl, 12.5 μL 2 M spermidine, bring the final volume to 50 ml with ddH₂O, and add 1 Roche Complete Protease Inhibitor EDTA-Free tablet.

Alternative: dissolve 1 tablet inhibitor into 2 ml ddH₂O.

Store @ 4°C for up to 1 week.

reagent	Final concentration	Volume for 50 ml	Volume for 10 ml	Volume for 5 ml	Volume for 2 ml	Volume for 1 ml
1 M HEPES pH 7.5	20 mM	1 mL	0.2 ml	100 μL	40 μL	20 μL
5 M NaCl	150 mM	1.5 mL	0.3 ml	150 μL	60 μL	30 μL
2 M spermidine	0.5 mM	12.5 μL	2.5 μL	1.25 μL	0.5 μL	0.25 μL
ddH ₂ O		45.4875 ml	9.0974 ml	4.54875 ml	1.81948 ml	910 μL
Protease inhibitor		2 ml	0.4 ml	200 μL	80 μL	40 μL

NP-40 Digitonin wash buffer (permeabilize):

0.01% NP40, 0.01% Digitonin in wash buffer.

prepare fresh each day and store @4°C.

Reagent	Final concentration	Volume for 2 ml	Volume for 1.5 ml	Volume for 1 ml	Volume for 0.5 ml
10% NP40	0.01%	2 μL	1.5 μL	1 μL	0.5 μL
1% Digitonin	0.01%	20 μL	15 μL	10 μL	5 μL
Wash buffer		989 x 2 μL	741.75 x2 μL	989 μL	494.5 μL

Digitonin is supplied at 2% in DMSO. Dilute 1:1 with ddH₂O to make a 1% stock solution. Avoid more than five freeze thaw cycles. Can be kept at -20°C up to 6 Months.

Antibody buffer:
antibody buffer mixes 8 μL 0.5 M EDTA and 6.7 μL 30% BSA with 2 ml NP40-Digitonin wash buffer and chill on ice. (Mix 4 μL 0.5 M EDTA and 10 μL 1% BSA with 1 ml NP40-DIGITONIN wash buffer and chill on ice) - prepare fresh each day and store @4°C.

Reagent	Final concentration	Volume for 1 ml	Volume for 0.5 ml	Volume for 0.25 ml
0.5 M EDTA	2 mM	4 μL	2 μL	1 μL
30% BSA	0.001%	3.35 μL	1.675 μL	0.8375 μL
NP40-Digitonin wash buffer	N/A	992.65 μL	496.325 μL	248.1625 μL

300 Wash buffer:
mix 1 ml 1 M HEPES pH 7.5, 3 ml 5 M NaCl, 12.5 μL 2 M spermidine, bring the final volume to 50 ml with dH₂O, and add 1 Roche Complete Protease Inhibitor EDTA-Free tablet.
Alternative: dissolve 1 tablet inhibitor into 2 ml dH₂O.
~Store @ 4°C for up to 1 week.

reagent	Final concentration	Volume for 50 ml	Volume for 10 ml	Volume for 2 ml	Volume for 1 ml
1 M HEPES pH 7.5	20 mM	1 mL	0.2 ml	40 μL	20 μL
5 M NaCl	300 mM	3 mL	0.6 ml	120 μL	60 μL
2 M spermidine	0.5 mM	12.5 μL	2.5 μL	0.5 μL	0.25 μL
dH ₂ O		43.9875 ml	8.7974 ml	880 x 2 μL	880 μL
Protease inhibitor		2 ml	0.4 ml	80 μL	40 μL

(segmentation solution:
Mix 1ml 300-wash buffer and 10 μL 1 M MgCl₂ (final 10 mM)
~Store @4°C for up to 1 week

Reagent	Final concentration	Volume for 14 ml	Volume for 12 ml	Volume for 6 ml	Volume for 1 ml
1 M MgCl ₂	10 mM	140 μL	120 μL	60 μL	10 μL
300-wash buffer		13.86 ml	11.88 μL	5.94 ml	990 μL

- 19

Discard reagent and Wash the tissue with

200 μL

 wash buffer II at

Room temperature

 for

00:05:00

 .

5m

20

Discard reagent and Wash the tissue with

200 μL

 NP40-Digitonin wash buffer at

Room temperature

 for

00:05:00

 .

5m

Primary antibody binding

21

Starting concentrations (can be further optimised, depending on the antibody)

1:50 primary antibody in antibody-binding buffer.

Final volume

100 μL

 per sample.

22

Discard reagent and add

100 μL

 primary antibody mix.

23

Incubate at

4 °C

Overnight

20m

Note

Make sure the liquid can cover the tissue region of interest.

protocols.io | <https://dx.doi.org/10.17504/protocols.io.q26g71n79gwz/v1> October 22, 2024 36/49



Nano-Tn5 binding

- 24 Discard reagent and Wash the tissue with  200 μL 300-wash buffer at  Room temperature for  00:05:00 . 5m
- 25 Starting concentrations of nano-Tn5s mix
1:100 nano-Tn5s in 300-wash buffer.
Final volume  100 μL per sample.
- 26 Discard reagent and add  100 μL nano-Tn5s mix at  Room temperature for  01:00:00 . 1h
- 27 Discard reagent and Wash the tissue with  200 μL 300-wash buffer at  Room temperature for  00:05:00 . 5m
- 28 Discard reagent and add  200 μL Tagmentation buffer.
- 29 Incubate for  01:00:00 at  37 °C . 1h
- 30 Discard reagent and stop the tagmentation by adding  200 μL of  40 millimolar (mM) EDTA at  Room temperature for  00:05:00 . 5m

Surface Protein Antibody staining

- 31 Discard reagent and add  200 μL of cell surface staining buffer at  Room temperature for  00:02:00 . 2m
- 32 Discard reagent and add  200 μL of 1:20 TruStain FcX™ in Cell Staining Buffer at  4 °C for  00:15:00 . 15m
- 33 Discard reagent and add  200 μL of cell surface staining buffer at  Room temperature for  00:01:00 . 1m
- 34 Discard reagent and add  200 μL of oligonucleotide-labeled antibody-derived tags (1:100) in cell surface staining buffer.



- 35 Incubate for 00:30:00 at 4 °C . 30m
- 36 Discard reagent and add 200 μL of cell surface staining buffer at Room temperature for 00:01:00 . 1m
- 37 Discard reagent and add 200 μL of Fab Fragment Goat Anti-Mouse IgG (1:25) in cell surface staining buffer.
- 38 Incubate for 00:15:00 at 4 °C . 15m

In situ reverse transcription

- 39 Discard reagent and tissue was fixed with 200 μL 2 % volume formaldehyde in PBS (with 0.05 U μL^{-1} RNase Inhibitor) at Room temperature for 00:10:00 . 10m
- 40 Discard reagent and quench with 1 mL 1.25 Molarity (M) glycine (with 0.05 U μL^{-1} RNase Inhibitor) for 00:05:00 . 5m
- 41 Discard reagent and add 200 μL of 0.5x PBS at Room temperature for 00:05:00 . 5m
- 42 Prepare the RT mixture

For one sample:

12.5 μL	5X RT Buffer
0.8 μL	RNase inhibitor
3.1 μL	dNTPs
6.1 μL	Maxima H Minus Reverse Transcriptase
20.3 μL	0.5X PBS
10 μL	RT primer
9.6 μL	ddH ₂ O



- 43 Discard reagent and add 62.5 μL RT mixture.
- 44 Incubate at Room temperature 00:30:00 . 30m
- 45 Incubate at 42 $^{\circ}\text{C}$ 01:30:00 . 1h 30m
- 46 Discard reagent and add 200 μL of 1X NEBBuffer 3.1 at Room temperature for 00:05:00 . 5m
- 47 Discard reagent, remove the clamp, and dip the tissue slide into ddH₂O.
- 48 Air dry the slide.

In-tissue ligation with barcode A

- 49 **Prepare buffer:** 5m
2X annealing buffer:
 2 millimolar (mM) EDTA, pH 8.0
 20 millimolar (mM) Tris-HCl, pH 7.5
 100 millimolar (mM) NaCl
- Prepare annealed barcode An**
 25 μL 2X annealing buffer
 12.5 μL 100 micromolar (μM) barcode An
 12.5 μL 100 micromolar (μM) ligation linker 1
- Prepare annealed barcode Bn**
 25 μL 2X annealing buffer
 12.5 μL 100 micromolar (μM) barcode Bn
 12.5 μL 100 micromolar (μM) ligation linker 2

Run the following program:



🌡️ 95 °C ⌚ 00:05:00 , and cool to 🌡️ 12 °C , -0.1 °C/sec.

Note

Annealed barcodes can be stored at 🌡️ -20 °C

Prepare spatial barcode mix, for each channel

🧴 27 µL 10X T4 DNA ligase buffer

🧴 11 µL T4 DNA ligase

🧴 5.4 µL 5% Triton X-100

🧴 11.58 µL 10X NEB buffer3.1

🧴 175 µL Nuclease free water

Prepare spatial barcode mix, for each channel:

🧴 4 µL ligation mix

🧴 1 µL annealed barcode An (such as A1)

- 50 Apply the microfluidic chip A to the tissue slide to cover the region of interest.

Note

Ensure the microfluidic chip is dust free by applying tape on the microfluid chip and tissue slide to gently remove dust.

- 51 Image full chip A on the tissue slide with no clamp using Keyence microscopy with an 10× objective lens.
- 52 Load the slide into a clamp. Secure the clamp onto the slide, and tighten the screws.
- 53 Load the spatial barcode mix A to the inlets, and attach a vacuum gasket over the outlets.



54 Flow the chip until the solution has reached the outlet wells for all channels.

55 Remove the vacuum gasket and incubate the tissue slide at  37 °C for  00:30:00 .

30m

56 Suck the solution out of the inlet and outlet.

57 Dip the slide in RNase free water and air dry the slide.

In-tissue ligation with barcode B

10m

58 Follow the steps from 50 to step 57, and only need to replace the barcode A with barcode B..

59 Image the tissue slide using Keyence microscopy with an 10× objective lens.

Lysis

60 Prepare the reverse crosslinking buffer (for one tissue slide)

 5 µL 1M Tris-HCl,pH 8.0

 10 µL 10mM EDTA,pH 8.0

 10 µL 10% SDS

 4 µL 5M NaCl

 2 µL 20mg/ml proteinase K

 69 µL Nuclese-free water

61 Attach the PDMS reservoir and 12mm diameter clamp to the tissue.

62 Add  100 μL reverse crosslinking buffer and cover the hole with Parafilm or tapes.

63 Incubate the tissue slide at  60 °C for  02:00:00 .

2h

64 Collect the lysate into a fresh PCR tube and Incubate the tissue slide at  60 °C for  Overnight .

Pause point Store lysate at -80 °C.

65 Purify the lysate with Zymo DNA Clean and concentrator kit in a 5:1 ratio (sample : buffer).

66 Elute the purified DNA with  100 μL RNase free water.

gDNA and cDNA separation

10m

67 Prepare buffer:

2x Wash Buffer (2xW&B)

reagents	Final concentration	volume	volume	volume
1M Tris-HCl, pH8.0	10mM	2.5ul	5ul	10ul
5M NaCl	2M	100ul	200ul	400ul
0.5M EDTA, pH8.0	1mM	0.5ul	1ul	2ul
Nuclear-free water	2x588ul	147ul	294ml	588ul
Total		0.25ml	0.5ml	1ml
SUPERase IN inhibitor 20U/ul	0.375U/ul	6.25ul	12.5ul	
		1x reaction	2x reactions	5x reactions

1x Binding Wash buffer-Tween (1x TBW)				
reagents	Final concentration	volume	volume	volume
1M Tris-HCl, pH8.0	5mM	5ul	25ul	50ul
5M NaCl	1M	200ul	1ml	2ml
0.5M EDTA, pH8.0	0.5mM	1ul	5ul	10ul
10% Tween 20	0.05%	5ul	25ul	50ul
Nuclear-free water		789ul	3.945ml	7.89ml
Total		1ml	5ml	10ml
SUPERase IN inhibitor 20U/ul	0.0275U/ul	1.38ul	6.875ul	

1x Tris-HCl Buffer (1xTris-HCl Buffer)				
reagents	Final concentration	volume	volume	volume
1M Tris-HCl, pH8.0	10mM	5ul	10ul	20ul
10% Tween 20	0.1%	5ul	10ul	20ul
Nuclear-free water		490ul	980ul	2x980ul
Total		500ul	1ml	2ml
SUPERase IN inhibitor 20U/ul	0.05U/ul	1.25ul	2.5ul	5ul

- 68 Wash  40 µL Dynabeads™ MyOne™ Streptavidin C1 beads with  800 µL 1XTBW buffer.
- 69 Repeat step 68 for 2 times.
- 70 Resuspend the beads with  100 µL 2xBW buffer and the purified DNA from step 66.
- 71 Agitate the DNA-beads at  Room temperature for  01:00:00 .
- 72 Use a magnetic rack to separate beads (containing cDNA/ADT) and supernatant (containing gDNA).

1h

RNA/ADT library construction



- 73 Add  400 μL of 1X TBW buffer to the seperated beads and resuspend.
- 74 Separate the beads from the supernatant on the magnetic rack for 1 min and remove the supernatant.
- 75 Repeate step 73-74.
- 76 Add  400 μL of 1XTris-HCl buffer to the seperated beads and resuspend.
- 77 Agitate beads for  00:05:00 at  Room temperature . 5m
- 78 Separate the beads from the supernatant on the magnetic rack for 1 min and remove the supernatant.
- 79 Add  400 μL of RNase free water to the seperated beads.
- 80 Prepare TSO solution
- | | |
|---|---|
|  22 μL | 10mM dNTPs |
|  44 μL | 5X Maxima RT buffer |
|  44 μL | 20% Ficoll PM-400 solution |
|  5.5 μL |  100 micromolar (μM) TSO primer |
|  11 μL | Maxima H Minus Reverse Transcriptase |
|  5.5 μL | RNase Inhibitor |
|  88 μL | RNase free water |
- 81 Separate the beads from the supernatant on the magnetic rack for 1 min and remove the supernatant.
- 82 Add  200 μL TSO solution and gently shake the beads at  Room temperature for  00:30:00 30m



- 83 Gently shake the beads at 42 °C for 01:30:00 1h 30m
- 84 Separate the beads from the supernatant on the magnetic rack for 1 min and remove the supernatant.
- 85 Add 400 µL of 1XTris-HCl buffer to the seperated beads and resuspend.
- 86 Agitate beads for 00:05:00 at Room temperature . 5m
- 87 Separate the beads from the supernatant on the magnetic rack for 1 min and remove the supernatant.
- 88 Add 400 µL of RNase free water to the seperated beads.
- 89 Prepare the PCR mixture
- | | |
|---------|----------------------------------|
| 110 µL | 2X KAPA HIFI HotStart Master Mix |
| 8.8 µL | 10 micromolar (µM) Primer 1 |
| 8.8 µL | 10 micromolar (µM) Primer 2 |
| 1 µL | 10 micromolar (µM) Primer 3 |
| 91.4 µL | Nuclease Free water |
- 90 Add 200 µL PCR mix to the separated beads. Splite the beads into 4 PCR tubes.
- 91 PCR thermocyling was performed as following:



Cycling conditions	
95°C	3 min
Then 5 cycles of:	
98°C	20 sec
65°C	45 sec
72°C	3 min
Hold at 4 °C	

92 Separate the beads from the supernatant on the magnetic rack for 1 min and transfer the supernatant to a fresh tube (1.5 ml).

93 Add  10 µL 20X EvaGreen and split into 5 qPCR tubes ( 46 µL for each tube).

94 Run with a qPCR machine with the following thermocycling conditions:

Cycling conditions	
95°C	3 min
Then 20 cycles of:	
98°C	20 sec
65°C	20 sec
72°C	3 min
72°C	5 min
Hold at 4 °C	

Once the qPCR signal began to plateau, reactions are stopped.

95 PCR reactions are purified using a 0.6x ratio of SPRI beads according to the manufacturer's protocol, and the cDNA is eluted with  20 µL Nuclease free water.

Note

Keep the supernatant for ADT library construction.

96 Confirm the quality and quantity of the cDNA library with the agilent tapestation.

**Note**

The sample can be stored at -20 °C

97 RNA library construction

Nextera XT DNA Library Preparation Kit is used for final RNA library construction.

98 ADT library construction

During SPRI cleanup (Step 95), the supernatant is saved and the short DNA derived from antibody oligos is purified with 2.0× SPRI beads:

- Transfer the supernatant to a new tube and add 1.4X SPRIselect reagent to obtain a final SPRI volume of 2X SPRI.
- purify the SPRI beads according to the manufacturer's protocol, and the beads are eluted with 50 µL Nuclease free water.
- Second time purification with 2.0X SPRI beads and the beads are eluted with 30 µL Nuclease free water.

ADT library amplification:

Prepare 100µl PCR reaction with purified ADT:

- 30 µL purified ADT fraction
- 50 µL 2x KAPA Hifi PCR Master Mix.
- 1.5 µL 10 micromolar (µM) Nextera DNA Index.
- 1.5 µL 10 micromolar (µM) P5 oligo.
- 17 µL Nuclease free water.

PCR program:

95°C 3 min

95°C 20 sec |

60°C 30 sec | ~ 6-10 cycles

72°C 20 sec |

72°C 5 min



PCR reactions are purified using a 1.6x ratio of SPRI beads according to the manufacturer's protocol, and the cDNA is eluted with  20 μ L Nuclease free water.

Confirm the quality and quantity of the cDNA library with the agilent tapestation.

The sample can be stored at  -20 °C

ATAC/nano-CUT&Tag library construction

- 99 During Step 72, the supernatant is saved and the gDNA is purified with Zymo DNA Clean and concentrator kit in a 5:1 ratio (sample : buffer).

The column is eluted with  45 μ L Nuclease free water.

ATAC/nano-CUT&Tag **library amplification:**

Prepare 100ul PCR reaction with purified ADT:

-  45 μ L purified ATAC/nano-CUT&Tag fraction
-  50 μ L 2X NEB Master Mix.
-  2.5 μ L  10 micromolar (μ M) Nextera DNA Index.
-  2.5 μ L  10 micromolar (μ M) P5 oligo.

PCR program:

58°C 5 min

72°C 5 min

98°C 30 sec

98°C 10 sec |

60°C 30 sec | ~ 13 cycles

72°C 5 min

PCR reactions are purified using a 1.2x ratio of SPRI beads according to the manufacturer's protocol, and the final gDNA library is eluted with  20 μ L Nuclease free water.

Confirm the quality and quantity of the library with the agilent tapestation.

The sample can be stored at  -20 °C

