

Jan 27, 2023

Spatial Transcriptomics for OCT using 10x Genomics Visium

 [Nature](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9d3opg3e/v1>

Andrew Houston¹, Siqi Chen¹, Feng Chen¹

¹Department of Medicine, Washington University in St Louis



Andrew Houston

Washington University in St. Louis

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9d3opg3e/v1>

External link: <https://doi.org/10.1038/s41586-024-08087-4>

Protocol Citation: Andrew Houston, Siqi Chen, Feng Chen 2023. Spatial Transcriptomics for OCT using 10x Genomics Visium. protocols.io <https://dx.doi.org/10.17504/protocols.io.x54v9d3opg3e/v1>

Manuscript citation:

Mo, CK., Liu, J., Chen, S. *et al.* Tumour evolution and microenvironment interactions in 2D and 3D space. *Nature* **634**, 1178–1186 (2024). <https://doi.org/10.1038/s41586-024-08087-4>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 25, 2023

Last Modified: January 27, 2023

Protocol Integer ID: 75822

Keywords: 10x genomics visium protocol, spatial transcriptomic, using 10x genomics visium, 10x genomics visium, 10x genomic, gene expression, frozen tissue sample, gene

Abstract

To detect gene expression spatially mapped across a fresh-frozen tissue sample adapted from 10x Genomics Visium protocol.



Visium Fresh-Frozen

- 1 Start [here](#)
If performing OCT tissue embedding for the first time be sure to consult pgs.5-9
If utilizing a new tissue type for the first time be sure to utilize the Tissue Optimization Guide which can be found [here](#)
- 2 RNA extraction and assessment is outlined in the above guide but not always necessary
Most OCT blocks have preserved RNA quality well due to their preservation/storage
- 3 Perform all fixation/H&E staining/imaging steps on the Visium Spatial Gene Expression Slide (PN-2000233)
Utilize the guide [here](#)
Typically it is not necessary to utilize a coverslip while imaging
- 4 Immediately proceed to the Visium Library Construction
Start on pg. 35 using this guide [here](#)
- 5 Sequence the libraries according to the requirements on pg. 59 of CG000239