

Feb 07, 2025

Specific (first) PCR of the Leeselab 2-step PCR protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1dmy7vr2/v1>



Marie-Thérèse Werner¹, Dominik Buchner¹

¹University of Duisburg-Essen, Aquatic Ecosystem Research

Aquatic Ecosystem Rese...



Dominik Buchner

University of Duisburg-Essen, Aquatic Ecosystem Research

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1dmy7vr2/v1>

Protocol Citation: Marie-Thérèse Werner, Dominik Buchner 2025. Specific (first) PCR of the Leeselab 2-step PCR protocol. protocols.io <https://dx.doi.org/10.17504/protocols.io.ewov1dmy7vr2/v1>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 29, 2025

Last Modified: February 07, 2025

Protocol Integer ID: 119269

Keywords: pcr, first PCR step, two-step PCR, pcr of the leeselab, first pcr step, product in the second pcr, step pcr protocol, step pcr protocol this protocol, specific first step of the leeselab, specific primer, step primer, second pcr, primer, pcr, universal primer, addition to the specific primer, well pcr plate, leeselab, specific first step, fwhr2n, target sequence, second step, bp tag

Abstract

This protocol describes how to perform the specific first step of the Leeselab 2-step PCR protocol in  10 µL reactions. In the first PCR step, universal primers (e.g. fwhF2/fwhR2n) are used to amplify the target sequence. There are up to 24 different first-step primers individually labeled with a 4 - 6 bp tag (see table in the material section). These are used to individually tag up to 24 different 96-Well PCR Plates. In addition to the specific primer and the tag, each of the primers contains a universal overhang, which is used to further amplify the product in the second PCR. In the second step of the 2-step-PCR, each well is tagged with an individual Index (1-96) enabling the identification of each sample (well) after library sequencing.

Guidelines

Follow general lab etiquette. Wear gloves to prevent contaminating the samples. Clean the workspace before starting with 80% EtOH.

Materials

Below all materials needed for the protocol are listed. Vendors and part numbers are listed but interchangeable.

Chemicals:

Multiplex Mastermix  5x Hot-Start Multiplex PCR-Mastermix **Bio&SELL Catalog #BS91.522.1000MPP**

Forward Primer (universal), see Excel spreadsheet below

Reverse Primer (universal), see Excel spreadsheet below

PCR-grade water

 Invitrogen UltraPure DNase/RNase-Free Distilled Water **Fisher Scientific Catalog #11538646**

Labware:

Micro-Twist-Top-Tube  Micro-Twist-Top-Tube **Sarstedt Catalog #72.693**

PCR plate Hard-Shell 96-well  Hard-Shell 96-well PCR plate **Bio-Rad Laboratories Catalog #HSP9601**

PCR foil, DNase-/RNase-free  PCR foil, DNase-/RNase-free **Sarstedt Catalog #95.1993**

Devices:

Vortexer

Centrifuge for SBS format microplates

Table centrifuge for tubes

Pipettes and Multichannel pipettes including pipette tips

Thermocycler compatible with low-profile full-skirted PCR plates

Documents:

 First_step_primer_Leeselab.xlsx 9.7KB

Protocol materials

 Hard-Shell 96-well PCR plate **Bio-Rad Laboratories Catalog #HSP9601**

 Micro-Twist-Top-Tube **Sarstedt Catalog #72.693**

 5x Hot-Start Multiplex PCR-Mastermix **Bio&SELL Catalog #BS91.522.1000MPP**

 UltraPure™ DNase/RNase-Free Distilled Water **Thermo Fisher Scientific Catalog #10977023**

 PCR foil, DNase-/RNase-free **Sarstedt Catalog #95.1993**

 Invitrogen UltraPure DNase/RNase-Free Distilled Water **Fisher Scientific Catalog #11538646**

Safety warnings

 Reagents are potentially damaging to the environment. Dispose waste responsibly.

Before start

Make sure all materials are available before starting. Make sure all materials and samples are unfrozen.

Note

The protocol described here is designed for the use of PCR plates, but can also be done in tubes, strips, or any sufficient reaction vessel. The recommended shaking speeds are adjusted to the plates mentioned in the materials.

Preparations

26m

- 1 Calculate the required volumes for the PCR Mastermix using the Excel spreadsheet.

1m

Note

You can download the Excel spreadsheet here:

 Mastermix_calculation.xlsx



- 2 Prepare the workplace to ensure contamination-free handling.


2m



Workspace prepared for a right-handed person to work from left (pipette tips, samples) to right (waste). Samples and (unused) labware are placed in the back.

- 2.1 Label the  Micro-Twist-Top-Tube **Sarstedt Catalog #72.693** (e.g., "MM") and  Hard-Shell 96-well PCR plate **Bio-Rad Laboratories Catalog #HSP9601** (e.g. project name, date, "1st_pcr", "A").

2m



Multiplex Mastermix (M), forward (fwd) and reverse (rev) primer, PCR-grade water and PCR Mastermix (MM).

- 3 Briefly vortex and shortly spin down all samples ( DNA extract) and required chemicals.

2m




- 4 Pipette the PCR Mastermix in a  Micro-Twist-Top-Tube **Sarstedt Catalog #72.693** containing the calculated volumes of

3m

 5x Hot-Start Multiplex PCR-Mastermix **Bio&SELL Catalog #BS91.522.1000MPP** , forward primer, reverse primer, and



 UltraPure™ DNase/RNase-Free Distilled Water **Thermo Fisher Scientific Catalog #10977023**



Preparation of the PCR Mastermix (MM).

Note

This step has to be performed separately per plate when using the tagging scheme mentioned above since per plate an individual combination of primers is used per plate.

- 4.1 Briefly vortex and shortly spin down the prepared PCR Mastermix.

1m



- 5 Distribute  9 μL per sample in each well used of the

8m

 Hard-Shell 96-well PCR plate **Bio-Rad Laboratories Catalog #HSP9601** . To speed up this step, a Multidispenser can be used.






Distribution of the PCR Mastermix (MM) in all required wells.


- 6 Pipette  1 μL  DNA extract per well. This step has to be done with a Multichannel pipette to ensure correct assignment and minimize the contamination risk.

4m






Always maintain the order of samples (here top to bottom, left to right) for pipette tips, source plate and target plate (see below).



Use 1 μ L of each sample (DNA extract).




Fill the DNA extracts in the respective wells (same order on each plate).

- 7 Close the  Hard-Shell 96-well PCR plate **Bio-Rad Laboratories Catalog #HSP9601** with  PCR foil, DNase-/RNase-free **Sarstedt Catalog #95.1993** . Make sure that each well is properly sealed.

1m





Sealed PCR plate.

- 7.1 Close the extraction plate with PCR foil. Make sure each well is properly sealed.

1m



- 8 Briefly vortex and shortly spin down the sealed

1m

Hard-Shell 96-well PCR plate **Bio-Rad Laboratories Catalog #HSP9601** .



PCR

2h 2m

- 9 Place the sealed

1m

Hard-Shell 96-well PCR plate **Bio-Rad Laboratories Catalog #HSP9601** in a Thermocycler and close it properly.



- 9.1 Depending on the Thermocycler, a lid or pad should be placed on the PCR plate before closing the Thermocycler.

1m



- 10 Enter/select the correct program (depending on used primers and chemicals) and start the polymerase chain reaction.

2h



Note

PCR program for fwhF2/fwhR2n:

Step 1: 95°C for 5 minutes (initial denaturation)

Step 2: 95°C for 30 seconds (denaturation)

Step 3: 58°C for 90 seconds (annealing)

Step 4: 72°C for 30 seconds (elongation)

--- repeat for 20 cycles ---

Step 5: 68°C for 10 minutes (final elongation)

Note

For samples with low template amounts (e.g. eDNA from water), also 30 cycles of 1st step PCR can be used.

Storage

1m

- 11 Vortex and spin down the



Hard-Shell 96-well PCR plate **Bio-Rad Laboratories Catalog #HSP9601** to use it

for further processing, e.g. gel electrophoresis, and/or store it at 4°C (short-term) or -20°C (long term).

1m

