



Feb 27, 2023

Splenocyte Preparation for Immune Reconstitution (Humanisation)

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7y6e9gwz/v1>

Sandra Petrus-Reurer¹, Kourosh Saeb-Parsy¹

¹University of Cambridge



Sandra Petrus-Reurer

University of Cambridge

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7y6e9gwz/v1>

Protocol Citation: Sandra Petrus-Reurer, Kourosh Saeb-Parsy 2023. Splenocyte Preparation for Immune Reconstitution (Humanisation) . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7y6e9gwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 27, 2023



Last Modified: February 27, 2023

Protocol Integer ID: 77679

Keywords: Mouse immune reconstitution, Humanisation, Splenocytes, splenocyte preparation for immune reconstitution, splenocyte preparation, immune reconstitution, mice

Funders Acknowledgements:

MRC UKRMP

Grant ID: G101985

Abstract

This protocol guides through the steps required to humanise mice with human splenocytes



- 1 Thaw vial (approx. 50×10^6 cells/vial) in water bath 37C for 2-3min (leave crystal)
- 2 Add cells into 10mL/vial of warm 50% RPMI media + 50% FBS
- 3 Centrifuge 7min at 300g
- 4 Combine all vials in 10mL RPMI media + 10% FBS
- 5 Centrifuge 7min at 300g
- 6 Resuspend all vials in 10mL RPMI media + 10% FBS + add DNase-I (typically 40uL/mL if DNase-I 1mg/mL) at 37C (incubator) for 15-20min (depending on the amount of clumps seen). Check and flick the tube every 5 min
- 7 Centrifuge 7min at 300g
- 8 Resuspend in 10mL RPMI media + 10% FBS and filter with 70um filter (white)
- 9 Count cells:
-If with Neubauer chamber: resuspend 1:2 (10uL cells + 10uL Trypan Blue); Calculation: average of cells per quadrant x 2 (dil) x 10^4 (cells/mL) x 10 (mL)
-If automatic counter: resuspend 1:10 (10uL cells + 90uL Trypan Blue); Calculation: average of two reads and take live cells x 5 (dil) x 10 (mL)
- 10 Resuspend pellet in right volume to: 50×10^6 cells/mL (10×10^6 cells/200uL=50,000 cells/uL) for NSG mice, or for NSG-dKO mice 100×10^6 cells/mL (20×10^6 cells/200uL=100,000 cells/uL) with PBS + 2% FBS in universal tubes/epps
- 11 Freeze leftover cells (max 50×10^6 cells/vial) with 10% DMSO in FBS (1mL per vial)
- 12 Take to the animal facility x2 the number of cells and volume needed per donor (in case there are losses due to clumping, dead volumes etc.)



- 13 Inject 200uL/mouse intraperitoneally

- 14 Follow humanization weekly by performing tail vein bleeds with heparin (applying red cell lysis buffer x3, 10min on ice) and flow cytometry for the immune markers: hCD45, mCD45, hCD19, hCD3, hCD8, hCD4, 7AAD