

Feb 27, 2019

Splitting 96 Well Plates for gDNA Extraction and Continuing Culture

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.x78frw

Celeste M M. Karch¹, Rita Martinez¹, Jacob Marsh¹


¹Washington University in St Louis

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Celeste M M. Karch

Washington University in St Louis



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x78frw>

Protocol Citation: Celeste M M. Karch, Rita Martinez, Jacob Marsh 2019. Splitting 96 Well Plates for gDNA Extraction and Continuing Culture. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x78frw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 16, 2019

Last Modified: February 27, 2019

Protocol Integer ID: 20448

Keywords: plates for gdna extraction, gdna extraction, plate

Attachments



Comprehensive

Genomi...


31KB

Guidelines

This protocols is part of the **Screening Edited iPSC Clones collection**.

Troubleshooting

Safety warnings











 Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

Before start



It takes approximately 1 week for iPSC picked into 96 well plates to be sufficiently confluent for freezing and screening. For screening purposes, a fraction of the cells picked into one well of a 96 well plate will be saved for DNA Extraction and the remaining will be kept in culture or frozen down.

Split cells upon reaching maximum 80% confluence and minimum 40%



- 1 Coat 96 well plate with  50 μ L Matrigel per well.
- 2 Incubate at  37 °C for  01:00:00 .
- 3 Prepare plate for expansion by aspirating Matrigel from plate.
- 4 Add  50 μ L mTesR1 + 10 uM Rock Inhibitor to appropriate wells.
- 5 Aspirate media from original plate.
- 6 Wash with  200 μ L PBS and aspirate.
- 7 Add  25 μ L of 0.05% Trypsin.
- 8 Incubate at  37 °C for  00:05:00
- 9 Tap to lift cells from plate.
- 10 Check under microscope to ensure that cells have detached from plate.
- 11 Add  50 μ L FBS and tap to mix.
- 12 Transfer  50 μ L to 96 well PCR plate, while maintaining the location of each sample (this plate will be used for gDNA extraction).
- 13 Transfer remaining cells (~30uL) to 96 well plate containing mTesR1.



- 14 Incubate at  37 °C .
- 15 After  24:00:00 , complete daily media changes with mTesR1.