

Sep 08, 2020

SPOT2 protocol

DOI

dx.doi.org/10.17504/protocols.io.bk28kyhw

Huimin Zhao¹, stlane2¹, Guanhua Xun¹

¹University of Illinois at Urbana-Champaign

XPRIZE Rapid Covid Tes...



stlane2

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bk28kyhw

Protocol Citation: Huimin Zhao, stlane2, Guanhua Xun 2020. SPOT2 protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bk28kyhw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 08, 2020

Last Modified: September 08, 2020

Protocol Integer ID: 41792



Materials

STEP MATERIALS

⊗ Magnesium Sulfate (MgSO₄) Solution - 6.0 ml **New England Biolabs Catalog #B1003S**

⊗ E gene primer mix **Catalog #N/A**

⊗ Non-primer oligos **Catalog #N/A**

⊗ Bst 2.0 WarmStart DNA Polymerase - 8,000 units **New England Biolabs Catalog #M0538L**

⊗ WarmStart RTx Reverse Transcriptase - 250 rxns **New England Biolabs Catalog #M0380L**

⊗ Isothermal Amplification Buffer - 6.0 ml **New England Biolabs Catalog #B0537S**

⊗ Nuclease-free Water - 25 ml **New England Biolabs Catalog #B1500S**

⊗ Deoxynucleotide Solution Mix - 8 umol of each **New England Biolabs Catalog #N0447S**

⊗ N gene primer mix **Catalog #N/A**

⊗ Saliva sample **Catalog #N/A**

⊗ Manganese(II) chloride tetrahydrate **Sigma Aldrich Catalog #M3634**

⊗ Reporter probe 1 **Catalog #N/A**

⊗ Reporter Probe 2 **Catalog #N/A**

⊗ Non-CRISPR nuclease **Catalog #N/A**



Protocol materials

⊗ Non-primer oligos **Catalog #N/A**

⊗ WarmStart RTx Reverse Transcriptase - 250 rxns **New England Biolabs Catalog #M0380L**

⊗ Nuclease-free Water - 25 ml **New England Biolabs Catalog #B1500S**

⊗ Reporter Probe 2 **Catalog #N/A**

⊗ N gene primer mix **Catalog #N/A**

⊗ Saliva sample **Catalog #N/A**

⊗ Magnesium Sulfate (MgSO₄) Solution - 6.0 ml **New England Biolabs Catalog #B1003S**

⊗ E gene primer mix **Catalog #N/A**

⊗ Manganese(II) chloride tetrahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M3634**

⊗ Reporter probe 1 **Catalog #N/A**

⊗ Bst 2.0 WarmStart DNA Polymerase - 8,000 units **New England Biolabs Catalog #M0538L**

⊗ Isothermal Amplification Buffer - 6.0 ml **New England Biolabs Catalog #B0537S**

⊗ Deoxynucleotide Solution Mix - 8 umol of each **New England Biolabs Catalog #N0447S**

⊗ Non-CRISPR nuclease **Catalog #N/A**

⊗ Non-CRISPR nuclease **Catalog #N/A**

⊗ Reporter Probe 2 **Catalog #N/A**

⊗ WarmStart RTx Reverse Transcriptase - 250 rxns **New England Biolabs Catalog #M0380L**

⊗ Isothermal Amplification Buffer - 6.0 ml **New England Biolabs Catalog #B0537S**

⊗ Saliva sample **Catalog #N/A**

⊗ Nuclease-free Water - 25 ml **New England Biolabs Catalog #B1500S**

⊗ Manganese(II) chloride tetrahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M3634**

⊗ Non-primer oligos **Catalog #N/A**

⊗ Reporter probe 1 **Catalog #N/A**

⊗ Magnesium Sulfate (MgSO₄) Solution - 6.0 ml **New England Biolabs Catalog #B1003S**

⊗ E gene primer mix **Catalog #N/A**

⊗ Deoxynucleotide Solution Mix - 8 umol of each **New England Biolabs Catalog #N0447S**

⊗ N gene primer mix **Catalog #N/A**

⊗ Bst 2.0 WarmStart DNA Polymerase - 8,000 units **New England Biolabs Catalog #M0538L**



- 1 Label test kit bag with patient's name and a unique sample number for the test run. Provide a 20 uL microcap and PCR tube and have patient collect a saliva sample and place it into then seal the PCR tube. Label PCR tube with sample number matching the test kit bag.
- 2 Place PCR tubes into floating rack and insert into sous vide water bath at 95C. Incubate for 5 minutes for sample pretreatment.

Equipment

Thermal Immersion Cooker

NAME

Sous Vide - low cost waterbath

TYPE

Wancle

BRAND

P00252

SKU

🌡 95 °C

🕒 00:05:00

- 2.1 During incubation, scan barcodes on test kits and enter patient identifying information into the SPOT2 device.

Equipment

SPOT2 Device

NAME

Data logging fluorometer

TYPE

University of Illinois

BRAND

N/A

SKU

Equipment

| | |
|---------------------------|-------|
| USB Laser Barcode Scanner | NAME |
| Barcode Scanner | TYPE |
| WoneNice | BRAND |
| WN3300 | SKU |

- After completion of pretreatment, use the second provided 20 uL microcap to transfer a small volume of saliva from the PCR tube to the capillary tube, containing SPOT assay mastermix. Place capillary tubes into floating rack and insert into sous vide water bath at 63C for 30 minutes.

🌡️ 63 °C ⌚ 00:30:00

SPOT assay mastermix:

| | Initial concentration | Final concentration | Amount (μL) |
|------------------------------------|-----------------------|---------------------|-------------|
| Upper compartment | | | |
| WarmStart® Bst 2.0 | 8000 units/mL | 160 units/mL | 2 |
| WarmStart® RTx | 15,000 units/mL | 150 units/mL | 1 |
| Isothermal amplification buffer | 10X | 0.5X | 8 |
| dNTPs | 10 mM | 0.7 mM | 5.6 |
| MgSO4 | 100 mM | 4 mM | 3.2 |
| N gene primer mix | 10X | 0.25X | 2 |
| E gene primer mix | 10X | 0.25X | 2 |
| Saliva samples | | | 5 |
| Non-CRISPR nuclease | 5 mg/mL or 55 μM | 1.375 uM | 2 |
| MnCl2 | 50 mM | 0.5 mM | 0.8 |
| Non-primer oligos (total 6 oligos) | 100 μM | 625 nM | 3 |



| | | | |
|---------------------|-------------|---------------|--------|
| Reporter probe 1 | 100 μ M | 156.25 nM | 0.125 |
| Reporter probe 2 | 100 μ M | 312.5 μ M | 0.25 |
| Nuclease-free water | | | 44.025 |
| Total | | | 80 |



Bst 2.0 WarmStart DNA Polymerase - 8,000 units **New England Biolabs Catalog #M0538L**



2 μ L



WarmStart RTx Reverse Transcriptase - 250 rxns **New England Biolabs Catalog #M0380L**



1 μ L



Isothermal Amplification Buffer - 6.0 ml **New England Biolabs Catalog #B0537S**



8 μ L



Deoxynucleotide Solution Mix - 8 μ mol of each **New England Biolabs Catalog #N0447S**



5.6 μ L

[M] 0.7 millimolar (mM)



Magnesium Sulfate (MgSO₄) Solution - 6.0 ml **New England Biolabs Catalog #B1003S**



3.2 μ L

[M] 4 millimolar (mM)



N gene primer mix **Catalog #N/A**



E gene primer mix **Catalog #N/A**



Saliva sample **Catalog #N/A**



Nuclease-free Water - 25 ml **New England Biolabs Catalog #B1500S**



Non-CRISPR nuclease **Catalog #N/A**



Manganese(II) chloride tetrahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M3634**





0.8 μ L

[M] 0.5 millimolar (mM)



Non-primer oligos **Catalog #N/A**

 Reporter probe 1 **Catalog #N/A**

 Reporter Probe 2 **Catalog #N/A**

- 4 Transfer floating capillary rack to sous vide at 95C for 5 minutes.

 95 °C  00:05:00

- 5 Individually transfer capillary tubes to SPOT2 device and press "Run" on the SPOT software. An integrated RFID chip on each capillary will automatically link test result to the patient information entered previously.

Equipment

| | |
|--------------------------|-------|
| SPOT2 Device | NAME |
| Data logging fluorometer | TYPE |
| University of Illinois | BRAND |
| N/A | SKU |