


Sep 01, 2023

🌐 Staining of GFAP and IBA1 in PDGFR-B/Td-tomato brain sections

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn36mml5d/v1>



Daniel Manrique-Castano¹

¹Université Laval



Daniel Manrique-Castano

Université Laval

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn36mml5d/v1>

Protocol Citation: Daniel Manrique-Castano 2023. Staining of GFAP and IBA1 in PDGFR-B/Td-tomato brain sections. protocols.io <https://dx.doi.org/10.17504/protocols.io.14egn36mml5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 01, 2023

Last Modified: September 01, 2023

Protocol Integer ID: 87251

Keywords: Staining, GFAP, IBA1, Immunofluorescence, iba1 in pdgfr, tomato brain sections this protocol, tomato brain section, mouse brain section, fixed mouse brain section, pdgfr, gfap, staining, iba1

Abstract

This protocol is suitable for the staining of GFAP and IBA1 in PDGFR-B/Td-Tomato fixed mouse brain sections.

Guidelines

Read the entire protocol before starting the procedure.

Note that this protocol uses 3 hours (room temperature) incubation for primary antibodies.

Do not submit the tissue to antigen retrieval or other acidic solutions to preserve the Td-tomato signal.

Works with the listed antibodies and dilutions. For other references, preliminary tests are highly recommended.

Materials

Equipment

ImmEdge® Hydrophobic Barrier PAP Pen

NAME

Vector

BRAND

H-4000

SKU

<https://vectorlabs.com/products/histology/immedge-hydrophobic-barrier-pen>

LINK

Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization

SPECIFICATIONS

Permeabilization solution = [M] 0.3 % (v/v) Triton

⊗ Triton X-100 Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML / = [M] 0.1 % (v/v)

⊗ TWEEN 20 Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7949 /

[M] 0.3 Molarity (M) ⊗ Glycine in 1x PBS

Antibody buffer = [M] 0.1 % (v/v) ⊗ TWEEN 20 Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7949 /

[M] 1 % (v/v) ⊗ Normal Donkey Serum Merck MilliporeSigma (Sigma-Aldrich) Catalog #D9663-10ml

[M] 1 Mass / % volume

⊗ Bovine Serum Albumin (BSA) Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7906

Primary antibodies:

⊗ GFAP Monoclonal Antibody (2.2B10) Thermo Fisher Scientific Catalog #13-0300

⊗ Anti Iba1 Rabbit antibody FUJIFILM Wako Pure Chemical Corporation Catalog #019-19741

Secondary antibodies:

⊗ Alexa Fluor® 647 AffiniPure Donkey Anti-Rat IgG (H L) Jackson ImmunoResearch Laboratories, Inc. Catalog #712-605-153

⊗ Anti-Rabbit Invitrogen - Thermo Fisher Catalog #A21206



⊗ DAPI (46-Diamidino-2-Phenylindole Dilactate) **Invitrogen - Thermo Fisher Catalog #D3571**

Mounting media:

⊗ Fluoromount-G **Electron Microscopy Sciences Catalog #17984-25**

Protocol materials

⊗ Normal Donkey Serum **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D9663-10ml**

⊗ Fluoromount-G **Electron Microscopy Sciences Catalog #17984-25**

⊗ Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML**

⊗ Anti-Rabbit **Invitrogen - Thermo Fisher Catalog #A21206**

⊗ DAPI (46-Diamidino-2-Phenylindole Dilactate) **Invitrogen - Thermo Fisher Catalog #D3571**

⊗ TWEEN 20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7949**

⊗ Bovine Serum Albumin (BSA) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7906**

⊗ Alexa Fluor® 647 AffiniPure Donkey Anti-Rat IgG (H L) **Jackson ImmunoResearch Laboratories, Inc. Catalog #712-605-153**

⊗ Anti Iba1 Rabbit antibody **FUJIFILM Wako Pure Chemical Corporation Catalog #019-19741**

⊗ GFAP Monoclonal Antibody (2.2B10) **Thermo Fisher Scientific Catalog #13-0300**

⊗ Glycine

Safety warnings

! DAPI is highly toxic. Handle it with care.



Tissue preparation and blocking

20m

- 1 Take out sections from -80 and place them in an incubator/plate for 00:20:00 at 37 °C . This step is performed to ensure tissue attachment to the crystal slides. 20m
- 2 Draw a hydrophobic barrier on each slide using **ImmEdge® Hydrophobic Barrier PAP Pen** and let dry for about 00:10:00 minutes. This will prevent buffers from escaping the tissue area. 10m

Note

There are other hydrophobic pens on the market. However, this one is recommended for its quality, consistency, and durability.

- 3 To initially rehydrate and permeabilize the tissue, place the slides in a slide jar containing the **permeabilization solution** with shaking for 00:30:00 30m

Note

The proposed permeabilization solution contains **glycine**, which is suitable for reducing PFA-related autofluorescence, especially in the 488 channel.

Antigen retrieval or any treatment with acid is not recommended given that it destroys the endogenous Td-tomato signal from PDGFR-B+ cells.

- 4 To prevent unspecific binding, incubate brain sections in **permeabilization solution** containing 5 % (v/v) 1h



Normal Donkey Serum **Merck MilliporeSigma (Sigma-Aldrich)** Catalog #D9663-10ml

and 1 Mass / % volume



Bovine Serum Albumin (BSA) **Merck MilliporeSigma (Sigma-Aldrich)** Catalog #A7906

for 01:00:00 at Room temperature .



Note

Please note that blocking serum must be chosen according to **secondary antibody species**. For the procedure depicted in this protocol, all secondary antibodies are raised in donkey. However, goat secondary antibodies will also be suitable.

Additionally, consider **avoiding BSA** when **primary antibodies** come from goat or sheep. BSA can generate unspecific background.

Antibody incubation

3h

- 5 When blocking is finished, decant the buffer (no washing is required) and incubate **primary antibodies** in **antibody buffer** for  03:00:00 at  Room temperature according to **table 1**.

3h



	A	B	C	D	E
	Antibody	Company	Reference	Specie	Dilution
	Gfap	Invitrogen	13-0300	Rat	1:500
	Iba1	Wako	019-19741	Rabbit	1:500

Table 1. Primary antibodies

Note

Please note that antibody buffers **do not contain Triton X-100**. This detergent tends to break the hydrophobic barrier and the sections may not be adequately incubated. In addition, given that the sections are already permeabilized, we do not contemplate the use of this detergent.

Please note that the proposed antibody incubation lasts only **3 hours**. Although in most cases staining protocols recommend overnight incubation for primary antibodies, the tests performed in our lab showed better efficiency at room temperature.

The indicated **Iba1 antibody** is superior to others available in the market.

6 When primary antibody incubation is finished, wash the sections  00:05:00 5x with  0.05 % (v/v)  TWEEN 20 Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7949 in **PBS**.

5m



7 Incubate **secondary antibodies** in  0.1 % (v/v)  TWEEN 20 Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7949 for  01:00:00 at  Room temperature according to **table 2**.

1h



	A	B	C	D	E
	Antibody	Company	Reference	Channel	Dilution
	Donkey anti-rat	Jackson	712-605-153	647	1:500
	Donkey anti-rabbit	Invitrogen	A21206	488	1:500
	Dapi	Invitrogen	D3571	405 (Dapi)	1:5000

Table 2. Secondary antibodies

8 When secondary antibody incubation is finished, wash the sections  00:05:00 3x with  0.05 % (v/v)  TWEEN 20 Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7949 in **PBS**. Follow this with  00:05:00 2x washes with **PBS** to remove all detergent traces.

10m



9 Clean the remaining buffer on the slides using absorbent tissue and mount the sections with a drop of  Fluoromount-G Electron Microscopy Sciences Catalog #17984-25 .