

Feb 08, 2024

Version 2

Standard DAB Staining for Free-floating Fixed NHP Brain Tissue V.2

DOI

dx.doi.org/10.17504/protocols.io.261ged767v47/v2

Andreea Bostan¹

¹Department of Neurobiology, University of Pittsburgh



Andreea Bostan

University of Pittsburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ged767v47/v2>

Protocol Citation: Andreea Bostan 2024. Standard DAB Staining for Free-floating Fixed NHP Brain Tissue . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ged767v47/v2> Version created by **Andreea Bostan**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 08, 2024



Last Modified: May 31, 2024

Protocol Integer ID: 94869

Keywords: ASAPCRN, Immunostaining, DAB, NHP Brain Tissue, fixed nhp brain tissue, nhp brain tissue, brain tissue, fixed brain tissue section, brain tissue section, biotin abc complex, avidin, tissue, standard dab



Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020519

Abstract

This protocol details the procedure for immunohistochemical 3,3'-Diaminobenzidine (DAB) staining of free-floating fixed brain tissue sections using the avidin/biotin ABC complex.

This protocol has been tested with free-floating non-human primate (NHP) and rodent (mouse, rat) brain tissue that has been fixed (10% formalin or 4% paraformaldehyde), cryoprotected (sucrose or glycerol gradients), and cryo-sectioned  20 μm -  50 μm .

Guidelines

When using 6 well tissue culture plates [Falcon, 353046] to react individual sections, you will need **2+ mL** solutions for ***each*** well plate.

When using circular staining nets [e.g., Brain Research Laboratories #4115] to react multiple series of sections, you will need **50 mL** solutions for ***each***.



Materials

Tissue:

Brain tissue sections (20 - 50 μm).

Materials/Equipment:

- Tissue culture plates or circular staining nets
- Orbital shaker
- Fume hood
- Nitrile Gloves
- Glass slides (charged or subbed)

Reagents:

- Phosphate-buffered saline (PBS)
- Hydrogen Peroxide: H_2O_2 (3% or 30%)
- Distilled water: dH_2O
- Primary Antibody
- Secondary Antibody (to match the host of the primary antibody)
- Normal Serum Blocking Solution (e.g., Normal Horse Serum, S-2000-20; Normal Goat Serum, S-1000-20)
- Vectastain Elite ABC Peroxidase Kit (Standard) (PK-6100) (Vector Laboratories)
- ABC-HRP Kit

Examples:

Vectastain ABC-HRP Kit, Peroxidase (Mouse IgG) (PK-4002, Vector Laboratories)

Vectastain ABC-HRP Kit, Peroxidase (Rabbit IgG) (PK-4001, Vector Laboratories)

- DAB Substrate Kit

Examples:

Peroxidase (HRP) with Nickel (3,3'-diaminobenzidine) (SK-4100) (Vector Laboratories)

ImmPACT DAB (SK-4105)

Troubleshooting



Safety warnings

- ⚠ Use appropriate care when using hydrogen peroxide (reactive, can cause skin/eye damage) and DAB (suspected carcinogen). Collect DAB solution for chemical waste disposal.



Part I (Day 1)

3h

- 1 Bring tissue to Room temperature in buffer (e.g., Phosphate buffered saline, PBS) on an orbital shaker for 30 minutes. 00:30:00 .

30m

- 2 Prepare **Peroxide Solution (0.3 - 3 % H₂O₂)** in dH₂O.

5m

E.g., for 10 mL 0.3% H₂O₂ use:

- 100 µL 30% H₂O₂
- 9900 µL dH₂O

- 3 Prepare **Blocking Serum Solution** (e.g. Normal Horse, Normal Goat Serum) using a serum that matches the **host of the secondary antibody** (e.g. Normal Horse Serum for a Horse anti-Mouse secondary, Normal Goat Serum for a Goat anti-Rabbit secondary).

5m

E.g., in 10 mL buffer (PBS) add:

- 150 µL normal serum (or 3 drops of normal serum if using an ABC kit, e.g. Vectastain ABC-HRP Kit, Peroxidase Mouse IgG PK-4002, Rabbit IgG PK-4001)

- 4 Prepare **Primary Antibody Solution** at the appropriate dilution in buffer (e.g., 1:1000 in PBS).

5m

- 5 **Rinse** in buffer (e.g. PBS) on a shaker at Room temperature : **3 × 3-5 minutes.**

15m

00:03:00 - 00:05:00

- 6 Quench endogenous peroxide in **Peroxide Solution (0.3 - 3 % H₂O₂)** on a shaker at

1h

Room temperature : **30 - 60 minutes.** 00:30:00 - 01:00:00

- 7 **Rinse** in buffer (e.g. PBS) on a shaker at Room temperature : **3 × 3-5 minutes.**

15m

00:03:00 - 00:05:00

- 8 Incubate in **Blocking Serum Solution** on a shaker at RT: **1 hour.**

1h

DO NOT RINSE after blocking serum.



- 9 Incubate in **Primary Antibody Solution** on a shaker at 4 °C Overnight , or 20h
longer (20 - 72 hours depending on the antibody).

Part II (Day 2)

4h

- 10 Bring tissue (in the **Primary Antibody Solution**) to Room temperature on a shaker 30m
(**30 - 60 minutes**). 00:30:00 - 01:00:00
- 11 Prepare **ABC Solution** in buffer (e.g. PBS) (**at least 30 minutes before use**). 5m
 00:30:00 .
- 12 Prepare **Secondary Antibody Solution (1:200)** in buffer (e.g. PBS). 5m
In 10 mL buffer add:
- 150 µL (= 3 drops of normal serum from a Vector Labs kit) of normal serum (matched to the host of your secondary antibody)
 - 50 µL (= 1 drop secondary antibody from a Vector Labs kit) of biotinylated secondary antibody (matched to the host of your primary Antibody)
- 13 **Rinse** in buffer (e.g. PBS) on a shaker at Room temperature : **3 × 3-5 minutes**. 15m
 00:03:00 - 00:05:00 .
- 14 Incubate in **Secondary Antibody Solution** on a shaker at Room temperature : **30 minutes**. 30m
 00:30:00 .
- 15 **Rinse** in buffer (e.g. PBS) on a shaker at Room temperature : **3 × 3-5 minutes**. 15m
 00:03:00 - 00:05:00 .
- 16 Incubate in **ABC Solution** on a shaker at Room temperature : **60 minutes**. 1h
 01:00:00 .



17 **Rinse** in buffer (e.g. PBS) on a shaker at Room temperature : **3 × 3-5 minutes.**

15m

00:03:00 - 00:05:00 .

18 Prepare **Peroxide Substrate Solution** in dH₂O.

5m

To use the Vector Labs DAB Peroxidase Substrate Kit (SK-4100):

In 5 mL dH₂O:

- 2 drops Reagent 1
- 4 drops Reagent 2
- 2 drops Reagent 3
- [optional] 2 drops of Reagent 4 (Nickel) if a black reaction product is desired

Note: Mix well before use. Use immediately.

19 Incubate in **Peroxide Substrate Solution on** a shaker at Room temperature :

6m

00:03:00 - 00:06:00 .

Note: Watch the tissue closely to avoid high background staining.

20 **Rinse** in buffer (e.g. PBS) on a shaker at Room temperature : **3 × 3-5 minutes.**

15m

00:03:00 - 00:05:00 .

21 Mount tissue on glass slides (subbed or charged) in 1:8 buffer in dH₂O and let air dry.

22 Rinse slides with dH₂O and let air dry (preferably in a hood).

23 Coverslip clean and dry slides with Cytoseal 60 (Thermo Fisher #830-16).

Protocol references

https://vectorlabs.com/productattachments/protocol/VL_SK-4100_UserGuide_LBL02267.pdf

<https://vectorlabs.com/products/vectastain-elite-abc-hrp-kit-standard>