

Mar 28, 2025

Standard Protocol for Magnetic Spheroid Formation

 [Methods and Protocols](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4z772vo5/v1>

Marine Simonneau¹

¹LabTAU



magali perier

INSERM U1032

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4z772vo5/v1>

External link: <https://doi.org/10.3390/mps8040075>

Protocol Citation: Marine Simonneau 2025. Standard Protocol for Magnetic Spheroid Formation . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.3byl4z772vo5/v1>

Manuscript citation:

Perier M, Wang L, Simonneau M, Ngo-Reymond J, Guillermet-Guibert J, Lafond M, Lafon C (2025) Formation and Characterization of Two Magnetic Three-Dimensional Spheroid Models of Murine Pancreatic Adenocarcinoma. Methods and Protocols 8(4). doi: [10.3390/mps8040075](https://doi.org/10.3390/mps8040075)



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 28, 2025

Last Modified: March 28, 2025

Protocol Integer ID: 125617

Keywords: pancreatic adenocarcinoma spheroid model, magnetic spheroid formation, standard protocol for magnetic spheroid formation, pancreatic cancer cell, spheroid formation, using magnetic aggregation, magnetic aggregation, spheroid, fibroblast

Abstract

Our study introduces three-dimensional (3D) pancreatic adenocarcinoma spheroid models using magnetic aggregation of pancreatic cancer cells and immortalized fibroblasts



Materials

- KPC A219 Cells
- iMEF Cells
- Nanoshuttles (Greiner Bio-One 657846)
- Complete DMEM/F-12 Culture Medium (10% FBS, 1% L-Glutamine, 1% Penicillin/Streptomycin) (Gibco 21331-020)
- DPBS (Gibco 14190-086)
- Trypsin (Gibco 25300-054)
- TGF β Solution (Transforming Growth Factor) (0.01 mg/ml) (PEPROTECH 100-21)
- 96-Magnet Plate (Greiner Bio-One 130188)
- 24-Magnet Plate
- Medium Aspirator
- Tube Holder
- Malassez/Thomas Hemocytometer
- Blue Trypan
- Micropipettes 20 μ L and 100 μ L (STAR LAB)
- Pipetboy
- T75 Culture Flask (Corning 430641U) / T25 Culture Flask (Thermo Fisher Scientific 156367)
- 96-Well Transparent Plate with Flat Bottom (Greiner



Bio-One 655970) / 96-Well Transparent Plate with Round Bottom (Greiner Bio-One 650970)

- 15ml Centrifuge Tubes (Fisher Scientific 431885) / 50ml Centrifuge Tubes (Fisher Scientific 431176)

Day -3 Cell Passage

- 1 ●
Warm the DMEM culture medium, DPBS, and trypsin in a water bath.
- 2 ●
Aspirate the medium from the flasks.
●
Add 5 mL of DPBS to the T75 flask and 3 mL of DPBS to the T25 flask to rinse the cells.
● Aspirate the DPBS.
- 3 ●
Add 2 mL of trypsin to the T75 flask and 1 mL of trypsin to the T25 flask.
●
Incubate KPC A219 cells for 6 minutes and iMEF cells for 5 minutes in the incubator (37°C, 5% CO₂).
●
Check if the cells detach properly from the flask walls.
●
Add 10 mL of DMEM to the T75 flask and 5 mL of DMEM to the T25 flask to inhibit the trypsin reaction.
- 4 ●
Take 20 µL of cell suspension and mix it with 20 µL of Blue Trypan.
●
Distribute the mixed solution on the Malassez/Thomas cell counter, then count the cells.
- 5 ●
Centrifuge the cell suspensions for 5 minutes at 300g in a centrifuge at room temperature.

Preparation of Cell Suspensions at 2×10^6 Cells/mL

- 6 ●
Calculate the volume of DMEM needed to dilute or concentrate the cell suspensions to a concentration of 2×10^6 cells/mL:

Volume of DMEM to add = Total number of cells / 2×10^6

- 7 ●
Aspirate the supernatant and add the calculated volume of DMEM to achieve a cell suspension at 2×10^6 cells/mL.
- Homogenize the cell suspension by pipetting.

Preparation of New T75 or T25 Flasks for KPC A219/iMEF Cell Culture

- 8 ●
To form 0-40 spheroids: use 1 T25 flask of KPC and 1 T25 flask of iMEF.
- To form 40-100 spheroids: use 1 T25 flask of KPC and 1 T75 flask of iMEF / 2 T25 flasks of iMEF.
- 9 ●
For the T75 flask: add 15 mL of DMEM + 300 μ L of cell suspension at 2×10^6 cells/mL.
- For the T25 flask: add 5 mL of DMEM + 100 μ L of cell suspension at 2×10^6 cells/mL.

Cell Incubation

- 10 ●
Observe the plate under the optical microscope.
- Incubate all flasks in the incubator (37°C, 5% CO₂).

Day -1: Addition of TGF Beta and Nanoshuttles

- 11 **Addition of TGF β to iMEF Flask 24 Hours Before Co-culture** ●
For the T75 iMEF flask, add 75 μ L of TGF β to achieve 50 ng/mL of TGF β in the culture medium.
- For the T25 iMEF flask, add 25 μ L of TGF β to achieve 50 ng/mL of TGF β in the culture medium.
- 12 **Addition of Nanoshuttles to KPC A219 and iMEF Flasks the Day Before Co-culture** ●
For KPC A219: add 100 μ L of Nanoshuttles to the T25 flask.



For iMEF: add 50 μL of Nanoshuttles to the T25 flask and 150 μL of Nanoshuttles to the T75 flask.

Day 0: Co-culture and Spheroid Formation

13 **Cell Passage and Preparation of Cell Suspensions at 2×10^6 Cells/mL Concentration**



Perform the cell passage and record the concentrations and cell counts.

14 **Cell Count**

Prepare the cell suspensions of KPC and iMEF at 2×10^6 cells/mL concentration and record the volumes of the cell solution.

15 **Preparation of KPC and iMEF Mixing Solution**

Calculate the volume of KPC suspension, iMEF suspension, and DMEM culture medium needed to form spheroids:

1 spheroid = 10,000 KPC + 20,000 iMEF



Mix the KPC suspension, iMEF suspension, and DMEM after the calculations.

16



Dispense 150 μL of the mixed solution into each well of the 96-well plate.

17 **Incubate**

the cells in the 96-magnet plate in the incubator (37°C, 5% CO₂) for 6 hours, then remove the magnet plate.