



Apr 03, 2025

Standard Protocol for Spheroid Formation in the 96-Well Plate by centrifugation

 Methods and Protocols

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v97m6qg3e/v1>

Jacqueline Ngo-Reymond¹

¹INSERM U1032



magali perier

INSERM U1032

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v97m6qg3e/v1>

External link: <https://doi.org/10.3390/mps8040075>

Protocol Citation: Jacqueline Ngo-Reymond 2025. Standard Protocol for Spheroid Formation in the 96-Well Plate by centrifugation. protocols.io <https://dx.doi.org/10.17504/protocols.io.x54v97m6qg3e/v1>

**Manuscript citation:**

Perier M, Wang L, Simonneau M, Ngo-Reymond J, Guillermet-Guibert J, Lafond M, Lafon C (2025) Formation and Characterization of Two Magnetic Three-Dimensional Spheroid Models of Murine Pancreatic Adenocarcinoma. *Methods and Protocols* 8(4). doi: [10.3390/mps8040075](https://doi.org/10.3390/mps8040075)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 03, 2025

Last Modified: April 03, 2025

Protocol Integer ID: 126094

Keywords: spheroid production by centrifugation, standard protocol for spheroid formation, spheroid formation, well plate by centrifugation, spheroid production, centrifugation, spheroid, well plate

Abstract

This protocol enables spheroid production by centrifugation.

Materials

KPC A219 Cells

iMEF Cells

KPC A219 and iMEF are cultured separately in T75 flask in complete DMEM/F12 in the incubator (37°C, 5% CO₂). They have two passages in a week.

Complete DMEM/F-12 Culture Medium (10% FBS, 1% L-Glutamine, 1% Penicillin/Streptomycin) (Gibco 21331-020)
DPBS (Gibco 14190-086)

Trypsin (Gibco 25300-054)

TGFβ Solution (Transforming Growth Factor) (0.01 mg/ml) (PEPROTECH 100-21)

Medium Aspirator

Tube Holder

Malassez/Thomas Hemocytometer

Blue Trypan

Micropipettes

20μL and 200μL (STAR LAB)

Pipetboy

T75 Culture Flask (Corning 430641U) / T25 Culture

Flask (Thermo Fisher Scientific 156367)

96-Well Transparent Plate with Flat Bottom (Greiner

Bio-One 655970) / 96-Well Transparent Plate with Round Bottom (Greiner Bio-One 650970)

15ml Centrifuge Tubes (Fisher Scientific 431885) / 50ml

Centrifuge Tubes (Fisher Scientific 431176)

Day -1 before cell passage : Addition of TGF Beta and Nanoshuttles

- 1 For the T75 iMEF flask, add 75 μ L of TGF β to achieve 50 ng/mL of TGF β in the culture medium.
For the T25 iMEF flask, add 25 μ L of TGF β to achieve 50 ng/mL of TGF β in the culture medium.
- 2 Incubate the flasks in the incubator (37°C, 5% CO₂).

Day 0: Co-culture and Spheroid Formation

3 Cell Passage and Preparation of Cell Suspensions at 2×10^6 Cells/mL Concentration

- 3.1 Warm the DMEM culture medium, DPBS, and trypsin in a water bath.
- 3.2 Aspirate the medium from the flasks.
Add 5 mL of DPBS to the T75 flask and 3 mL of DPBS to the T25 flask to rinse the cells.
Aspirate the DPBS.
- 3.3 Add 2 mL of trypsin to the T75 flask and 1 mL of trypsin to the T25 flask.
Incubate KPC A219 cells for 6 minutes and iMEF cells for 5 minutes in the incubator (37°C, 5% CO₂).
Check if the cells detach properly from the flask walls.
Add 10 mL of DMEM to the T75 flask and 5 mL of DMEM to the T25 flask to inhibit the trypsin reaction.
- 3.4 Take 20 μ L of cell suspension and mix it with 20 μ L of Blue Trypan.
Distribute the mixed solution on the Malassez/Thomas cell counter, then count the cells.
- 3.5 Centrifuge the cell suspensions for 5 minutes at 300g in a centrifuge at room temperature.
- 3.6 Calculate the volume of DMEM needed to dilute or concentrate the cell suspensions to a concentration of 2×10^6 cells/mL:
Volume of DMEM to add = Total number of cells / 2×10^6
Aspirate the supernatant and add the calculated volume of DMEM to achieve a cell suspension at 2×10^6 cells/mL.
Homogenize the cell suspension by pipetting.



Preparation of KPC and iMEF Mixing Solution

- 4 Calculate the volume of KPC suspension, iMEF suspension, and DMEM culture medium needed to form spheroids:
1 spheroid = 10,000 KPC + 20,000 iMEF
Mix the KPC suspension, iMEF suspension, and DMEM after the calculations.

Dispensing the Mixed Solution into the 96-Well Plate

- 5 Dispense 150 μ L of the mixed solution into each well of the 96-well plate round bottom.

Plate centrifugation

- 6 Change the rotor adapted to plate
Equilibrate with another plate with the same volume
Centrifuge at room temperature 4 minutes, at 300G
Rotate the plate in 180°
Centrifuge again at room temperature 4 minutes, at 300G

Cell Incubation

- 7 Incubate the cells in the incubator (37°C, 5% CO₂)
Perform treatment 4 days after