

Jun 16, 2022

## Static insulin secretion analysis of isolated islets

DOI

[dx.doi.org/10.17504/protocols.io.bp2l61dxkvqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l61dxkvqe/v1)

Julien Ghislain<sup>1</sup>, Vincent Poitout<sup>2</sup>, Caroline CT Tremblay<sup>1</sup>, Marine Croze<sup>1</sup>

<sup>1</sup>Montreal Diabetes Research Center, CRCHUM, Montréal, QC, Canada;

<sup>2</sup>Montreal Diabetes Research Center, CRCHUM, and Department of Medicine, Université de Montréal, Montréal, QC, Canada

Laboratory of Vincent Poitout

Tech. support email: [julien.ghislain.chum@ssss.gouv.qc.ca](mailto:julien.ghislain.chum@ssss.gouv.qc.ca)



Julien Ghislain

CRCHUM

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l61dxkvqe/v1>

**Protocol Citation:** Julien Ghislain, Vincent Poitout, Caroline CT Tremblay, Marine Croze 2022. Static insulin secretion analysis of isolated islets . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l61dxkvqe/v1>

**Manuscript citation:**

Free fatty acid receptor 4 inhibitory signaling in delta cells regulates islet hormone secretion in mice. Croze ML, Flisher MF, Guillaume A, Tremblay C, Noguchi GM, Granziera S, Vivot K, Castillo VC, Campbell SA, Ghislain J, Huising MO, Poitout V. Mol Metab. 2021 Mar;45:101166. doi: 10.1016/j.molmet.2021.101166. PMID: 33484949

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 29, 2022

**Last Modified:** June 16, 2022

**Protocol Integer ID:** 60048

**Keywords:** pancreatic islet, somatostatin, insulin, glucose-stimulated insulin secretion, static insulin secretion analysis, hour assay from isolated pancreatic islet, insulin secretion, secreted insulin, intracellular insulin content, isolated pancreatic islet, isolated islet, insulin, islet number, suitable for islet, islet, assessing sst secretion, islet number per well, glucose, sst secretion, hour assay

**Funders Acknowledgements:**

**NIH**

Grant ID: 5R01DK058096-15

**CIHR**

Grant ID: 0406014428

**NSERC**

Grant ID: 51047

## Abstract

This protocol describes the steps to measure insulin secretion in a static, 1-hour assay from isolated pancreatic islets. It is suitable for islets isolated from both rodents and humans. We routinely apply this protocol to assess beta-cell function in response to glucose but can be easily adapted to interrogate the response to a variety of secretagogues (eg. fatty acids, hormones). Briefly batches of 10 islets are pre-incubated in triplicate in KRB solution at 2.8 mM glucose twice for 20 min followed by incubation in either 2.8 mM or 16.7 mM glucose for 1 hour. Secreted insulin is measured in the supernatant and intracellular insulin content, after acid-alcohol extraction, by radioimmunoassay. This protocol is also suitable for assessing SST secretion, however we recommend increasing the islet number per well from 10 to at least 20 due to the relative lower levels of SST compared to insulin.



## Materials

☒ Sodium Chloride **Fisher Scientific Catalog #S271**

☒ Calcium Chloride Dihydrate **Fisher Scientific Catalog #C79**

☒ Potassium Phosphate dibasic (KH<sub>2</sub>PO<sub>4</sub>) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9791**

☒ Potassium Chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #746436**

☒ Magnesium sulfate heptahydrate (MgSO<sub>4</sub>) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #63138**

☒ Sodium bicarbonate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S6014**

☒ HEPES **Merck MilliporeSigma (Sigma-Aldrich) Catalog #H6147**

☒ D-(+)-Glucose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G-7528**

☒ Fatty Acid Free Heat Shock Bovine Serum Albumin Powder **equitech bio, inc. Catalog #BAH66-0500**

☒ Ethanol (100%, Molecular Biology Grade) **Fisher Scientific Catalog #BP2818500**

☒ RPMI 1640 Medium **Thermo Fisher Scientific Catalog #11875093**

☒ FCS (Fetal Calf Serum) **Life Technologies**

☒ Rat insulin RIA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #RI-13K**

☒ Human Insulin-Specific RIA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #HI-14K**

## Troubleshooting



## Preparation of KRB solution and plates

### 1 Prepare KRB stocks

KRB stock I: weigh 27.7 g NaCl bring to 1 L with milliQ water

KRB stock II: weigh 1.494 g  $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$  and bring to 1 L with milliQ water

KRB stock III: weigh 0.648 g  $\text{KH}_2\text{PO}_4$  and bring to 1 L with milliQ water

KRB stock IV: weigh 1.415 g KCl, 1.17 g  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$  and 8.52 g  $\text{NaHCO}_3$  and bring to 1 L with milliQ water

1 Molarity (M) Glucose: weigh 18.02 g and bring to 100 mL with milliQ water. Filter sterilize using a 0.2  $\mu\text{m}$  filter.

### 2 Prepare KRB solution

Determine the number of static conditions for the assay in order to prepare a sufficient volume of KRB. Remember that you will have two pre-incubation steps and the picking, along with extra media to wash between these steps. In a beaker combine equal volumes of the four KRB stock solutions to achieve the desired volume. Add 2.38 mg/mL

HEPES powder and swirl to dissolve. Then add 1 mg/mL BSA (fatty acid free), but do not mix as the BSA will stick to the sides. Cover with plastic wrap (put holes in top) and place in the 37 °C incubator for > 01:00:00. Adjust the solution to

7.35 using 1 Molarity (M) NaOH. The solution should start at ~ 7.2.

### 3 2.8 mM Glucose condition and islet picking and washing

Calculate the required volume of 2.8 millimolar (mM) Glucose in KRB for the two islet pre-incubations and static ( 1 mL /well) plus 1 mL /well and about

70 mL for the islet picking and wash steps. Add 2.8  $\mu\text{L}$  of 1 Molarity (M) Glucose/ml of KRB.

#### 16.7 mM Glucose condition

Calculate the required volume of 16.7 millimolar (mM) Glucose in KRB for the static ( 1 mL /well) and add 16.7  $\mu\text{L}$  of 1 Molarity (M) Glucose/ml of KRB needed.



#### Note





Often there are other reagents to be added to the final static conditions, such as fatty acids, inhibitors or agonists. These additional components may require that separate KRB solutions be prepared. In the case of fatty acids addition prepare the KRB solution without BSA.




1h



#### 4 Prepare plates







Prepare islet picking plates by adding  1 mL of  2.8 millimolar (mM) Glucose in KRB to three wells of a 24-well plate for each static sample.

Prepare pre-incubation plates by adding  1 mL of  2.8 millimolar (mM) Glucose in KRB to three wells of a 24-well plate for each static sample. Repeat this step for a second pre-incubation plate. Place these plates in an incubator with  5 % volume CO<sub>2</sub> at  37 °C .



Prepare static incubation plate by adding  1 mL of experimental KRB to three wells of a 24-well plate for each static sample. Place these plates in the incubator with  5 % volume CO<sub>2</sub> at  37 °C .



### Islet Picking and Incubations



2h 40m






- 5 Following isolation, the islets should be allowed to recover in recovery medium (RPMI /  10 % (v/v) serum /  11.1 millimolar (mM) glucose) for  01:00:00 at  37 °C . Wash the islets in a petri dish containing  20 mL of  2.8 millimolar (mM) Glucose in KRB.


2h 50m






Pick the islets (in triplicate batches of 10) into the islet picking plate wells. Using a pipette transfer the islets from the picking plate to the first pre-incubation plate. Incubate at  37 °C for  00:20:00 .

Then transfer the islets to the second pre-incubation plate. Incubate at  37 °C for  00:20:00 .




Then transfer the islets to the incubation plate and incubate at  37 °C for  01:00:00 .

During the incubation, label two  1.5 mL tubes per sample for collection of the KRB media containing secreted insulin. Prepare and label  1.5 mL tubes filled with  1 mL acidified ethanol (  75 % (v/v) ethanol /  1.5 % (v/v) HCl) for insulin content analysis.

At the end of the static incubation, collect the islets and transfer them to the tubes pre-filled with acidified ethanol. Cap the vials and store at  -20 °C overnight.

Then, transfer the media from each well of the static incubation into a  1.5 mL tube, centrifuge at  4000 rpm, 4°C,  00:05:00 , transfer the supernatant to a new  1.5 mL tube and store at  -20 °C until ready to complete the insulin assay.



The next day retrieve the insulin content analysis tubes, vortex and centrifuge at  4000 rpm, 4°C, 00:05:00 . Transfer the supernatant to labelled  1.5 mL tubes. Store the insulin content samples at  -20 °C , until ready to complete the insulin assay.

## Radioimmunoassay

- 6 Radioimmunoassay kits are used to measure insulin levels. Several kits are available from MilliporeSigma. For the protocol please refer to the manufacturer's instruction.