

Jun 23, 2025

Stereo-seq for FF-OCT and FFPE

DOI

<https://dx.doi.org/10.17504/protocols.io.14egnyp3pv5d/v1>

Joanna Bi¹

¹Stanford University

Human BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu



Joanna Bi

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egnyp3pv5d/v1>

Protocol Citation: Joanna Bi 2025. Stereo-seq for FF-OCT and FFPE . protocols.io

<https://dx.doi.org/10.17504/protocols.io.14egnyp3pv5d/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 23, 2025

Last Modified: June 23, 2025

Protocol Integer ID: 220775

Keywords: seq for ff, seq protocol, ffpe sample, ffpe, stereo, seq

Abstract

Stereo-seq protocol for both FF-OCT and FFPE samples

FF-OCT

- 1 FF-OCT samples were sectioned into 10 μm sections using a cryostat (Leica, CM1950) and mounted onto 1 cm x 1 cm Stereo-seq T chips (Complete Genomics, 211ST114).
- 2 Mounted chips were processed through the Version 1.2 fresh-frozen Stereo-Seq workflow using kit components (Complete Genomics, 111KT114) and according to the manufacturer's protocol.
 - 2.1 Briefly, chips were dried at 37°C for 4 minutes, followed by methanol fixation for 30 minutes at -20°C, Qubit ssDNA staining (Invitrogen, Q10212), and FITC imaging (Zeiss, Axioscan Z1).
 - 2.2 Tissues were then permeabilized for 24 minutes, followed by reverse transcription for 3 hours at 42°C, tissue removal, and cDNA collection.
 - 2.3 cDNA were then purified (Vazyme, N411-01) and PCR amplified using kit components.
- 3 Library preparation was performed using a commercial kit (Complete Genomics, 111KL114) according to manufacturer's directions.
- 4 DNA Nanoballs (DNBs) were prepared using the OneStep DNB kit (MGI, 1000020563) and sequenced on a DNBSeq-T7 Sequencer (Complete Genomics, 900-000698-00) using the PE100 sequencing set (Complete Genomics 940-001889-00, with custom primer set 940-000037-00).
 - 4.1 Sequencing was done using 50 cycles for Read 1 (dark cycles 26-40), 100 cycles for Read 2, and 10 cycles of index read.
- 5 Primary analysis was performed using STOmics' Standard Analysis Workflow pipeline (SAW, Version 7), aligning to the human reference genome Hg38.

FFPE

- 6 FFPE samples were sectioned into 5 μm sections using a microtome (Dakewe, MT1) and mounted onto Stereo-seq N chips (STOmics, 210CN114-EA).
- 7 The chips were processed through the Stereo-Seq FFPE early-access workflow using the kit components (STOmics, 211KN114-EA) according to the manufacturer's protocol.



- 7.1 Briefly, mounted chips were dried at 42°C for 3 hours, 37°C for 16 hours, and 60°C for 1 hour.
- 7.2 Tissue was cleared in Sub-X (Leica, 3803670) twice for 20 minutes each, followed by rehydration with serial incubations in 100% ethanol, 96% ethanol, 90% ethanol, 80% ethanol, 70% ethanol, and nuclease-free water.
- 7.3 Decrosslinking was performed for 30 minutes at 85°C, followed by methanol fixation for 20 minutes at -20°C.
- 7.4 Staining, imaging, and permeabilization (30 minutes) were performed in the same way as for FF-OCT samples.
- 7.5 Following permeabilization, chips were processed through dimerization with kit reagents for 1 hour at room temperature.
- 8 Subsequent reverse transcription, tissue removal, cDNA collection, cDNA amplification, library preparation, and DNB preparation were conducted according to the manufacturer's protocol.
- 9 Sequencing was performed on a DNBSeg-T7 Sequencer using custom MDA primer 5' CTGCTGACGTACTGAGAGGC 3' and read 2 primer 5' GAGACGTTCTCGACTCAGCAGAGGG 3'.
- 9.1 The sequencing protocol included 91 cycles of Read 1 (dark cycles 26-85), 100 cycles of Read 2 (dark cycles 1-25), and 10 cycles of index read.
- 10 Primary analysis was performed using STOmics' Standard Analysis Workflow pipeline (SAW FFPE EA, with cell segmentation added from SAW V7 modified by expansion to 10 pixels).