



Feb 03, 2023

Stereologic Cell Population Estimates in Cleared Human Brain Tissue

DOI

dx.doi.org/10.17504/protocols.io.4r3l273dgg1y/v1

Bridget Wicinski¹, Patrick.hof¹

¹Icahn School of Medicine at Mount Sinai



Bridget Wicinski

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l273dgg1y/v1>

Protocol Citation: Bridget Wicinski, Patrick.hof 2023. Stereologic Cell Population Estimates in Cleared Human Brain Tissue. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l273dgg1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 28, 2023

Last Modified: February 03, 2023

Protocol Integer ID: 75991

Keywords: stereologic cell population estimate, step protocol for stereologic assessment, stereologic assessment, stereologic analysis, stereologic tool, mbf bioscience stereo investigator cleared tissue software, region of cortex, analysis for distinct neuronal population, numerical quantities of various morphologic parameter, cleared human brain tissue, neuronal density, human brain tissue, contouring of cortical layer, neuronal densities in each layer, human cerebral cortex, cleared postmortem specimen, systematic sampling, various morphologic parameter, surface area, distinct neuronal population, volume of reference, neuronal population, laminar surface area, contouring, object throughout the analysis, sampling, cortical layer, cortex, sampling design, large specimen, object under study, encompassing broca, sampling parameter, slice from large specimen, brain, defined region, specific cytoplasmic marker, counting frame size

Funders Acknowledgements:

Patrick Hof

Grant ID: U01 MH117023-01

Abstract

We provide a step-by-step protocol for stereologic assessment in the human cerebral cortex using cleared postmortem specimens with a laminar level of analysis for distinct neuronal populations identified by immunolabeling for specific cytoplasmic markers. Stereologic tools provide a set of simple rules and formulas to estimate numerical quantities of various morphologic parameters in tissue, with precision, accuracy and a design free of observer-induced bias. Parameters such a number, density, volume, surface area, or length can be easily estimated based on the systematic sampling of a region of interest (e.g., a layer of a cytoarchitecturally defined region of cortex, defined as a “volume of reference”), with an observer-independent random design, that gives each object under study (e.g., a neuronal population identified by a specific protein marker) the same probability to be sampled once and only once in its volume of reference, based on strict sampling criteria that are kept constant for a given object throughout the analysis. Stereologic analysis was performed on cleared slice from large specimens encompassing Broca’s area (Brodmann’s area 44/45), up to 4×4×2 cm, imaged at 3×3×3 μm pixel dimensions, and immunostained for NeuN, calretinin, and somatostatin, using the MBF Bioscience Stereo Investigator Cleared Tissue software (version 2022.1.1) with an Optical Fractionator design. Data were recorded as estimates of total number of neurons as well as neuronal densities in each layer. The protocol provides the detail of the sampling parameters including counting frame size, grid size, and disector height, definition and contouring of cortical layers and use of their boundaries to estimate laminar surface areas and volume. Full details on software operation and sampling design are explained.

Attachments



[Protocol_SI-CTE.pdf](#)

797KB

Troubleshooting

