

Jan 11, 2024

Version 1

Stereotaxic Injections and Implants V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5xdw4g1b/v1>

Jonathan Tang^{1,2,3}

¹Department of Neuroscience, Zuckerman Mind Brain Behavior Institute, Columbia University, New York, NY 10027, USA;

²Seattle Children's Research Institute, Center for Integrative Brain Research, Seattle, WA, USA;

³University of Washington School of Medicine, Department of Pediatrics, Seattle, WA, USA



Christine Weber-Schmidt

Allen Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5xdw4g1b/v1>

Protocol Citation: Jonathan Tang 2024. Stereotaxic Injections and Implants . protocols.io

<https://dx.doi.org/10.17504/protocols.io.8epv5xdw4g1b/v1> Version created by **Christine Weber-Schmidt**

Manuscript citation:

Tang, J.C.Y., Paixao, V., Carvalho, F. *et al.* Dynamic behaviour restructuring mediates dopamine-dependent credit assignment. *Nature* (2023). <https://doi.org/10.1038/s41586-023-06941-5>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 22, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 92645

Keywords: ASAPCRN, mouse, brain, AAV, stereotaxic injection, methods for stereotaxic injection, implants from tang et al, implant

Funders Acknowledgements:

Life Sciences Research Fellowship

NINDS K99/R00 Award

Grant ID: 1K99NS112575

NIH

Grant ID: 5U19NS104649

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-020551

Abstract

Methods for Stereotaxic injections and implants from Tang et al 2023.

Materials

Please see references for similar materials to what was used.

VTA AAV Viral Steriotaxic Injections

- 1 Using a standard mouse surgical suite (see references for example), 750nl of rAAV.EF1a.DIO.hChR2(H134R).eYFP or rAAV.EF1a.DIO.eYFP ($3\text{--}4 \times 10^{12}$ vg/ml, AAV5, University of North Carolina Vector Core; $1\text{--}2 \times 10^{13}$ vg/ml, AAV1, Addgene, 27056-AAV1 and 20298-AAV1) were injected into each hemisphere of the VTA of 3–4-month-old DAT-Cre mice.
 - 1.1 For viral injections, the coordinates are AP - 3.52 mm, ML - $\pm$ 0.35 mm, DV - 4.3 mm. Injections were made at 0.2 Hz pulses. Each pulse injects 4.6 nl volume.
 - 1.2 Injected needles were kept in place in the injection site for ~15 minutes before withdrawal.

Dual Optic Fiber Cannula Implant

- 2 For each mouse, a dual optic fiber cannula (200/240 μ m diameter, 6 mm length, 0.7 mm center-to-center FLT, 0.22 NA; Doric, DFC_200/240-0.22_6mm_DF0.7_FLT) was placed 200 μ m above the injection site and fixed to the skull.
 - 3 Next, a 4-position receptacle connector (Harwin Inc., M52-5000445) was fixed anteriorly to the dual optic fiber cannula, with its posterior edge set at -0.6 mm.
 - 3.1 Skull implants are then fixed with dental cement.
 - 3.2 A 4-position connector (Harwin Inc., M52-040023V0445) with pins removed from one end was used to cap the receptacle connector.

Photometry Experiments

- 4 The conditions used for VTA injections and implants were as above, and 3–5-month-old DAT-Cre males were used.
 - 4.1 Additionally, 1 μ l and 500 nl of AAV9-hSyn-GRAB-rDA1m (2×10^{13} vg/ml; RRID:Addgene_140556) were injected into the dorsal striatum (AP 0.5 mm, ML +2.1 (right), DV 2.3 (from brain surface)) and ventral striatum (AP 1.15mm, ML +1.65 (right), DV 4.2 (from Bregma)), respectively.



- 5 For photometry fiber implants, mono fiberoptic cannula were used (400/430 μ m diameter, 4 mm length (dorsal striatum) and 6 mm length(ventral striatum), 0.37 NA, 1.25 diameter ferrule, flat; Doric, MFC_400/430-0.37_6mm_MF1.25_FLT (ventral striatum) and MFC_400/430-0.37_4mm_MF1.25_FLT (dorsal striatum)). Implants were inserted at a 22 degrees angle.

Protocol references

Paxinos, George & Franklin, K. B. J. The mouse brain in stereotaxic coordinates / George 1431 Paxinos, Keith B.J. Franklin. (2001).

Allen Institute for Brain Science 2023. Stereotaxic Injection by Nanoject

Protocol. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bp2l6nr7kgge/v4>Version created by **Allen Institute**