

Nov 21, 2020

stLFR library construction

 [GigaByte](#)

 In 2 collections

DOI

dx.doi.org/10.17504/protocols.io.bpwxmpfn

Chunhua Li¹, Xianwei Yang¹, Libin Shao¹, Rui Zhang¹, QunLiu¹, Mengqi Zhang¹, Shanshan Liu¹, Shanshan Pan¹, Weizhen Xue¹, Congyan Wang¹, Chunyan Mao¹, He Zhang², Guangyi Fan¹

¹BGI-Qingdao, BGI-Shenzhen, Qingdao 266555, China.;

²BGI-Qingdao, BGI-Shenzhen, Qingdao 266555, China, Department of Biology, Hong Kong Baptist University, Hong Kong, China.

GigaScience Press

BGI



lichunhua



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bpwxmpfn

External link: <https://doi.org/10.46471/gigabyte.32>

Protocol Citation: Chunhua Li, Xianwei Yang, Libin Shao, Rui Zhang, QunLiu, Mengqi Zhang, Shanshan Liu, Shanshan Pan, Weizhen Xue, Congyan Wang, Chunyan Mao, He Zhang, Guangyi Fan 2020. stLFR library construction. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bpwxmpfn>

Manuscript citation:

Chunhua Li, Xianwei Yang, Libin Shao, Rui Zhang, Qun Liu, Mengqi Zhang, Shanshan Liu, Shanshan Pan, Weizhen Xue, Congyan Wang, Chunyan Mao, He Zhang, Guangyi Fan, Bicolor angelfish (*Centropyge bicolor*) provides the first chromosome-level genome of the Pomacanthidae family, *Gigabyte*, 2021 <https://doi.org/10.46471/gigabyte.32>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: November 20, 2020

Last Modified: November 21, 2020











Protocol Integer ID: 44727

Keywords: Bicolor Angelfish, genome, stLFR



Abstract

This protocol is used to clarify the process of stLFR library preparation for Bicolor Angelfish (*Centropyge bicolor*) genome.




- 1 Sample preparation
Remove the long-fragment gDNA from  4 °C and store on ice.
- 2 Transposon Insertion
 - 2.1 Transfer  10 ng long-fragment gDNA gently into a 0.2ml PCR tube. Without mixing, add Nuclease Free Water to a total volume of  36.8 µL. Collect and dispense long fragment DNA slowly (the process should take >  00:00:10) each time when pipetting. 10s
 - 2.2 Prepare the stLFR_SamBarTIE and TE Buffer. Dilute the stLFR_SamBarTIE 16-fold with TE Buffer. Add  6 µL TE Buffer into a new 0.2 mL PCR tube and then transfer  2 µL stLFR_SamBarTIE to the tube. Vortex intermittently for 4 times (2s each) to mix. Label it as "4× dilution stLFR_SamBarTIE".
 - 2.3 Add  18 µL TE Buffer into a new 0.2 mL PCR tube and transfer  6 µL of the 4× dilution stLFR_SamBarTIE into the tube. Vortex intermittently for 4 times (2s each) to mix the tube. Label it as "stLFR_SamBarTIE (Working Mix)." stLFR_SamBarTIE (Working Mix) can be used for 12 reactions.
 - 2.4 According to the requirements of the sample sequencing pooling strategy, mix the stLFR_SamBarTIE Working Mix.
 - 2.5 Prepare the Transposon Insertion Reaction Mix on ice as shown. Vortex intermittently for 4 times (2s each) to mix.
 - 2.6 Transfer  13.2 µL Transposon Insertion Reaction Mix to the DNA sample from step 2.1 (total volume is  50 µL). Mix by very gently pipetting 10 times using wide-bore tip. Briefly centrifuge the tube. Transfer the tube to the thermocycler and start the Transposon Insertion Reaction Program.
 - 2.7 Store the tubes on ice after the transposon insertion step has finished.





2.8 Pipette  42.5 μL TE buffer to a new 0.2 mL tube then transfer  7.5 μL transposon-inserted product (from step 2.6) to this new 0.2 mL tube. Mix by inverting the tube very gently and collect liquid to the bottom of the tube by briefly centrifuging (1 second) on a microcentrifuge. Label this new tube as the sample tube.


3 Capture






3.1 Vortex Capture Beads to mix thoroughly before use. Pipette  30 μL Capture Beads per sample to the tube.

3.2 Place the tube on a magnetic separation rack. Once the liquid is clear, carefully collect and discard the supernatant.

3.3 Pipette  50 μL Wash Buffer I per sample into the 0.2 mL PCR tube or 1.5 mL EP tube. Ensure that Wash Buffer I can cover all of the Capture Beads.

3.4 Rotate the tube 180 degrees within the rack such that the beads are forced to pass through the Wash Buffer I. Repeat the tube rotation once. Carefully remove and discard the supernatant once the liquid in the tube is clear. Pipette  50 μL Capture Buffer per sample to resuspend the Capture Beads.

3.5 Transfer  50 μL Capture Beads from step 3.4 to the sample tube from step 2.7 and mix thoroughly by inverting gently at least 10 times.



3.6 Centrifuge the product from step 3.5 for 1 second and place on the rotator in the incubator. Immediately start rotating the sample. In this experiment, After the first  00:10:00 of incubation at  60 $^{\circ}\text{C}$, switch the temperature of the incubator to  45 $^{\circ}\text{C}$. Open the door of the incubator to accelerate cooling, then close the door of the incubator and start the  00:50:00 countdown once the temperature drops to  48 $^{\circ}\text{C}$.




1h


4 Ligation Reaction 1

4.1 Centrifuge the product from step 3.6 and allow the product to cool to room temperature.






4.2 Ensure the product has cooled to room temperature, then transfer  30 μL Ligation Reaction 1 Mix to the  100 μL sample. Mix by inverting the tubes gently at least 10 times then briefly centrifuge (1s). Place the sample tube on the rotator and turn it on.

4.3 Perform Ligation Reaction 1 with the incubation in  Room temperature (20°C to 25°C)  10m
 00:10:00 . After incubation, centrifuge the sample and place it on the magnetic separation rack. Carefully remove and discard the supernatant once the liquid is clear.



4.4 Pipette  180 μL of Wash Buffer II into the sample tube. Rotate the tube 180 degrees while on the magnetic separation rack to let the beads move through the Wash Buffer II. Repeat the tube rotation once. Carefully remove and discard the supernatant once the liquid is clear. Keep the Capture Beads in the Wash Buffer II until the Digestion Reaction Mix 1 is ready.

5 Digestion Reaction 1

5.1 Transfer  100 μL Digestion Reaction 1 Mix to the sample tube from step 4.4. Mix by inverting the tube gently at least 10 times followed by an instantaneous centrifugation (1s).

5.2 Place the sample tube on the rotator and turn it on. Perform the Digestion Reaction 1 incubation in  37 °C  00:10:00 . Remove the sample tube from the incubator at the end of the reaction. Immediately add TIS Buffer.

6 Termination Reaction

6.1 Remove the sample from  37 °C , centrifuge briefly, and store at room temperature. Immediately add  11 μL TIS Buffer to each sample from step 5.2.

6.2 Ensure the sample tube is sealed tightly. Mix the sample tube by vortexing at medium speed for 3 to 5 seconds to make sure the beads are fully resuspended. Centrifuge the



tube for 00:00:01 and place on the rotator. Start the rotator and perform the incubation in Room temperature (20°C to 25°C) 10 minutes.

- 6.3 After incubation, centrifuge the sample and place on the magnetic separation rack. Carefully remove and discard the supernatant once the liquid is clear. Keep tube on the magnetic separation rack and pipette 150 μL Wash Buffer II into the sample tube. Mix the beads by vortexing for 00:00:05 . Centrifuge briefly and place the tube back onto the magnetic rack for 00:02:00 . Carefully remove and discard the supernatant once the liquid becomes clear. Repeat this step twice.

2m 5s

7 Pre-Ligation Reaction 2

- 7.1 Pipette 24 μL Pre-Ligation 2 Reaction mix into the sample tube from step 6.3. Mix the sample tube by vortexing for 3-5 seconds to ensure the beads are fully resuspended. Centrifuge the tube for 00:00:01 and place on the rotator stored in the 37 °C incubator. Perform the incubation in 37 °C 00:10:00 .

10m 1s

- 7.2 When the Pre-Ligation 2 Reaction is complete, remove the product from incubator immediately and centrifuge briefly. Keep it at Room temperature and proceed to next step.

8 Ligation Reaction 2

- 8.1 Pipette all 76 μL Ligation Reaction 2 mix into the sample tube for a total volume of 100 μL . Mix the sample tube by vortexing for 00:00:10 . Centrifuge the tube for 00:00:01 to make sure beads are fully resuspended.

11s


- 8.2 Place the tube on rotator and perform the incubation in Room temperature (20°C to 25°C) 02:00:00 .

2h



- 8.3 After incubation, centrifuge the sample and pipette 80 μL Wash Buffer II into the sample tube. Place the sample tube on the magnetic separation rack for 00:02:00 .

2m







Carefully remove and discard the supernatant when the liquid becomes clear.

- 8.4 Keep the sample tube on the magnetic separation rack and pipette  180 μL Wash Buffer II into the sample tube. Rotate the tube 180 degrees within the magnetic separation rack to allow the beads to move through Wash Buffer II. Repeat once. Carefully remove and discard the supernatant when the solution becomes clear. Make sure Wash Buffer II is completely removed.











9 PCR

- 9.1 Pipette  150 μL of the PCR mix into sample tube. Use the pipet to mix the beads until fully resuspend. Transfer  75 μL of the sample to a different 0.2 mL tube.
- 9.2 Place all samples on the thermocycler. Make sure the beads are fully resuspended before starting.
- 9.3 Centrifuge the sample and place on the magnetic separation rack. Transfer all the supernatant of the two PCR tube from the same sample into a new 1.5 mL EP tube and mix together.
- 9.4 After confirming complete recovery of the supernatant, discard the original PCR tube.

10 PCR Product Purification

- 10.1 Remove the DNA Clean beads from  4 $^{\circ}\text{C}$ and equilibrate to  Room temperature for at least  00:30:00 before use. Vortex at full speed for  00:00:10 to ensure the beads are completely resuspended. 30m 10s
- 10.2 Measure the volume of the PCR product from step 9.3. Add 0.7-fold DNA Clean beads to the PCR product. Vortex the tube to mix the beads with the sample.
- 10.3 Incubate the sample at  Room temperature for  00:10:00 . 10m



- 10.4  75 μL Centrifuge the tube and place on the magnetic separation rack. Wait for  00:02:00 or until the solution is clear. Discard the supernatant. 2m
- 10.5 Keep the sample tube on the magnetic separation rack and transfer  500 μL 80% (v/v) ethanol into the tube. Let stand for  00:00:30 . Carefully remove and discard the supernatant. Repeat this step once more. 30s
- 10.6 Keep the sample tube on the magnetic separation rack, open the cap of tube and air-dry the beads for 3 to 5 minutes until no wetness is observed (the surface of the beads will dim). Do not over dry the beads as this will significantly decrease the elution efficiency (cracks can be observed on pellet).
- 10.7 Remove the sample tube from the magnet and add  33 μL TE Buffer for elution. Vortex for  00:00:03 to resuspend the beads, then briefly centrifuge. 3s
- 10.8 Incubate the sample at  Room temperature for  00:05:00 . 5m
- 10.9 Place the sample tube on the magnetic separation rack and wait for  00:02:00 or until the supernatant is clear. Transfer  31 μL of supernatant from the sample tube to a new 1.5 mL EP tube. Do not disturb or pipette the beads. 2m