

Dec 22, 2025

STN optotagging

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzen3rvx1/v1>

Asa Mackenzie¹

¹Lund University



Åsa Wallén-Mackenzie

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzen3rvx1/v1>

Protocol Citation: Asa Mackenzie 2025. STN optotagging. protocols.io
<https://dx.doi.org/10.17504/protocols.io.rm7vzen3rvx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 22, 2025



Last Modified: April 21, 2026

Protocol Integer ID: 235674

Keywords: ASAPCRN, optotagging stn

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020600

Abstract

STN optotagging

STN optotagging

- 1 A custom-made optrode was made with an optic fiber (100 μ m diameter, Thorlabs) connected to a laser (MBL-III-473 nm-100 mW laser, CNI Lasers, Changchun, China) mounted and glued on the recording glass micropipette which was filled with 2% pontamine sky blue in 0.5 mM sodium acetate (tip diameter 1–2 mm, resistance 10–15 MU).
- 2 The distance between the end of the optic fiber and the tip of the recording pipette varied between 650 nm and 950 nm.
- 3 Extracellular action potentials were recorded and amplified with an Axoclamp-2B and filtered (300 Hz/0.5 kHz).
- 4 Single extracellular spikes were collected online (CED 1401, SPIKE2; Cambridge Electronic Design).
- 5 The laser power was measured before starting each experiment using a power meter (PM100D; Thorlabs). The baseline was recorded for 100 s for each neuron before starting any light stimulation protocols which were set and triggered with Spike2 software.
- 6 Light protocol consisted in a peristimulus time histogram (PSTH, 0.5 Hz, 5 ms pulse duration, 5–8 mW) for at least 100 s followed, after returned to baseline, by a “behavioral” protocol corresponding to the parameters used for behavioral experiments (20 Hz, 5 ms pulse duration, 5–8 mW).