

Jun 13, 2025

Stony Coral Single Cell Isolation

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzqx72vx1/v1>

Maria Ingersoll¹, Sarah Davies¹

¹Boston University



Maria Ingersoll

Boston University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzqx72vx1/v1>

Protocol Citation: Maria Ingersoll, Sarah Davies 2025. Stony Coral Single Cell Isolation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzqx72vx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 05, 2024



Last Modified: June 13, 2025

Protocol Integer ID: 114326

Keywords: Coral, Single Cell Isolation, Cnidaria, symbiotic stony coral *Oculina*, cell rna, live cell, live cells for downstream use, cell, rna

Funders Acknowledgements:

National Science Foundation

Grant ID: IOS-1937650

National Science Foundation

Grant ID: Graduate Research Fellowship

National Science Foundation Research Traineeship

Grant ID: DGE 1735087

Boston University Marine Program

Grant ID: Warren McLeod Marine Fellowship

Abstract

This protocol was optimized in the Davies lab at Boston University for use with the facultatively symbiotic stony coral *Oculina arbuscula*. The protocol is designed to isolate live cells for downstream use in 10X single-cell RNA-sequencing.



Materials

Reagents

Sterile DI Water

0.2 micron filter-sterilized Calcium-Free Seawater (CfFSW)

Fetal Bovine Serum (FBS)

10X PBS

1X PBS

1M HEPES

Protector RNase Inhibitor (Sigma Aldrich 3335399001)

0.4% Trypan Blue Solution (filtered through 0.2 micron filter to remove aggregates)

RNase Zap or equivalent

Materials

Heat-bath or hot-plate capable of 56°C

6-well plate

Sterile 200µl pipette tips

Forceps/tweezers

Sterile 1.5 ml tubes

Sterile 15 ml conical tubes

Sterile 70 micron cell strainers (e.g., Fisher Scientific 08-771-2)

Sterile 40 micron cell strainers (e.g., Fisher Scientific 08-771-1)

Microcentrifuge capable of 4°C

Hemocytometer

Light Microscope

Safety warnings

- ❗ 1. As much as possible, all steps should be performed on ice
- 2. Do not vortex cell suspension at any point during the protocol
- 3. Limit bubbles in cell suspension
- 4. Cells should be processed promptly following dissociation to maintain viability (viability decreases within an hour following isolation)



- 1 Heat inactivate FBS for 30 minutes at 56°C
- 2 Sterilize tools (i.e., forceps, pipettes, etc.) and bench space with RNase Zap
- 3 Prepare cell media¹ and keep on ice
- 3.1 To make 20 ml of cell media:
 - 6.6 ml 10X PBS (final concentration of 3.3X)
 - 0.4 ml heat-inactivated FBS (final concentration of 2% by volume)
 - 0.4 ml 1M HEPES (final concentration of 20mM)
 - Fill to 20 ml with sterile DI H₂O
- 4 Add 10 mls of CafFSW² to two wells of a 6-well plate, and add 10 mls of cell media to a third well. Keep plate on ice.
- 5 Obtain live coral frag (this protocol has been optimized for frags of symbiotic or aposymbiotic branching coral ~1.5 cm in length and ~1 cm in width) and place in first well of CafFSW
- 6 Incubate for 1 minute (on ice)
- 7 With sterilized forceps, move frag to second well of CafFSW, and incubate for 1 minute (on ice)
- 8 With sterilized forceps, move frag to well containing cell media. Using forceps and sterile 200 µl pipette tips, carefully scrape all tissue from the coral skeleton and from within the polyps³. This step should last no more than 10 minutes. If necessary, this step can be performed on the bench, not on ice.

TIP: Dig the pipette tip or forcep tip into a polyp and make small circles to extract cells from within the polyp



- 9 Filter full cell suspension through one 70 micron cell strainer into a 15 ml conical tube.


NOTE: for filtration steps, pipette 1 ml aliquots through strainer by holding tip against mesh, depressing the plunger slowly and assuring all liquid passes through the strainer. Keep plunger depressed while removing tip from mesh. Cell suspension may be mucus-y at this point, so some volume (and cells) will inevitably be lost.



- 10 Filter 70-micron-filtered suspension through one 40 micron cell strainer into a new 15 ml conical tube
- 11 Filter suspension through a second 40 micron cell strainer into a third 15 ml conical tube
- 12 Centrifuge 1 ml aliquot(s) of filtered cell suspension at 300xg for 10 minutes at 4°C



- 13 Discard supernatant and resuspend pellet in 100 μ l of cell media with 0.2 U/L Protector RNase Inhibitor
 - 13.1 Right before use, add 0.5 μ l of stock RNase Inhibitor (at 40 U/ μ l) per 100 μ l cell media
- 14 Check cell viability and perform cell counts (using hemocytometer)
 - 14.1 Combine 2.5 μ l of 0.4% Trypan Blue with 12.5 μ l of 1X PBS and 10 μ l of final cell resuspension (for final concentration of 0.04% Trypan)
 - 14.2 For 10X, cell viability should be over 90% (dead cells are completely permeated with blue dye)
 - 14.3 Symbiotic cells (algal cells) should be counted in cell count and viability calculations



Example of cell suspension from symbiotic *Oculina arbuscula* with 0.04% Trypan Blue visualized at 40X magnification

- 15 Dilute cells to desired concentration in cell media with RNase Inhibitor (see step 13.1 for concentration)
- 15.1 For 10X at the Boston University Single Cell Core, cells were provided at concentrations of 700-1,200 cells/ μ l



Protocol references

1. Rosental B, Kozhekbaeva Z, Fernhoff N, Tsai JM, Traylor-Knowles N. Coral cell separation and isolation by fluorescence-activated cell sorting (FACS). *BMC Cell Biology*. 2017 Aug 29;18(1):30.
2. Roger LM, Reich HG, Lawrence E, Li S, Vizgaudis W, Brenner N, et al. Applying model approaches in non-model systems: A review and case study on coral cell culture. *PLoS ONE* 2021 Apr 8;16(4):e0248953.
3. Levy S, Elek A, Grau-Bové X, Menéndez-Bravo S, Iglesias M, Tanay A, et al. A stony coral cell atlas illuminates the molecular and cellular basis of coral symbiosis, calcification, and immunity. *Cell*. 2021 May 27;184(11):2973–2987.e18.