

May 07, 2020

## Storage and Processing of Tissue for bulk RNA Isolation

DOI

[dx.doi.org/10.17504/protocols.io.bf52jq8e](https://dx.doi.org/10.17504/protocols.io.bf52jq8e)

Aaron Horning<sup>1</sup>

<sup>1</sup>Stanford University

NCIHTAN

Human BioMolecular Atl...



Aaron Horning

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bf52jq8e>

**Protocol Citation:** Aaron Horning 2020. Storage and Processing of Tissue for bulk RNA Isolation. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bf52jq8e>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** May 07, 2020

**Last Modified:** May 07, 2020

**Protocol Integer ID:** 36762

**Keywords:** bulk rna isolation, rna, processing of tissue, tissue, storage

## Troubleshooting



## Storing Tissue for RNA Isolation

- 1 After resection, tissue is immediately flash frozen in 50 ml Falcon tube or 2 ml cryovial tubes by placing tissue into tube then simply dropping filled tube into liquid nitrogen.
- 2 The tubes are stored on dry ice or -80C until they can be stored in a liquid nitrogen tank

## Bulk RNA Isolation

- 3 The large tissue pieces are smashed in a freezing cold mortar and pestle with liquid nitrogen. The desired amount of tissue (30mg - 60mg) is then pulverized and swept into a 1.5 ml centrifuge tube sitting on dry ice to prevent thawing of the tissue.

### Note

Its critically important for RNA integrity to not let the tissue thaw. Thawing will activate the endogenous RNases within the tissues thus depleting the chemical we wish to assay.

- 4 Follow the Qiagen All Prep Kit protocol found here: <https://www.qiagen.com/us/resources/resourcedetail?id=bbd50261-3b80-4657-ad58-6a5a97b88821&lang=en>. Complete the RNA isolation steps.
- 5 Briefly, be sure to lyse and homogenize the tissue with RLT buffer with beta-mercaptoethanol with an 18 gauge needle until the needle runs smoothly. Additionally, if once the lysed tissue in RLT appears to be at all yellow and not clear, dilute the lysed tissue in RLT in more RLT buffer. If the RLT is too dense, it will not run through the columns and you'll loose material.
- 6 Elute in water and freeze at -80C.