



Apr 28, 2023

Strategies for optimizing the isolation and expansion of sensitive patient-derived duodenoids, ileoids and colonoids .pdf

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx9d6og8j/v1>

Katlynn Bugda Gwilt¹, Jay Thiagarajah¹

¹Boston Children's Hospital



Katlynn Bugda Gwilt

Boston Children's Hospital, Harvard Medical School

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx9d6og8j/v1>

Protocol Citation: Katlynn Bugda Gwilt, Jay Thiagarajah 2023. Strategies for optimizing the isolation and expansion of sensitive patient-derived duodenoids, ileoids and colonoids .pdf. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kxygx9d6og8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working as described

Created: April 28, 2023

Last Modified: April 28, 2023

Protocol Integer ID: 81156

Keywords: colonoids from human sample, pediatric patients with genetic intestinal epithelial disorder, genetic intestinal epithelial disorder, colonoids from these patient, enteroids from patient, enteroid growth, colonoid, derived duodenoid, ileoid, intestinal disease, culturing enteroid, onset inflammatory bowel disease, congenital diarrhea, primary duodenoid, enteroid, duodenoid, patients without any intestinal disease, enteropathy, ileum, initial culture time after cell isolation, colon, cellular phenotype, novel disease pathophysiology, endoscopic biopsies from patient, endoscopic biopsy, defects in proliferation, reproducible phenotype, cell isolation

Funders Acknowledgements:

NIDDK

Grant ID: RC2 DK118640

NIDDK

Grant ID: F32 DK130267

Abstract

These protocols are an optimization of previous protocols and the original work described in Sato et al., 2009 which provided a foundation for the development of primary duodenoids, ileoids and colonoids from human samples. Here, we optimized methods specifically for use with samples from pediatric patients with genetic intestinal epithelial disorders. These protocols apply specifically to enteroids generated from endoscopic biopsies from patients with Congenital Diarrheas and Enteropathies (CoDE) or Very-Early-Onset Inflammatory Bowel Disease (VEO-IBD) or age-matched patients without any intestinal disease and include samples from the duodenum (duodenoids), ileum (ileoids), and colon (colonoids). The duodenoids, ileoids and colonoids from these patients exhibit a range of cellular phenotypes that potentially impact enteroid growth and maintenance including increased apoptosis, defects in proliferation, polarity and vesicular trafficking. These characteristics make it potentially challenging for long-term culture and expansion, limiting the availability of cells for functional experiments. Using a modified culture media, an expanded initial culture time after cell isolation from patient biopsies, and gentle passaging techniques, we were able to successfully culture and expand several patient lines over multiple passages and maintain the cultures for >5 years. The media conditions and protocols described here allow for reproducible phenotypes as well as scaling for larger functional studies on patient lines. These protocols also provide a useful starting point for further optimization for generating and culturing enteroids from patients with novel disease pathophysiology.

Attachments



Strategies for optim...

237KB

Materials

	Key resources table	
	Reagent	Supplier
	TrypLE™ Express Enzyme (1X), phenol red	Thermo Fisher Scientific
	Advanced DMEM/F12	Gibco
	L-WRN Conditioned Media	ATCC (Cite Clevers)
	Wnt Conditioned Media	ATCC
	R-spondin-1 conditioned Media	n/a
	Noggin Conditioned Media	n/a
	Recombinant Noggin	Peperotech
	Glutamax	Gibco
	HEPES	Gibco
	Primocin	Invivogen
	Normocin	Invivogen
	B27	Gibco
	N2	Gibco
	Nicotinamide	Sigma-Aldrich
	N-acetyl-cysteine	Sigma-Aldrich
	A83-01	Sigma-Aldrich
	SB202190	Sigma-Aldrich
	EGF	Peprtech
	Gastrin	Sigma-Aldrich
	Y27632	Sigma-Aldrich
	Prostaglandin E2	Sigma-Aldrich
	CHIR99021	Sigma-Aldrich
	Corning® Matrigel® Growth Factor Reduced	Corning



	(GFR) Basement Membrane Matrix, Phenol Red-free, LDEV-free, 10 mL		
	Cultrex	R&D Technologies	BME001-10
	Cultrex Organoid Harvesting Solution, 100 mL	R&D Technologies	3700-100-01
	Corning® Cell Recovery Solution, 100 mL	Corning	Cat#354253
	Collagen IV from human placenta	Sigma-Aldrich	Cat#C5533
	Collagenase, Type 1	Stem Cell	Catalog # 07416
	6.5 mm Transwell with 0.4 µm pore polyester membrane insert, TC-treated, sterile, 48/cs	Corning	Cat#CLS3470-48EA
	24 well plate	Costar	Cat#3526

Before start

Dissolve all reagents to the appropriate concentration according to manufacturer's protocols. A * indicates that there were no modifications to the reagent.



Protocol references

Sato, T., Vries, R. G., Snippert, H. J., van de Wetering, M., Barker, N., Stange, D. E., van Es, J. H., Abo, A., Kujala, P., Peters, P. J., & Clevers, H. (2009). Single Lgr5 stem cells build crypt-villus structures in vitro without a mesenchymal niche. *Nature*, 459(7244), 262–265. <https://doi.org/10.1038/nature07935>