

Oct 16, 2023

Strep pull down assay

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4jmpolo5/v1>

Minghao Chen¹, Xuefeng Ren¹

¹Department of Molecular and Cell Biology, University of California, Berkeley, Berkeley, CA 94720, USA



Minghao Chen

University of California, Berkeley

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4jmpolo5/v1>

Protocol Citation: Minghao Chen, Xuefeng Ren 2023. Strep pull down assay. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.3byl4jmpolo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 05, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 82917

Keywords: ASAPCRN, assay strep, gst, mutant, assay

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) initiative

Grant ID: ASAP-000350

Abstract

Strep pull down assay for FIP200(NTD)-TSF WT or mutants were co-transfected with GST-ATG13(363-517) and/or ULK1(MIT)-MBP



Expression

2d 0h 40m

1 Transfect 10 ml HEK GNTi cells at concentration of 2×10^6 cells/ml

2 Dilute PEI with Warm Hybridoma-SFM(1X)

3 In a separate tube, dilute DNA with Hybridoma-SFM(1X)

4 Add PEI to DNA dilution. Incubate mixture for 00:30:00 at 37 °C

30m

5 Add mixture to cells. Let cells grow for 48:00:00

2d

6 Harvest Cells 500 rpm , 4 °C , 00:10:00

10m

7 Wash pellet with cold PBS. Store pellet at -80 °C until purification.

Pull down

2h 50m

8 Homogenize the pellets in 0.5 ml of lysis buffer/protease inhibitors/1% TritonX-100, and clarified after 40000 x g , 00:15:00 .

15m

9 Incubate the lysate with 30 µl Strep-Tactin Sepharose resin (IBA-Lifesciences) at 4 °C for 03:00:00

3h

10 Wash the beads four times, and eluted in 50 µl lysis buffer/4 mM desthiobiotin.



- 11 Mix 18 μ l eluent with lithium dodecylsulfate (LDS)/BME buffer, heated at 60°C for  00:05:00 and subjected to SDS/PAGE gel. The gel was then stained with Coomassie brilliant blue G250.

5m