



Aug 09, 2024

Striatal oxidative stress quantification (OxiSelect™ In Vitro ROS/RNS Assay Kit (Green Fluorescence))

DOI

dx.doi.org/10.17504/protocols.io.4r3l2qbqjl1y/v1

Livia Hecke Morais¹

¹California Institute of Technology



Livia Hecke Morais

California Institute of Technology

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2qbqjl1y/v1>

Protocol Citation: Livia Hecke Morais 2024. Striatal oxidative stress quantification (OxiSelect™ In Vitro ROS/RNS Assay Kit (Green Fluorescence)). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l2qbqjl1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: August 09, 2024

Last Modified: August 09, 2024

Protocol Integer ID: 105019

Keywords: striatal oxidative stress quantification, vitro ro, oxiselect, rns assay kit, green fluorescence, assay

Funders Acknowledgements:

ASAP

Abstract

Striatal oxidative stress quantification (OxiSelect™ In Vitro ROS/RNS Assay Kit (Green Fluorescence))

Attachments



STA-347-in-vitro-ros...

321KB

Troubleshooting

Oxidative stress quantification(OxiSelect™ In Vitro ROS/RNS Assay Kit (Green Fluorescence))

- 1 H_2O_2 is quantified using OxiSelect™ In Vitro ROS/RNS Assay Kit (Green Fluorescence) from CellBio Labs (Cat # STA-347-5) according to the manufacturer's instructions.

Sample preparation

- 2 Following dissection, whole striatal tissue is homogenized in 100 μL PBS using a tissue homogenizer.
- 3 Samples are diluted (1:25) in PBS. Different dilutions should be tested for different tissues and experimental conditions. Prepare enough for 110 microliters of diluted samples (50 μL x duplicates + spare volume).

Standards Preparation (please refer to manufacturer protocol attached)

- 4 Due to light-induced auto-oxidation, the stock DCFH solutions must be protected from light.
- 5 To prepare the Hydrogen Peroxide standards, first perform a 1:4400 dilution of the stock Hydrogen Peroxide in deionized water. Use only enough for immediate applications (e.g., Add 5 μL of Hydrogen Peroxide to 22 mL deionized water). This solution has a concentration of 2 mM.
- 6 Use the 2 mM H_2O_2 solution to prepare standards in the concentration range of 0 μM – 20 μM by further diluting in PBS (see Table 3 in attached manual). H_2O_2 diluted solutions and standards should be prepared fresh. Use the table below as a reference guide only. The volumes and concentrations of the standard may be adjusted by the user.

Prepare other reagents (please refer to manufacturer protocol attached)

- 7 1X Stabilization Solution: Dilute the 10X Stabilization Solution 1:10 by adding 1.5 mL of solution to 13.5 mL of deionized water. Stir or vortex to homogeneity. Store the solution at 4°C.
- 8 1X Catalyst: Prior to use, dilute the 250X Catalyst 1:250 in PBS. Vortex thoroughly. Prepare only enough for immediate applications (e.g., add 10 μL of Catalyst to 2.49 mL PBS for 50 wells).
- 9 DCFH Solution: Prepare only enough DCFH Solution for immediate applications in an amber tube or aluminum foil covered tube. Prepare DCFH Solution by diluting the stock



solution of DCF-DiOxyQ 1:5 with Priming Reagent (e.g., for 50 assays, add 25 μ L DCF-DiOxyQ to 100 μ L Priming Reagent). Vortex to homogeneity. Incubate the solution for 30 minutes at room temperature. Next, dilute the reaction 1:40 with 1X Stabilization Solution (e.g., for 50 assays, add 125 μ L DCF-DiOxyQ/ Priming Reagent reaction to 4.875 mL of Stabilization Solution). Vortex to homogeneity. Protect the solution from light. This solution is now stable in the DCFH form and ready to use. The solution may be stored at -20°C for up to one week when protected from light.

Assay (please refer to manufacturer protocol attached)

- 10 Prepare and mix all reagents thoroughly before use. Each sample, including unknown(s) and standard(s), should be assayed in duplicate or triplicate.
- 11 Add 50 μ L of unknown sample or hydrogen peroxide standard to wells of a 96-well plate suitable for fluorescence measurement
- 12 Add 50 μ L of Catalyst to each well. Mix well and incubate 5 minutes at room temperature.
- 13 Add 100 μ L of DCFH solution to each well. Cover the plate reaction wells to protect them from light and incubate at room temperature for 15-45 minutes.
- 14 Read the fluorescence with a fluorescence plate reader at 480 nm excitation / 530 nm emission.