

Nov 10, 2021

## Sub-ARTIC Illumina SARS-CoV-2 Spike sequencing protocol (LoCost) V3.2

DOI

<https://dx.doi.org/10.17504/protocols.io.btpjnmkn>

Gavin Horsburgh<sup>1</sup>, Kathryn Maher<sup>1</sup>, Steve Paterson<sup>2</sup>, Terry Burke<sup>1</sup>, Paul J Parsons<sup>3</sup>

<sup>1</sup>NERC Environmental Omics Facility, Ecology & Evolutionary Biology, School of Biosciences, University of Sheffield, Sheffield, UK;

<sup>2</sup>NERC Environmental Omics Facility, Department of Evolution, Ecology and Behaviour, Institute of Infection, Veterinary and Ecological Sciences, University of Liverpool, Liverpool, UK;

<sup>3</sup>University of Sheffield



Gavin J Horsburgh

University of Sheffield

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.btpjnmkn>

**Protocol Citation:** Gavin Horsburgh, Kathryn Maher, Steve Paterson, Terry Burke, Paul J Parsons 2021. Sub-ARTIC Illumina SARS-CoV-2 Spike sequencing protocol (LoCost) V3.2. [protocols.io https://dx.doi.org/10.17504/protocols.io.btpjnmkn](https://dx.doi.org/10.17504/protocols.io.btpjnmkn)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 26, 2021

**Last Modified:** November 10, 2021

**Protocol Integer ID:** 48587

**Keywords:** SARS-CoV-2, COVID, variant, sequencing, Illumina, ARTIC , SubARTIC, amplicon, sequencing v3, sequencing protocol, clinical rna sample, spike region of sar, method from the neof liverpool illumina artic protocol, neof liverpool illumina artic protocol, wastewater sample, short amplicon, using short amplicon, spike region

## Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Abstract

This protocol describes a procedure for sequencing the Spike region of SARS-CoV-2 using short amplicons (146-208bp). The method has proved to be successful with both clinical RNA samples and degraded wastewater samples. The primers are unique to this method. The library prep procedure has been heavily adapted from the ncov-2019 sequencing v3 (ARTIC) protocol by Josh Quick (<https://www.protocols.io/view/ncov-2019-sequencing-protocol-v3-locost-bh42j8ye>) and the "low cost" method from the NEOF Liverpool Illumina ARTIC protocol.



## Materials

Primers are listed in the attached  SubArtic primers v3-2 091121.csv

LunaScript<sup>®</sup> RT SuperMix Kit, New England BioLabs, Cat# E3010L

Q5<sup>®</sup> Hot Start High-Fidelity 2X Master Mix, New England BioLabs, Cat# M0494L

NEBNext<sup>®</sup> Ultra<sup>™</sup> II DNA Library Prep Kit for Illumina<sup>®</sup>, New England BioLabs, Cat# E7645L

(contains NEBNext Ultra II End Prep Mix and buffer, Ultra II Q5 Master Mix, Ligation Enhancer, Ultra II Ligation Master Mix)

NEBNext<sup>®</sup> Multiplex Oligos for Illumina<sup>®</sup> (Dual Index Primers Set 1), NewEngland BioLabs, Cat# E7600S (contains i5/i7 indexes, Adaptor for Illumina and User enzyme)

NEBNext<sup>®</sup> Multiplex Oligos for Illumina<sup>®</sup> (Dual Index Primers Set 2), NewEngland BioLabs, Cat# E7780S (contains i5/i7 indexes, Adaptor for Illumina and User enzyme)

AMPure XP<sup>®</sup> paramagnetic beads for PCR purification, Beckman Coulter, Cat# A63881

Kapa Library Quantification Kit for Illumina<sup>®</sup>, Complete Kit (ROX Low), Roche, Cat# KK4873 (for Quantstudio 12k)

## Before start

Before starting, generate the "odd" and "even" primer pools as follows:

1. Fully resuspend lyophilised oligonucleotides in 1x TE to a concentration of 100 micromolar ( $\mu\text{M}$ ) , vortex thoroughly and spin down
2. Sort the odd and even primer sets into separate batches and label two 1.5-ml tubes
3. Starting with the even primer set add the volume ( $\mu\text{l}$ ) given overleaf to the pooled 1.5-ml tube (between 7.5 - 15). Repeat with the odds. Vortex and spin down. These are your pooled stocks.
4. Dilute the pools one in ten across several aliquots with molecular grade water. Vortex and spin down. These are your working primer pools.



SubArtic primers v3-2 091121.csv



## cDNA prep

1h 25m

- 1 In a freshly bleached Pre-PCR hood (ideally in an isolated clean room), that has been subjected to UV irradiation for  00:40:00 , add LunaScript and RNA sample to a tube/well as follows (making sure to include negative controls):

1h

|  | A                    | B         |
|--|----------------------|-----------|
|  | Component            | Vol/Rxn   |
|  | LunaScript Super Mix | 2 $\mu$ l |
|  | RNA sample           | 8 $\mu$ l |

- 1.1 Mix thoroughly by pipetting up and down several times, seal plate, and centrifuge briefly.
- 1.2 Incubate the reaction as follows (with heated lid):

10m

15m

|  | A               | B      |
|--|-----------------|--------|
|  | Temperatur<br>e | Time   |
|  | 25°C            | 2 min  |
|  | 55°C            | 10 min |
|  | 95°C            | 1 min  |
|  | 4°C             | Hold   |

## Multiplex PCR

4h 10m

- 2 Primers are separated into two pools, odd and even, depending on where they sit across the Spike region. See the guidelines section for further details. In a clean pre-PCR hood set up two PCR reactions, one per pool as follows:

30m



| A                                | B           | C           |
|----------------------------------|-------------|-------------|
| Component                        | Rxn 1       | Rxn 2       |
| Pool Even                        | 0 µl        | 1.75 µl     |
| Pool Odd                         | 1.75 µl     | 0 µl        |
| Q5 Hot Start Hi-Fi 2x Master Mix | 6.25 µl     | 6.25 µl     |
| <b>Total</b>                     | <b>8 µl</b> | <b>8 µl</b> |

- 2.1 Add  4.5 µL of cDNA to respective wells to give a total volume of  12.5 µL . Include negative controls. Mix by pipetting, seal plate and centrifuge briefly.

10m

- 2.2 Perform PCR using the following program (with heated lid):

3h 30m

| A                     | B           | C     | D      |
|-----------------------|-------------|-------|--------|
| Step                  | Temperature | Time  | Cycles |
| Initial Denaturation  | 98°C        | 30 s  | 1      |
| Denaturation          | 98°C        | 15 s  | 35     |
| Annealing & Extension | 60°C        | 5 min |        |
| Hold                  | 4°C         | Hold  |        |

#### Note

Fewer PCR cycles can be used for samples of higher concentration.

## PCR pooling

30m

- 2.3 In a freshly bleached post-PCR hood that has been subjected to UV for  00:40:00 , combine Pool Even and Pool Odd to give  20 µL in each tube/well. Add  20 µL nuclease-free water to dilute products 1 in 2.

20m

**Note**

At least a 1 in 5 dilution is advised when using highly concentrated clinical samples subjected to a 35x cycle PCR.

- 2.4 Check the quality and concentration of negatives and a selection of samples using a fluorometer and/or Agilent TapeStation.

10m

**Note**

Samples can be normalised at this point if even barcode representation is required, but at the cost of time.

**NEBNext Ultra II End Prep**

50m

- 3 Prepare a master mix of the reagents as below by multiplying volumes by the number of samples; add 10% to allow for pipetting error.

5m

| A                                         | B           |
|-------------------------------------------|-------------|
| Component                                 | Vol/PCR Rxn |
| NEBNext Ultra II End Prep Enzyme Mix      | 0.6 µl      |
| NEBNext Ultra II End Prep Reaction Buffer | 1.4 µl      |
| <b>Total</b>                              | <b>2 µl</b> |

- 3.1 Mix well by pipetting and centrifuge briefly.

5m

- 3.2 Combine  2 µL of End Prep master mix with  10 µL of amplified cDNA and mix by pipetting. Spin down, seal, place in a thermocycler and run the following program:

40m



| A                  | B           |
|--------------------|-------------|
| <b>Temperature</b> | <b>Time</b> |
| 20°C               | 15 min      |
| 65°C               | 15 min      |
| 4°C                | Hold        |

## Adapter Ligation

45m

- 4 Prepare adapter ligation master mix by adding volumes as detailed below for each sample; adding 10% to allow for pipetting error.

10m

| A                                    | B                  |
|--------------------------------------|--------------------|
| <b>Component</b>                     | <b>Vol/PCR rxn</b> |
| NEBNext Ultra II Ligation Master Mix | 6 µl               |
| NEBNext Ligation Enhancer            | 0.2 µl             |
| NEBNext Adapter                      | 0.5 µl             |
| <b>Total</b>                         | <b>6.7 µl</b>      |

- 4.1 Add  6.7 µL ligation mix to each  12 µL amplified cDNA/mastermix and mix by pipetting. Incubate at  20 °C for  00:15:00 .

25m

- 4.2 Add  1 µL of USER enzyme to the ligation mixture. Mix by pipetting, centrifuge briefly and incubate at  37 °C for  00:15:00 .

10m

## Magnetic Bead clean up

1h 15m

- 5 Increase volume of sample from  19.7 µL to  25 µL by the addition of nuclease-free water.

5m

- 5.1 Perform a 0.9x Ampure XP bead clean by adding  22.5 µL of Ampure XP beads and mix by pipetting. Incubate for  00:05:00 at  Room temperature . 10m
- 5.2 Transfer tube/plate to magnet and allow beads to clump for  00:05:00 . Remove supernatant taking care not to disturb pellet, and discard. 15m
- 5.3 Wash the beads with  100 µL 80% ethanol for  00:00:30 then remove ethanol by pipetting. Repeat one more time making sure to remove any residual ethanol. 10m
- 5.4 Allow beads to air dry for around  00:05:00 ; ensure that all the ethanol has evaporated. Do not let beads dry to the point of cracking, as this could affect the amount of DNA recovered. 5m
- 5.5 Add  17 µL of nuclease-free water to each well and mix by pipetting. Allow to stand for  00:05:00 then transfer to magnet and allow beads to clump for  00:05:00 . Remove  15 µL into a fresh low-bind tube. 20m
- 5.6 Perform Qubit assay to determine concentration. 10m

## Addition of NEBNext indexes

1h 5m

- 6 Add the following components to separate wells of a PCR plate. Make sure to use the indexes in appropriate wells so that each sample has unique i5 and i7 indexes. If using less than 12 samples, the i5 primer can be replaced by the Universal PCR primer. 20m

| A                              | B       |
|--------------------------------|---------|
| Component                      | Vol/Rxn |
| Adapter ligated cDNA fragments | 7.5 µl  |
| NEBNext Ultra II Q5 Master Mix | 12.5 µl |
| Index Primer i7                | 2.5 µl  |



| A                | B            |
|------------------|--------------|
| Index Primer i5* | 2.5 µl       |
| <b>Total</b>     | <b>25 µl</b> |

\* if using 12 samples or fewer replace index primer i5 with the Universal primer.

6.1 Mix thoroughly by pipetting and centrifuge briefly.

5m

6.2 Place on thermocycler and perform PCR using the following program:

40m

| A                    | B           | C      | D       |
|----------------------|-------------|--------|---------|
| Cycle Step           | Temperature | Time   | Cycles  |
| Initial Denaturation | 98°C        | 30 s   | 1       |
| Denaturation         | 98°C        | 10 s   | 3 - 15* |
| Annealing/Extension  | 65°C        | 75 s   |         |
| Final Extension      | 65°C        | 5 mins | 1       |
| Hold                 | 4°C         | Hold   |         |

\*Refer to NEBNext protocol for number of cycles required. (Examples: 100 ng = ~3 cycles, 50 ng = ~3–4 cycles, 10 ng = ~6–7 cycles)

## Pool

1h 15m

7 Pool  5 µL of each sample together in a fresh low-binding microfuge tube and mix. This will give a final volume of  60 µL for twelve samples. If larger numbers of samples are used then pool  5 µL of each and mix. Split into aliquots of about  120 µL .

10m

7.1 Perform a 0.7x bead clean by adding  42 µL or  84 µL magnetic beads, depending on sample quantity. Leave for  00:05:00 to allow beads to bind DNA.

10m



- 7.2 Place on a magnet and allow beads to clump for 00:05:00 . Remove supernatant and discard. 10m
- 7.3 Wash pellet with 100  $\mu\text{L}$  80% ethanol for 00:00:30 and remove. Repeat one more time. Allow the ethanol to dry off for about 00:05:00 , making sure the pelleted beads do not dry out, which can reduce DNA recovery. 10m
- 7.4 Add 30  $\mu\text{L}$  of nuclease-free water and mix by pipetting. Leave for 00:05:00 at Room temperature to allow DNA to elute from the beads. 10m
- 7.5 Place the tube/plate on the magnet and let the beads clump for 00:05:00 . Elute DNA into a fresh low-binding tube. 10m
- 7.6 Perform TapeStation assay to determine if primer dimers are present. If present, perform another 0.7x bead clean and check again. 15m

## qPCR

2h

- 8 Perform qPCR by using a Kapa Library Quant kit for Illumina from Roche. Reactions are scaled down to 10  $\mu\text{L}$  . Use the appropriate Kapa kit for the type of qPCR machine used, see <https://pim-eservices.roche.com/eLD/api/downloads/ca670ceb-fb38-eb11-0291-005056a71a5d?countryIsoCode=pi> 1h 30m
- 8.1 Enter the values obtained from qPCR into the KAPA Library Quantification Data Analysis Template. 15m
- 8.2 Dilute sample to 4 nanomolar (nM) . 15m

## Library dilution

14m

- 9 Take 5  $\mu\text{L}$  of 4 nanomolar (nM) dilution of library and add 5  $\mu\text{L}$  of 0.1 Mass Percent NaOH to denature the DNA and incubate at 5m



🌡 Room temperature for ⌚ 00:05:00 .

9.1 Neutralise with 🧪 5  $\mu\text{L}$  of [M] 200 millimolar (mM) Tris-HCl (pH 7.0) for ⌚ 00:01:00 . 5m

9.2 Add 🧪 985  $\mu\text{L}$  cold HT1 (hybridisation buffer) to give [M] 20 picomolar (pM) denatured library 2m

9.3 Take 🧪 35  $\mu\text{L}$  of [M] 20 picomolar (pM) sample and add 🧪 465  $\mu\text{L}$  of HT1. This will give 🧪 500  $\mu\text{L}$  of [M] 1.4 picomolar (pM) sequencing mix. 2m

## PhiX control

18m

10 Thaw a tube of [M] 10 nanomolar (nM) PhiX stock . 5m

10.1 Combine the following volumes in a microcentrifuge tube:  
[M] 10 nanomolar (nM) PhiX (10  $\mu\text{L}$ ) and RSB (Resuspension Buffer) (15  $\mu\text{L}$ ). The total volume is 25  $\mu\text{L}$  at [M] 4 nanomolar (nM) .  
Vortex briefly and then pulse centrifuge. The [M] 4 nanomolar (nM) PhiX can be stored frozen for 3 months. 5m

10.2 Combine 🧪 5  $\mu\text{L}$  of [M] 4 nanomolar (nM) PhiX with 5  $\mu\text{L}$  of [M] 0.1 Mass Percent NaOH. Vortex briefly and pulse centrifuge. Incubate at 🌡 Room temperature for ⌚ 00:05:00 5m

10.3 Add 🧪 5  $\mu\text{L}$  of [M] 200 nanomolar (nM) Tris-HCl (pH 7.0 ), vortex briefly, pulse centrifuge and incubate at 🌡 Room temperature for ⌚ 00:01:00 1m

10.4 Add 🧪 985  $\mu\text{L}$  of HT1 (hybridisation buffer) to give a concentration of [M] 20 picomolar (pM) . Vortex briefly and pulse centrifuge. 1m



- 10.5 PhiX control at [1M] 20 picomolar (pM) can be frozen at -20°C for up to two weeks. After this time cluster numbers tend to decrease.

1m

## MiniSeq Sequencing

5m

- 11 Combine your library with PhiX control to give a final dilution of PhiX of 5%

5m