



Jul 21, 2018

Sudan Black staining of zebrafish larvae

 [PLOS Pathogens](#)

DOI

dx.doi.org/10.17504/protocols.io.rced2te

Emily Rosowski¹

¹University of Wisconsin - Madison



Emily Rosowski

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rced2te>

External link: <https://doi.org/10.1371/journal.ppat.1007229>

Protocol Citation: Emily Rosowski 2018. Sudan Black staining of zebrafish larvae. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rced2te>

Manuscript citation:

Rosowski EE, Raffa N, Knox BP, Golenberg N, Keller NP, Huttenlocher A (2018) Macrophages inhibit *Aspergillus fumigatus* germination and neutrophil-mediated fungal killing. PLoS Pathog 14(8): e1007229. doi: [10.1371/journal.ppat.1007229](https://doi.org/10.1371/journal.ppat.1007229)



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 27, 2018

Last Modified: July 21, 2018

Protocol Integer ID: 13414

Keywords: neutrophil staining, zebrafish larvae stain, neutrophil granule, sudan black, larvae

Abstract

Stains zebrafish larvae with Sudan Black B, which marks neutrophil granules.

Materials

MATERIALS

⊗ Phenol

⊗ Ethanol 100%

⊗ Sudan black B **Merck MilliporeSigma (Sigma-Aldrich) Catalog #199664**

⊗ Tween20

⊗ 16% Formaldehyde, methanol-free **Fisher Scientific Catalog #NC1040701**

⊗ Phosphate buffered saline powder, pH 7.4, for preparing 1 L solutions **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P3813**

⊗ Potassium hydroxide pellets

⊗ Hydrogen Peroxide, 30% **Fisher Scientific Catalog #H325-500**

Troubleshooting



Before start

You will need a rocker/nutator/rotator at 4°C and at RT.

Make both 1x and 10x stocks of PBS from powder packets, in water.

Make 70% Ethanol solution.

Make PBS-Tween = 1x PBS + 0.1% (v/v) Tween20.

Make 10% (w/v) Potassium hydroxide (KOH) in water.

Make Sudan Black B stock and working solution:

Stock: 0.18% (w/v) Sudan Black B in 70% Ethanol (store at RT, good for 1+ yr)

Working solution: 10 ml stock + 40 ml 70% Ethanol + 50 µl Phenol (store at RT, good for 6+ mo)

Anesthetize larvae

1. Anesthetize larvae with Tricaine.
2. Transfer to 1.5 ml microcentrifuge tube.

Fix larvae

2. 1. Make **fresh** 4% formaldehyde in PBS:
 - 1 ml 16% formaldehyde
 - 400 μ l 10x PBS
 - 2.6 ml H₂O
2. Pipet as much liquid off of larvae in tubes as possible.
3. Add 1 ml 4% formaldehyde per tube.
4. Cover in aluminum foil, rock at 4°C overnight.

Stain larvae

3. 1. Wash larvae in 1 ml 1x PBS, >5 min, rocking at 4°C.
2. Repeat for a total of 3 washes.
3. Pipet as much PBS as possible off of larvae, add 1 ml working solution of Sudan Black B (see Guidelines- Before you start).
4. Cover in aluminum foil, rock at RT for 1 hour.

Wash and Depigmentation

4. 1. Wash larvae in 1 ml 70% ethanol, 5 min, rocking at RT.
2. Repeat for a total of 3 washes.
3. Rehydrate larvae by washing in PBS-Tween, 5 min, rocking at RT.
4. Pipet as much PBS-Tween as possible off larvae, add 1 ml Depigmentation solution:
 - 400 μ l 10% KOH
 - 133 μ l 30% H₂O₂
 - 3.467 ml H₂O
5. Rock at RT for 5 minutes.
6. Take out a few larvae from tube, check depigmentation under a zoomscope. Pigment should be mostly lost, and sudan black signal should remain.
7. Wash larvae in 1 ml PBS-Tween, 5 min, rocking at RT.
8. Repeat for a total of 2 washes.
9. Store at 4°C. For longer-term storage change liquid to 1x PBS (no Tween).