



May 10, 2024

Version 2

Sugar Beet Tissue Collection for Genome Assembly and Annotation with Long Read Sequencing V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzbk92vx1/v2>

Olivia Todd¹, kevin.dorn¹

¹USDA-ARS



Olivia Todd

USDA-ARS

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzbk92vx1/v2>

Protocol Citation: Olivia Todd, kevin.dorn 2024. Sugar Beet Tissue Collection for Genome Assembly and Annotation with Long Read Sequencing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzbk92vx1/v2> Version created by **Olivia Todd**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: May 10, 2024

Last Modified: May 10, 2024

Protocol Integer ID: 99517

Keywords: RNA, OMNI-C, High molecular weight DNA, Sugar Beet, sugar beet tissue collection for genome assembly, sugar beet tissue collection, ars dorn lab protocol for beta vulgaris crop, wild relative genome sequencing, annotation with long read sequencing, beta vulgaris crop, sequencing, genome assembly, sugar beet, long read sequencing, tissue collection, various types of tissue collection, beta vulgari, annotation pipeline, ars dorn lab protocol, annotation with long read, annotation

Abstract

This protocol outlines the various types of tissue collection needed to follow the USDA-ARS Dorn lab protocol for *Beta vulgaris* crop wild relative genome sequencing. The annotation pipeline is in another protocol, and can be found by searching the author profile.

Materials

See materials list in the kit for each specified extraction.

QIAGEN RNeasy Plant MiniKit

PacBio Nanobind Kit

Safety warnings

! The high molecular weight DNA protocol uses betamercaptoethanol. Complete this extraction in a fume hood and dispose of hazardous waste properly.

Before start

It is best to plant out 20+ individuals of the line of interest, for all sequencing needed for this experiment. Be sure to keep the plants that you extract high molecular weight DNA from, as you can use that plant for the tissue collection for Omni-C later. It is best to have a highly homozygous population.



Harvesting High Molecular Weight Tissue for OMNIC sequencing

- 1 Place plant(s) of interest at 6 leaf stage or equivalent in a dark cabinet, drawer or room for 24-48 hours. Plants can be kept in the dark until tissue is fully etiolated or the consistency needed for successful extraction.
- 2 After 24-48 hours, harvest 5g young tissue into a 50ml falcon tube, removing the mid vein is optional. Place tube into -80C freezer overnight or flash freeze.
- 3 Submit tissue to sequencing company of choice.

Harvesting High Molecular Weight DNA for HiFi Sequencing

- 4 Place plant(s) of interest at 6 leaf stage or equivalent in a dark cabinet, drawer or room for 24-48 hours. Plants can be kept in the dark until tissue is fully etiolated or the consistency needed for successful extraction.
- 5 After 24-48 hours, harvest 1-5g young tissue into a 50ml falcon tube, removing the mid vein is optional. Place tube into -80C freezer overnight or flash freeze.
- 6 Follow the protocol below as written to isolate plant nuclei from 1-5g of dark treated tissue.
https://15a13b02-7dac-4315-baa5-b3ced1ea969d.filesusr.com/ugd/5518db_90f751986a4a4e7bbc8e648d467507d2.pdf?index=true
- 7 Repeat step 2 as necessary, until 1-5g of tissue has been collected. Additionally, you may choose to wait until the plant(s) of interest have more biomass. More than 5 grams of tissue may be needed, split over multiple reactions.
- 8 Place plants back in greenhouse to accumulate additional biomass if needed later.
- 9 Check for quality and presence of necessary DNA size for sequencing (>15kb) using a Pippin Pulse, running a 1% gel with Pre-set protocol example 3. Submit sample as specified by sequencing company of choice

Pippin Pulse Manual with protocols:

<https://sagescience.com/wp-content/uploads/2021/06/Pippin-Pulse-User-Manual-460022-Rev-F-1.pdf>

Tissue Collection and RNA Extraction for Iso-Seq

- 10 Determine plant treatments and which tissue types will be isolated for equimolar pooling of RNA.

e.g. To maximize transcript expression in sugar beet, two to three seedlings, 4 leaf stage plants and post-vernalized plants were harvested. Root and shoot tissue for adult plants were separated, ground with liquid N, and RNA extraction was completed, keeping root and shoot tissue separated for downstream pooling.
- 11 Complete RNA extraction with the RNeasy Plant Mini Kit from Qiagen.
- 12 Quality check RNA using Agilent Tapestation or gel electrophoresis 30S and 50S ribosome ratios. Determine approximate concentration with nanodrop.
- 13 Using sample submission data for sequencing company of choice and quality scores, determine how many individual RNA samples will be pooled into one tube for sequencing.
- 14 Uniformly dilute individual RNA samples such that the final concentration of the submission tube equals the amount required by the sequencing company.

e.g. Adding 5 uL of 5 ng/uL RNA for 10 samples will yield 250 ng of RNA in 50 uL (12.5 ug).