

Jun 20, 2025

Version 1

Sulfide Measurement Protocol V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvrq1oolmk/v1>

Gage R. Coon¹, Bosak Lab Protocols²

¹Method developed for Coon et al., in prep; ²Massachusetts Institute of Technology

Bosak Lab



Bosak Lab Protocols

Massachusetts Institute of Technology

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvrq1oolmk/v1>

Protocol Citation: Gage R. Coon, Bosak Lab Protocols 2025. Sulfide Measurement Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.6qpvrq1oolmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 20, 2025

Last Modified: June 20, 2025

Protocol Integer ID: 220658

Keywords: sulfide, cline assay, sulfide measurement protocol photometric determination of sulfide, sulfide measurement protocol photometric determination, hydrogen sulfide, h₂s concentration, sulfide level, water sample, form of hydrogen sulfide, proportional to the h₂s concentration, absorbance of the methylene blue color, yielding methylene blue, h₂, methylene blue, methylene blue color, absorbance, more concentrated sample

Abstract

Photometric determination of sulfide (S²⁻) in the form of hydrogen sulfide (H₂S) in water samples: Based on diamine and H₂S yielding methylene blue when in the presence of Fe³⁺. Absorbance of the methylene blue color at 670 nm is proportional to the H₂S concentration. This protocol is optimized for sulfide levels of roughly 0 to 10 mM. You can use more concentrated samples as long as absorbance is below 1 so it still follows Beer's Law. Samples can be stored in the fridge and just add diamine directly later. If unsure about dilutions, save two samples or extra volume.

Materials

All water should be deoxygenated by flushing with N₂ for 15+ minutes. Store diamine reagent in a dark place at 4°C if possible.

- Diamine reagent: Add 3.75 g N,N-dimethyl-p-phenylenediamine sulfate and 5.625 g ferric chloride to a beaker of 250 mL 50/50 HCl/deoxygenated water, store in a dark container.
- 0.05M zinc acetate
- 1 mM working sodium sulfide standard: Dilute concentrated standard from fridge to 0.05 M zinc acetate for a final dilution. Do this quickly and keep the concentrated standard anoxic!



Procedure

- 1 Standard Curve
 - 1.1 Create the following standards in Eppendorf tubes in the following format: Standard concentration: volume of 0.05M zinc acetate, volume of 1 mM sodium sulfide standard.
 - 1.2 0 mM (blank): 1.00 mL, 0 μ L
 - 1.3 0.2 mM: 0.98 mL, 20 μ L
 - 1.4 0.5 mM: 0.95 mL, 50 μ L
 - 1.5 0.8 mM: 0.92 mL, 80 μ L
 - 1.6 1 mM: 0.90 mL, 100 μ L
- 2 Prepare 1 mL samples diluted 100x in 0.05M zinc acetate (10 μ L sample + 990 μ L 0.05M zinc acetate).
- 3 Add 33 μ L of the diamine reagent to each Eppendorf tube.
- 4 Vortex and store in the dark for 20+ minutes.
- 5 Measure absorbance at 670 nm.



Protocol references

```
@article{CLINE1969,  
  title = {SPECTROPHOTOMETRIC DETERMINATION OF HYDROGEN SULFIDE IN NATURAL WATERS1},  
  volume = {14},  
  ISSN = {1939-5590},  
  url = {http://dx.doi.org/10.4319/lo.1969.14.3.0454},  
  DOI = {10.4319/lo.1969.14.3.0454},  
  number = {3},  
  journal = {Limnology and Oceanography},  
  publisher = {Wiley},  
  author = {CLINE, JOEL D.},  
  year = {1969},  
  month = may,  
  pages = {454-458}  
}
```