



Oct 31, 2023

SureSelect XT HS2 DNA to prepare libraries for single-cell Whole Genome Sequencing (scWGS) after single-cell Whole Genome Amplification (scWGA)

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9p3qzg3e/v1>

Ester Kalef-Ezra^{1,2}, Ben Harvey³, Katherine Roper³, Christos Proukakis^{1,2}

¹Department of Clinical and Movement Neurosciences, UCL Queen Square Institute of Neurology, London, UK;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815;

³Diagnostics and Genomics Group, Agilent Technologies LDA UK Ltd

University College London



Ester Kalef-Ezra

University College London



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9p3qzg3e/v1>



Collection Citation: Ester Kalef-Ezra, Ben Harvey, Katherine Roper, Christos Proukakis 2023. SureSelect XT HS2 DNA to prepare libraries for single-cell Whole Genome Sequencing (scWGS) after single-cell Whole Genome Amplification (scWGA). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v9p3qzg3e/v1>

License: This is an open access collection distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this collection and it's working.

Created: October 09, 2023

Last Modified: May 31, 2024

Collection Integer ID: 89822

Keywords: ASAPCRN, sureselect xt hs2 dna protocol, sureselect xt hs2 dna, cell whole genome sequencing, whole genome sequencing, cell whole genome amplification, genomic dna sample, whole genome, genome, library preparation for illumina, bravo automated liquid handling platform, dna

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research (MJFF)
Grant ID: 000430

Disclaimer

Acknowledgements and Funders:

We thank ICH.ZCR UCL Genomics Sequencing and Agilent for technical support.

This research was funded in part by Aligning Science Across Parkinson's [Grant ID: 000430] through the Michael J. Fox Foundation for Parkinson's Research (MJFF).

Note

ALL IMAGES CAN BE ENLARGED BY CLICKING AT THEM

Critical notes!

Please follow Good Laboratory Practices.

- To prevent DNA contamination, clean all surfaces and equipment before use with DNA AWAY Surface Decontaminant, separate the lab in pre- and post-PCR areas, and use filtered sterile pipette tips.
- For each protocol step that requires theremoval of tube cap strips, reseal the tubes with a fresh strip of caps.
- In all master mixes to be used manually, include 5-10% excess volume in the calculations.

Abstract

We adapted the SureSelect XT HS2 DNA protocol to prepare libraries using as input material samples after single-cell Whole Genome Amplification (scWGA), instead of genomic DNA samples, for library preparation for Illumina single-cell Whole Genome Sequencing (scWGS). This protocol can be employed either manually (Section 2: Option A) or in combination with automation using the Bravo Automated Liquid Handling Platform (Section 2: Option B).

Attachments



[861-2221.pdf](#)

1.1MB



[2022_SureSelect_Auto...](#)

5MB



[SureSelect Enzymatic...](#)

393KB



[SureSelect XT HS2 DN...](#)

1.1MB


Guidelines

We have used this protocol for DNA products after single-cell Whole Genome Amplification (scWGA) by either PicoPLEX (Takara), ResolveDNA (BioSkryB) and droplet MDA (Samplicx). We recommend to bead purify the PicoPLEX and ResolveDNA products prior library preparation. In all cases, the amplicon quantity and size should be assessed using Qubit and TapeStation or other methods prior to SureSelect XT HS2 DNA library preparation.


Bravo error/notification messages

Two error/notification messages are encountered when VWorks Bravo software is initializing:

1. If you encounter the G-axis error message shown below, select Ignore and Continue, leaving device in current state.



2. If you encounter the W-axis error message shown below, select Retry.



Materials

Commercial Reagents:

1. Equipment and consumables for DNA quality control:

DNA quantification. Here we used Qubit (Thermo Fisher Scientific):

Equipment	
Qubit Fluorometer	NAME
Fluorometer	TYPE
Invitrogen	BRAND
Q33238	SKU
https://www.thermofisher.com/order/catalog/product/Q33238#/Q33238 ^{LINK}	

-  Qubit™ dsDNA BR Assay Kit Thermo Fisher Scientific Catalog #Q32853
-  Qubit® dsDNA HS Assay Kit Thermo Fisher Scientific Catalog #Q32854
-  Qubit™ Assay Tubes Invitrogen - Thermo Fisher Catalog #Q32856

2. Nucleic acid analysis platforms. Here we used TapeStation (Agilent):

Equipment	
4150	NAME
TapeStation System	TYPE
Agilent	BRAND
G2992AA	SKU
https://www.agilent.com/store/en_US/Prod-G2992AA/G2992AA ^{LINK}	



-  Optical tube strips (8x Strip) **Agilent Technologies Catalog #401428**
-  Optical tube strip caps (8x strip) **Agilent Technologies Catalog #401425**
-  D1000 ScreenTape **Agilent Technologies Catalog #5067-5582**
-  D1000 Reagents **Agilent Technologies Catalog #5067-5583**
-  High Sensitivity D1000 ScreenTape **Agilent Technologies Catalog #5067-5584**
-  High Sensitivity D1000 Reagents **Agilent Technologies Catalog #5067-5585**

SureSelect XT HS2 performed manually:

-  SureSelect XT HS enzymatic fragmentation kit, 16 reactions **Agilent Technologies Catalog #5191-4079**
or
-  SureSelect Enzymatic Fragmentation Kit, 96 reactions **Agilent Technologies Catalog #5191-4080**
- SureSelect XT HS2 DNA Library Preparation Kit with Index Primer Pairs (Choose one of the following: Agilent G9981A (index 1-16), G9985A (index 1-96), G9985B (index 97-192), G9985C (index 193-288) or G9985D (index 289-384))

SureSelect XT HS2 with automation (Bravo):

-  SureSelect Enzymatic Fragmentation Kit, 96 Reactions, Auto **Agilent Technologies Catalog #5191-6764**
- SureSelect XT HS2 DNA Library Preparation Kit with Index Primer Pairs (Choose one of the following: Agilent G9985A (index 1-96), G9985B (index 97-192), G9985C (index 193-288) or G9985D (index 289- 384))

Consumables:

-  UltraPure[®] DNase/RNase-Free Distilled Water **Thermo Fisher Catalog #10977049**
-  Low TE Buffer **Invitrogen - Thermo Fisher Catalog #12090-015**
-  Ethanol (molecular biology grade, ≥99.8%) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #51976-500ML-F**
- Plasticware compatible with the selected thermal cycler. Here we used PCR 8-tube strips w/o caps
 8-tube strips for qPCR/PCR, without caps **VWR International (Avantor) Catalog #732-1517**
- Low binding filtered tips (sterile)
- Low-binding 1.5ml DNase/RNase-free (sterile). Here we used DNA LoBind Tubes,
 DNA LoBind[®] Tubes **Eppendorf Catalog #0030108418**
- Powder-free gloves
-  DNA AWAY[®] Surface Decontaminant, Surface decontaminant; 8.5 oz. (250mL) **Thermo Fisher Catalog #7010PK**
- Cleaning wipes for surface cleaning. Here we used Conti Washcloth Dry Brosch Direct PH5959

Consumables for SureSelect XT HS2 with automation:



- Processing plate:

- ⊗ Armadillo PCR Plate, 96-well, clear, clear wells, barcoded **Thermo Fisher Scientific Catalog #BC2396**

- PCR plate:

- ⊗ PCR plates, 96-well, Armadillo™ Standard Semi-skirted **VWR International (Avantor) Catalog #731-0649**

- Deep well plate for AMPure beads and master mixes:

- ⊗ Nunc® DeepWell® Plates with Shared-Wall Technology, 96U, 1mL, sterile **Thermo Fisher Catalog #260251**

- Waste Plate:
 - ⊗ Axygen™ Storage Microplates **Fisher Scientific Catalog #P-2ML-SQ-C**

- - ⊗ VersiCap Mat, 96-well, domed cap strips **Thermo Fisher Scientific Catalog #AB1810**

- - ⊗ Agilent Bravo 250uL Pipette Tips, Sterile, Filtered, for 96LT head **Catalog #19477-022**

- - ⊗ Reservoir, single cavity, polypropylene, 300 mL, 96 pyramids base geometry, 44 mm height, 25/pk **Agilent Technologies Catalog #201244-100**

Equipment:

- Thermal Cycler with 96-well, 0.2 ml block. Here we used Multigene Optimax Gradient Thermal cycler with 96 well block (Labnet MB0520)
- Manual version: Magnetic separator compatible with 0.2 ml tubes, strips of plates. Here we used either Magnetic Separator - PCR Strip (Takara 635011) or ResolveDNA Dual Volume Strip Tube Magnet (BioSkray PN100226)
- Plate or strip tube centrifuge. Here we used SciSpin Mini microfuge, blue, 7000rpm (SciQuip SS-6050) and plate centrifuge (Starlab N2631-0008)
- 1.5 ml tube centrifuge. Here we used Mini-Centrifuge (Fisherbrand 16617645)
- PCR cooler. Here we used PCR-Cooler, Blue, Capacity: 0.2 mL (Eppendorf 1019228)
- General lab pipettes (single and multi-channel)
- Nuclease-free filtered pipette tips sterile, if possible, use low binding tips
- 96-well plate mixer. Here we used Vortex mixer. Here we used WhirliMixer (Fisons Scientific 1993-520) and iSwix VT Digital Vortex Mixer (Appleton Woods ST6000)
- Ice bucket
- For automated version: Bravo NGS Option A robot (Agilent G5573A)

- ⊗ SMARTer-Seq™ Magnetic Separator - PCR Strip **Takara Bio Inc. Catalog #635011**

Equipment

SciSpin MINI Microfuge BLUE, inc 2 rotors of SS-6058 & SS-6059^{NAME}

SciQuip^{BRAND}

SS-6050^{SKU}

<https://sciquip.co.uk/scispin-mini-microfuge.html>^{LINK}

Equipment

Plate Centrifuge^{NAME}

StarLab^{BRAND}

N2631-0008^{SKU}

<https://www.starlabgroup.com/en/product/plate-centrifuge-n2631-0008.html>^{LINK}

Equipment

Mini-Centrifuge 100-240V, 50/60Hz Universal Plug, Grey^{NAME}

minicentrifuge^{TYPE}

Fisherbrand™^{BRAND}

16617645^{SKU}

<https://www.fishersci.co.uk/shop/products/fisherbrand-standard-mini-centrifuge/16617645>^{LINK}

Equipment

iSwix VT Digital Vortex Mixer

NAME

Appleton

BRAND

ST6000

SKU

<https://www.appletonwoods.co.uk/product/iswix-vt-digital-vortex-mixer-parent/>^{LINK}

Safety warnings



Safety Warnings:

Please follow the Safety Data Sheets (SDS) for all reagents for safe handling and safety hazards.

Disclaimer

Acknowledgements and Funders:

We thank ICH.ZCR UCL Genomics Sequencing and Agilent for technical support.

This research was funded in part by Aligning Science Across Parkinson's [Grant ID: 000430] through the Michael J. Fox Foundation for Parkinson's Research (MJFF).

Note

ALL IMAGES CAN BE ENLARGED BY CLICKING AT THEM

Critical notes!

Please follow Good Laboratory Practices.

- To prevent DNA contamination, clean all surfaces and equipment before use with DNA AWAY Surface Decontaminant, separate the lab in pre- and post-PCR areas, and use filtered sterile pipette tips.
- For each protocol step that requires theremoval of tube cap strips, reseal the tubes with a fresh strip of caps.
- In all master mixes to be used manually, include 5-10% excess volume in the calculations.

Attachments



861-2221.pdf
1.1MB



2022_SureSe
lect_Auto...
5MB



SureSelect
Enzymatic...
393KB



SureSelect
XT HS2 DN...
1.1MB

Files

 SEARCH

Protocol

NAME
Section 1: Enzymatic DNA Fragmentation (Manually)

VERSION 1

CREATED BY



Ester Kalef-Ezra
University College London

OPEN →

Protocol

NAME
Section 2: NGS library preparation for sequencing

VERSION 1

CREATED BY



Ester Kalef-Ezra
University College London

OPEN →

Protocol

NAME
Section 3: Libraries quality control (QC)

VERSION 1

CREATED BY



Ester Kalef-Ezra
University College London

OPEN →



Protocol references

This protocol was adapted with minor modifications from the following Agilent protocols:

SureSelect Enzymatic Fragmentation Kit

<https://www.agilent.com/cs/library/usermanuals/public/G9702-90050.pdf>

Manual version: SureSelect XT HS2 DNA System

<https://www.agilent.com/cs/library/usermanuals/public/G9985-90000.pdf>

Automated version: SureSelect XT HS2 DNA System Automated using Agilent NGS Bravo Option

A <https://www.agilent.com/cs/library/usermanuals/public/G9985-90020.pdf>