



Jun 22, 2023

Surgical procedures for simultaneous optogenetic manipulation in the SNc and acute fiber photometry in the dorsal striatum

DOI

dx.doi.org/10.17504/protocols.io.e6nvwd6qdlmk/v1

Qianzi Elena He¹

¹Northwestern University, Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815



Qianzi He

Northwestern University, Aligning Science Across Parkinson's...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwd6qdlmk/v1>

Protocol Citation: Qianzi Elena He 2023. Surgical procedures for simultaneous optogenetic manipulation in the SNc and acute fiber photometry in the dorsal striatum. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwd6qdlmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 22, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 83886

Keywords: ASAPCRN, acute fiber photometry in the dorsal striatum, optogenetic stimulation, simultaneous optogenetic stimulation, surgical procedures for simultaneous optogenetic manipulation, simultaneous optogenetic manipulation, craniotomy above the dorsal striatum, dorsal striatum, acute fiber photometry, optic cannula implantation into the snc, optic cannula implantation, surgical procedure, craniotomy, green dopamine sensor, viral injection

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research (MJFF)

Grant ID: ASAP-020600

Disclaimer

Appropriate ethics approval must be obtained before undertaking these experiments.

Abstract

This protocol describes the three surgical procedures required for **simultaneous optogenetic stimulation** in the Substantia Nigra pars compacta (SNc) and **acute fiber photometry** in the dorsal striatum:

1. **Viral injection** of a red-shifted channelrhodopsin, ChRmine, in the SNc, and a green dopamine sensor, GRAB-DA3m, in the dorsal striatum.
2. **Fiber-optic cannula implantation** into the SNc for optogenetic stimulation.
3. **Craniotomy** above the dorsal striatum for acute fiber photometry.

Troubleshooting

Viral injection for ChRmine and GRAB-DA3m

- 1 Perform unilateral ChRmine and GRAB-DA3m injection following the previous viral injection protocol (<https://www.protocols.io/view/midbrain-viral-injections-for-striatal-fiber-photo-5qpvr8zxv4o/v1?step=10> step 1-37).
- 2 In short, Anesthetize the mouse with isoflurane (5% for induction, 1-2% for maintenance with oxygen flowing at 0.8 liter per minute) and administer analgesics (0.02ml slow-release Buprenorphine 1mg/kg, and 1ml 0.1% lidocaine, both administered s.c.)
- 3 Secure the mouse on the stereotaxic frame and remove the skin on top of the skull.
- 4 Perform a 1mm circular craniotomy on either the left or the right SNc (coordinates from Bregma: AP -3.1, ML \pm 1.55) and another 1mm craniotomy above the ipsilateral striatum (coordinates from Bregma: AP +0.5mm, ML \pm 1.8mm).
- 5 Inject AAV-ChRmine($1E+12$ vg/ml) in a glass micropipette connected to a syringe and a pressure meter in the SNc craniotomy site at the four following depths (DV -3.8mm, -4.1mm, -4.4mm, -4.7mm). Inject 0.1ul of AAV-ChRmine at each depth, and 0.4ul in total.
- 6 Inject GRAB-DA3m ($2E+12$ vg/ml) with the same method in the striatum craniotomy site at the three following depths (DV -1.8mm, -2.2mm, -2.6mm). Inject 0.1ul at each depth, and 0.3ul in total.

Cannula and headplate implantation

- 7 Immediately after the viral injection, implant a borosilicate mono fiber-optic cannula (diameter = 400um, length = 4mm, NA = 0.66; Doric) in the same craniotomy site as the SNc injection.
- 8 First, secure the cannula with a stereotaxic cannula holder (Doric) on the stereotaxic frame. Place a metal headplate on the skull before lowering the cannula.
- 9 Lower the cannula till it is just above the cortex, adjusting the orientation so that the fiber is centered on the craniotomy, and the implant not obstructing the striatum craniotomy. Zero the z-coordinates on the stereotaxis frame.
- 10 Slowly lower the fiber implant till the bottom of the implant is firmly touching the skull (or the z-coordinate at -3.6, whichever comes first). Adjust the headplate to be perpendicular to the skull.



- 11 Secure the headplate and the implant with metabond and let it completely dry before loosening the cannula holder and raising it with the stereotaxic frame.
- 12 Finally, with a glass micropipette dipped in betadine, navigate to the striatum craniotomy coordinates, and gently mark it on the metabond. Drill a small circular dent with a hand drill, and mark the site with a marker for future craniotomy.
- 13 Remove the mouse from the stereotaxic frame and monitor till it recovers from anesthesia.
- 14 Wait for one week before training and at least two weeks before experiments.

Craniotomy for acute fiber photometry

- 15 One day prior to the experiment, when the mouse is already trained, perform craniotomy surgery for the acute fiber photometry following the same procedure described in previous fiber photometry protocol (<https://www.protocols.io/view/acute-striatal-or-midbrain-fiber-photometry-in-hea-4r3l27yj4g1y/v1>, step 23-28)
- 16 In short, after anesthetizing and administering analgesics to the mouse, place it on the stereotaxis frame.
- 17 Locate the previously marked dent on the metabond, corresponding to the striatum craniotomy site. With a hand drill, enlarge the dent till it is ~2mm in diameter. Slowly create a 2mm craniotomy, and make sure the dura mater is removed so it does not break the acute optic fiber during recording.
- 18 Stop any bleeding with PBS. Dry the site carefully with kimwipe and pointed cotton applicators. Seal the site with KWIK-SIL (WPI) to prevent the brain from being exposed.
- 19 Remove the animal from the stereotaxic frame and monitor till recovered from anesthesia.
- 20 Wait at least 24 hours before recording.



Protocol references

Maite Azcorra2023.Acute striatal or midbrain fiber photometry in head-fixed mice.**protocols.io**<https://dx.doi.org/10.17504/protocols.io.4r3l27yj4g1y/v1>

Maite Azcorra2023.Midbrain viral injections for striatal fiber photometry in mice.**protocols.io**<https://dx.doi.org/10.17504/protocols.io.5qpvor8zxv4o/v1>