

Nov 04, 2022

SVF Isolation and Immunophenotyping using Flow Cytometry in Leprosy Patients

DOI

dx.doi.org/10.17504/protocols.io.q26g7yym9gwz/v1

Sondang P. Sirait¹, Kusmarinah Bramono¹, Sri Linuwih Menaldi¹, Jeanne Adiwinata Pawitan^{2,3,4}, Wresti Indriatmi¹, Tiara Aninditha⁵

¹Dermatovenereology Department, Faculty of Medicine, Universitas Indonesia/Dr. Cipto Mangunkusumo General Hospital, Jakarta, Indonesia;

²Department of Histology, Faculty of Medicine, Universitas Indonesia, Jakarta, Indonesia;

³Stem Cells Medical Technology Integrated Service Unit, Dr. Cipto Mangunkusumo General Hospital/ Faculty of Medicine, Universitas Indonesia, Jakarta, Indonesia;

⁴Stem Cells and Tissue Engineering Research Center, Indonesia Medical Education and Research Institute (IMERI), Faculty of Medicine, Universitas Indonesia, Jakarta, Indonesia;

⁵Neurology Department, Faculty of Medicine, Universitas Indonesia/Dr. Cipto Mangunkusumo General Hospital, Jakarta, Indonesia



drsondangpenelitian

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7yym9gwz/v1>

Protocol Citation: Sondang P. Sirait, Kusmarinah Bramono, Sri Linuwih Menaldi, Jeanne Adiwinata Pawitan, Wresti Indriatmi, Tiara Aninditha 2022. SVF Isolation and Immunophenotyping using Flow Cytometry in Leprosy Patients. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7yym9gwz/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 02, 2022

Last Modified: November 04, 2022

Protocol Integer ID: 72207

Keywords: flow cytometry in leprosy patients adipose, using flow cytometry, flow cytometry, flow cytometry characterization, leprosy patients adipose, peripheral neuropathy, traditional flow cytometer, leprosy patient, immunophenotyping, neuropathy treatment, svf isolation, regenerative medicine, leprosy, ideal candidate for regenerative medicine, appropriate amount of mesenchymal cell, cd105 cell, mesenchymal stromal cell, liposuction, mononuclear cell

Funders Acknowledgements:

Universitas Indonesia

Grant ID: Hibah TADOK 2019

Abstract

Adipose derived stromal vascular fraction (SVF) contains a heterogeneous population of mononuclear cells, progenitor cells and about 1-10% are mesenchymal stromal cells. These cells are an ideal candidate for regenerative medicine for peripheral neuropathy. Leprosy is a disabling disorder with neuropathy, usually with consequences of permanent disability of the extremities. Before we could inject the SVF, we need to learn whether the SVF is appropriate to be used for neuropathy treatment using the flow cytometry characterization. Liposuction was done to four leprosy patients and adipose tissue was processed into SVF with a closed system. Following the process, the SVF were characterized and counted the CD73, CD90, CD34, CD45, Lin-Negatif, and CD105 cells using flow cytometry, cell counting, sterility and presence of mycobacteria. The process was done by ProStem using a traditional flow cytometer, without blocking and negative isotypes. The expected results from these protocols are the appropriate amount of mesenchymal cells needed to meet the IFATS and ISCT standards for SVF.

Troubleshooting

Isolation of Abdominal Adipose Tissue

- 1 The following procedure is performed after patient identification and obtaining patient consent.
- 2 The area was marked by using skin marker. After aseptic and antiseptic procedures the area was covered with a sterile surgical drape.
- 3 The local anesthesia was done with pehacain 2% and a small skin incision (1–3 mm) was made with number 11 blade afterwards.
- 4 The tumescent solution (TLA) was injected using a blunt tip cannula in the form of a fan-shaped. The tumescent solution is made based on Lilis formula with 1.000 mL NaCl 0.9%, 50 mL lidocaine 2%, 10 mL sodium bicarbonate 8.4%, and 1 mL epinephrine 1:1000.
- 5 The adipose tissue was collected using Mercedes aspiration cannula that connected with the canister and suction machine. The cannula, which is 3 mm in diameter, was directed in a radial, fan-shaped pattern. The fat above Scarpa's fascia was sucked out through the cannula with criss cross technique. Up to 100 mL of adipose tissue was collected.
- 6 The remaining tumescent solution is pushed caudally by pressing the operation area while the patient is standing.
- 7 The incision area was sutured and covered with sofratulle, gauze, sanitary pad, and Hypafix®. A compression bandage was applied using a large sized Hypafix®.

Isolation of Stromal Vascular Fraction

- 8 The tissue collected was processed into SVF using Sepax2® method.
- 9 The liposuction sample was transferred to a sterile collection bag without contact with air using a luer lock to luer lock connector. Then 0.15% collagenase solution (Serva) was added and incubated at 37 °C for 60 minutes.
- 10 Using CS-900.2 kit the incubated adipose tissue was connected to a Sepax2® device. Press "Start Procedure" to start a fully automated cellular separation.
- 11 For cell collection, the system diluted the residual pellets and extracted the content into the final bag, and automatically rinsed the chamber. After the procedure is completed, a



text "Remove Bags and Air Filter" will appear on the device. About twelve mL of SVF was collected in a sterile 20 mL syringe.

- 12 Samples were taken for count of live cells using Trypan Blue method, flow cytometric analysis, endotoxin testing using Lonza® kit, contamination assessments using Mycoplasma kit, and blood agar for bacterial culture.

Flow cytometry and cell counting

- 13 Cell count and viability are analyzed with **Automated Cell Counter Countess II** for enumerating the cells and cell counting.
- 14 After counting, 100.000 cells in a 100 µL medium were exported into the test tube.
- 15 The cells were added 20 µL positive antibody cocktail (contains CD90-FITC, CD73-APC, CD105-PerCP.Cy5.5, and 20 µL negative antibody cocktails (contains CD45-PE, CD34-PE, CD11b-PE, CD19-PE, and HLA-DR-PE.
- 16 Samples were incubated for 30 minutes and washed with Phosphate Buffer Saline (PBS) afterwards.
- 17 Samples in PBS were analyzed in Flow cytometer BD FACSCanto II using negative isotypes for the beads and positive isotypes for the CDs with enumeration control methods.
- 18 CD73+, CD90+, CD105+, Lineage Negative, CD34+, CD34+/CD45+ are counted and plotted into graphs using the methods.
- 19 The process was done by 1 clinical laboratory scientist, verified by another clinical laboratory scientist, and authorized by 1 clinical pathologist.