



Jan 09, 2026

## Synaptic vesicle preparation from mouse brain

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkyqndg5r/v1>

Akio Mori<sup>1</sup>, Robert Edwards<sup>1</sup>

<sup>1</sup>University of California San Francisco



Akio Mori

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.j8nlkyqndg5r/v1>

**Protocol Citation:** Akio Mori, Robert Edwards 2026. Synaptic vesicle preparation from mouse brain. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.j8nlkyqndg5r/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** November 20, 2025



**Last Modified:** January 12, 2026

**Protocol Integer ID:** 233001

**Keywords:** synaptic vesicle preparation from mouse brain, purification of synaptic vesicle, synaptic vesicle preparation, synaptic vesicle, brains of mice, mouse brain, glycerol gradient centrifugation, brain tissue, purification, gradient centrifugation

## Abstract

This protocol describes the preparation and purification of synaptic vesicles (SVs) from the brains of mice. Brain tissue is homogenized in sucrose-based SHT buffer, followed by sequential centrifugation steps to obtain the LP2 fraction. Synaptic vesicles are subsequently purified by glycerol gradient centrifugation, and SV-enriched fractions are collected for downstream biochemical analyses.

## Materials

- **Sucrose** ( $\geq 99.5\%$ ) (Sigma-Aldrich)
- **Mg-EGTA**
- **HEPES-Tris**, pH 7.4
- **cOmplete™ Protease Inhibitor Cocktail (11836153001, Roche)**
- **EGTA**
- **MgCl<sub>2</sub>**
- **Glycerol** (molecular biology grade)

## Buffers

- **SHT buffer:** 0.3 M sucrose, 1 mM Mg-EGTA, 10 mM HEPES-Tris, pH 7.4
- **Osmotic lysis buffer:** 8 mM HEPES-Tris, pH 7.4
- **Glycerol gradient buffer (5–25%):** 150 mM NaCl, 10 mM HEPES-Tris (pH 7.4), 1 mM EGTA, 0.1 mM MgCl<sub>2</sub>, 5–25% (v/v) glycerol (formed as a continuous gradient)



## Tissue homogenization

- 1 Dissect mouse brains and place them in ice-cold SHT buffer supplemented with protease inhibitor.
- 2 Homogenize the tissue in 10 mL SHT buffer using a glass–Teflon homogenizer ( 9 strokes at 900 rpm).

## Initial fractionation

- 3 Centrifuge homogenate at **871 × g, 10 min, 4°C**. Collect the supernatant (S1). 
- 4 Centrifuge S1 at **12,000 × g, 15 min, 4°C** to obtain the P2 pellet (crude synaptosomes). 

## Synaptosome washing and osmotic lysis

- 5 Resuspend the P2 pellet in SHT buffer and centrifuge at **14,500 × g, 15 min, 4°C**. Discard the supernatant 
- 6 Resuspend the pellet in 8 mM HEPES-Tris (pH 7.4) and incubate on ice for **45 min** to induce osmotic lysis. 
- 7 Homogenize briefly to complete disruption using a glass–Teflon homogenizer (five strokes at 3,000 rpm) .

## Preparation of LP2 fraction

- 8 Centrifuge the lysate at **32,500 × g, 20 min, 4°C** to obtain supernatant LS1. 
- 9 Centrifuge LS1 at **280,000 × g, 2 h, 4°C** to obtain the LP2 pellet. 
- 10 Resuspend the LP2 pellet in SHT buffer diluted to 40 mM sucrose (prepared by diluting standard SHT buffer: 0.3 M sucrose, 1 mM Mg-EGTA, 10 mM HEPES-Tris, pH 7.4).



- 10.1 Resuspension
  - a. Use a **200  $\mu$ L pipette tip cut at  $\sim 1/4$  from the distal end** and pipette up and down **100 times**.
  - b. Switch to a **200  $\mu$ L tip cut at  $\sim 1/8$  from the distal end** and pipette **100 times**.
  - c. Finally, use a **standard (uncut) 200  $\mu$ L tip** and pipette **100 times** to complete resuspension.
- 11 Aliquot the resuspended LP2, flash-freeze in liquid nitrogen, and store at  **$-80^{\circ}\text{C}$** .

## Glycerol gradient purification

- 12 Prepare a **5–25% glycerol gradient** using the glycerol gradient buffer using a gradient maker to form a continuous gradient.
- 13 Layer the LP2 sample onto the gradient.
- 14 Centrifuge at **270,000  $\times$  g for 1 h 4 min at  $4^{\circ}\text{C}$** .
- 15 Collect 0.5 mL fractions from the top using a gradient fractionator or connected fraction collector, ensuring smooth and continuous recovery of the gradient.
- 16 Use fractions 5–12, which usually contain purified synaptic vesicles, for subsequent experiments.

