

Oct 15, 2021

Synthesis of in vitro transcribed RNA from whole bacterial transcriptome

 [Nucleic Acids Research](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgb7r7yvpk/v1>

Bhargava Reddy Morampalli¹, Olin Silander¹

¹Massey University



Bhargava Reddy Morampalli

Massey University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgb7r7yvpk/v1>

External link: <https://doi.org/10.1093/nar/gkae601>

Protocol Citation: Bhargava Reddy Morampalli, Olin Silander 2021. Synthesis of in vitro transcribed RNA from whole bacterial transcriptome. [protocols.io](#) <https://dx.doi.org/10.17504/protocols.io.81wgb7r7yvpk/v1>

Manuscript citation:

Tan L, Guo Z, Shao Y, Ye L, Wang M, Deng X, Chen S, Li R (2024) Analysis of bacterial transcriptome and epitranscriptome using nanopore direct RNA sequencing. *Nucleic Acids Research* 52(15). doi: [10.1093/nar/gkae601](https://doi.org/10.1093/nar/gkae601)



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 10, 2021

Last Modified: May 18, 2026

Protocol Integer ID: 46157

Keywords: in vitro transcription, direct RNA sequencing, RNA modifications, cDNA for whole bacterial transcriptome, RNA modifications across the transcriptome, identifying RNA modification, bacterial transcriptome, whole bacterial transcriptome, RNA without any modification, native RNA, direct RNA, transcribed RNA, RNA, signals from native RNA, different ribosomal RNA, IVT RNA for entire transcriptome, ribosomal RNA, possible by direct RNA, nucleotides with epigenetic modification, trivial since bacterial transcriptome, entire transcriptome, transcriptome, different ribosomal RNAs for each type, getting IVT RNA, cDNA synthesis, using nanopore, epigenetic modification, vitro transcription, transcription, rRNA, nucleotide, Oxford Nanopore technology, bacteria, deep learning model, transcribed

Abstract

Identifying RNA modifications across the transcriptome has been made possible by direct RNA sequencing using nanopore sequencing from Oxford Nanopore Technologies. This is due to voltage changes happening when each nucleotide passes through the pore. Hypothesis is that since each nucleotide passing through the pore has a distinct change in voltage, nucleotides with epigenetic modifications (e.g. methylation) also have distinct voltage changes when passing through. This is identified by errors in basecalling.

Several software have been developed to identify locations of modifications using machine learning, deep learning models. In order to make highly accurate predictions, there is the necessity of comparing the signals from native RNA sequencing with RNA without any modifications. A few recent research publications get RNA through *in vitro* transcription but only for ribosomal RNAs (rRNA). This is easier since there are only 7 different ribosomal RNAs for each type i.e., 16S, 23S and 5S.

Getting IVT RNA for entire transcriptome is not trivial since bacterial transcriptome is complex and depends on the condition the bacteria has been grown in. We use the method of strand switching during cDNA synthesis to get cDNA for whole bacterial transcriptome which can then be used as a template for IVT.

Guidelines

Great precautions should be taken when handling RNA.



Materials

HiScribe™ T7 Quick High Yield RNA Synthesis Kit - E2050S

RNA Clean & Concentrator-5 - R1013

RNaseOUT

Primers: oligo-d(T)-22 primer, T7-Strand Switch Primer, PR2/T7 promoter primer (primer sequences commented in the protocol)

10 mM dNTP solution

Nuclease-free water (Thermofisher - #AM9937)

0.2 mL thin-walled PCR tubes

1.5 mL Eppendorf LoBind tubes

Maxima H Minus Reverse Transcriptase (200 U/uL) with 5x RT buffer (Thermofisher - #EP0751)

Pre-chilled freezer block at -20 C

LongAmp Taq 2X MasterMix (NEB - #M0287)

RNase Cocktail Enzyme Mix (Thermofisher - #AM2286)

RNAclean Agencourt XP beads

Before start

Clean the bench surfaces and pipettes using RNase ZAP to inhibit all RNases.

Always use RNase free centrifuge tubes, PCR tubes and pipette tips for all steps wherever required.



Addition of poly(A) tail to E. coli RNA

- 1 A reaction is setup in a PCR tube using the following components

	A	B	C	D	E	F	G	H
	Components	Volume (μL)						
	RNA	1 - 10 μg						
	10x E. coli poly(A) polymerase reaction buffer	2						
	ATP (10 mM)	2						
	RNase inhibitor	1						
	E. coli poly(A) polymerase	1						
	Water	upto 25 μL						
	Total	25						

- 2 Incubate the reaction mixture for 00:15:00 at 37 °C in a thermocycler followed by 00:20:00 at 65 °C for inactivation.

35m

- 3 Transfer the reaction mixture into a RNase free tube for cleanup step. Add 45 μL (1.8 volumes of reaction) of RNAClean AMPure XP beads and 105 μL (4.2 volumes) of 100 % ethanol to the reaction mixture

- 4 Flick the tube for mixing and incubate at Room temperature for 00:05:00 .

5m

- 5 Keep the tube on a magnetic rack and allow the beads to pellet towards the magnet

- 6 Prepare 500 μL of 80% ethanol during this time. Pipette the solution out of the tube taking care to not disturb the pellet



- 7 Wash the beads by pipetting  150 μL of 80% ethanol and repeat the step. Remove the ethanol, spin and put the tube back in the magnetic rack to remove any left over solution
- 8 Allow the beads to air dry for  00:00:30 making sure not to over dry. Add  20 μL of RNase free water to the beads and incubate at  Room temperature for  00:05:00 to elute RNA.
- 9 Place the tube back on the magnetic rack allowing the beads to separate and pipette the poly(A) tailed RNA in water into a new RNase free tube and keep it  On ice

5m 30s

cDNA first strand synthesis

1h 45m

- 10 Prepare a reaction with the following components in a PCR tube

A	B
Reagent	Volume
poly A+ RNA	1 μg (x μL)
anchored oligo d(T)22	1 μL
10 mM dNTPs	1 μL
RNase free water	9 - x μL
Total	11 μL

- 11 Incubate the reaction mixture at  65 $^{\circ}\text{C}$ for  00:05:00 in a thermocycler and then snap cool on a pre-chilled freezer block
- 12 In a separate tube, mix together the following:

5m



A	B
Reagent	Volume
5x RT Buffer	4 µl
RNaseOUT	1 µl
Nuclease-free water	2.5 µl
T7 - Strand-Switching Primer (SSP)	0.5 µL
Total	8 µl

13 Mix gently by *flicking the tube*, and *spin down*.

14 Add the strand-switching buffer to the snap-cooled, annealed mRNA, mix by *flicking the tube* and *spin down*

15 Incubate at  42 °C for  00:02:00

2m

16 Add 1 µl of Maxima H Minus Reverse Transcriptase. The total volume is now 20 µl.

17 Mix gently by *flicking the tube*, and *spin down*.

18 Incubate using the following protocol:

1h 35m

A	B	C	D
Cycle Step	Temperature	Time	No. of Cycles
Reverse transcription and strand-switching	42° C	90 mins	1
Heat inactivation	85° C	5 mins	1
Hold	4° C	∞	



RNA strand degradation and second strand synthesis

58m

- 19 Add  1 μL of RNase Cocktail Enzyme Mix (ThermoFisher, AM2286) to the reverse transcription reaction.
- 20 Incubate the reaction for  00:10:00 at  37 °C . 10m
- 21 Resuspend the AMPure XP beads by vortexing.
- 22 Transfer the sample to a clean 1.5 ml Eppendorf DNA LoBind tube.
- 23 Add  17 μL of resuspended AMPure XP beads to the reaction and mix by *flicking the tube*.
- 24 Incubate on a *Hula mixer* (rotator mixer) for  00:05:00 at  Room temperature . 5m
- 25 Prepare  500 μL of fresh 70% ethanol in nuclease-free water.
- 26 *Spin down* the sample and *pellet on a magnet*. Keep the tube on the magnet, and pipette off the supernatant.
- 27 Keep the tube on the magnet and wash the beads with  200 μL of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.
- 28 Repeat the previous step.
- 29 *Spin down* and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~  00:00:30 , but do not dry the pellet to the point of cracking. 30s



- 30 Remove the tube from the magnetic rack and resuspend pellet in  20 μL nuclease-free water.
- 31 Incubate on a *Hula mixer* (rotator mixer) for  00:10:00 at  Room temperature . 10m
- 32 *Pellet beads* on magnet until the eluate is clear and colourless.
- 33 Remove and retain  20 μL of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.
- 34 Prepare the following reaction in a 0.2 ml thin-walled PCR tube:

A	B
Reagent	Volume
2x LongAmp Taq Master Mix	25 μl
PR2 Primer (PR2) - T7 promoter primer (IDT)	3.75 μl
Reverse-transcribed sample from above	20 μl
Nuclease-free water	1.25 μl
Total	50 μl

- 35 Incubate using the following protocol: 17m

A	B	C
Temperature	Time	Cycle
94 °C	1 mins	1
50 °C	1 mins	1
65 °C	15 mins	1



	A	B	C
	4 °C	∞	

- 36 Resuspend the AMPure XP beads by vortexing.
- 37 Transfer the sample to a clean 1.5 ml Eppendorf DNA LoBind tube.
- 38 Add  40 μL of resuspended AMPure XP beads to the reaction and mix by *flicking the tube*.
- 39 Incubate on a *Hula mixer* (rotator mixer) for  00:05:00 at  Room temperature . 5m
- 40 Prepare  500 μL of fresh 70% ethanol in nuclease-free water.
- 41 *Spin down* the sample and *pellet on a magnet*. Keep the tube on the magnet, and pipette off the supernatant.
- 42 Keep the tube on the magnet and wash the beads with  200 μL of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.
- 43 Repeat the previous step.
- 44 *Spin down* and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~  00:00:30 , but do not dry the pellet to the point of cracking. 30s
- 45 Remove the tube from the magnetic rack and resuspend pellet in 21 μL nuclease-free water.
- 46 Incubate on a *Hula mixer* (rotator mixer) for  00:10:00 at  Room temperature . 10m

- 47 Pellet beads on magnet until the eluate is clear and colourless.
- 48 Remove and retain  21 µL of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.
- 49 Measure the concentration of the double stranded cDNA obtained after this step.

in vitro transcription

2h 15m

- 50 100 ng of double stranded cDNA synthesized is used for *in vitro* transcription.
- 51 New England Biolab's HiScribe™ T7 Quick High Yield RNA Synthesis Kit (E2050s) is used for RNA synthesis using cDNA as a template.
- 52 Thaw the necessary kit components, mix and pulse-spin in microfuge to collect solutions to the bottoms of tubes. Keep on ice.
- 53 Assemble the reaction at room temperature in the following order:

A	B
Components	Volume
Nuclease-free water	used to make up the volume to 20 µL
NTP Buffer Mix	10 µL (10 mM each NTP final)
Template DNA (cDNA)	X µL (100 ng - 1 µg)
T7 RNA Polymerase Mix	2 µL
Total Reaction Volume	20 µL

- 54 Mix thoroughly by pipetting and pulse-spin in a microfuge.



- 55 Incubate at 37°C for 02:00:00 . Incubation over 2 hours will no negatively affect the RNA synthesis 2h
- 56 $30\ \mu\text{L}$ of Nuclease-free water to the reaction mixture and $2\ \mu\text{L}$ of DNase I (included in the kit) is added after that.
- 57 Mix the reaction mixture by pipetting and incubate for 00:15:00 at 37°C 15m
- 58 Following the incubation, RNA is purified using any RNA cleanup kit. I used the RNA Clean & Concentrator-5 kit (R1013) and finally eluted in $20\ \mu\text{L}$ of RNase-free water.
- 59 Concentration is checked using Qubit RNA kit.