



Feb 01, 2024

α -Synuclein aggregation monitored by thioflavin-T (ThT) fluorescence in a plate reader

DOI

dx.doi.org/10.17504/protocols.io.8epv5x87dg1b/v1

Patricia Yuste-Checa¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5x87dg1b/v1>

Protocol Citation: Patricia Yuste-Checa, F Ulrich Hartl 2024. α -Synuclein aggregation monitored by thioflavin-T (ThT) fluorescence in a plate reader. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv5x87dg1b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 25, 2024

Last Modified: June 01, 2024

Protocol Integer ID: 94477

Keywords: ASAPCRN, synuclein aggregation, synuclein, plate reader, thioflavin, fluorescence

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol details how to efficiently monitor α -Synuclein aggregation by thioflavin T fluorescence in a plate reader.

Attachments



[α-Synuclein agg moni...](#)

96KB

Materials

Buffers:

▪ α -Synuclein fibril buffer:

	A	B
	NaN3	0.05%
	ThT	10 μ M
	KCl	150 mM
	Tris-HCl pH 7.6	50 mM


- [M] 1 millimolar (mM) Thioflavin T (ThT). Can be stored at -20°C .
- 6.5% NaN₃

96 well half-area plate of black polystyrene with a clear bottom **Corning Catalog #3881**

Troubleshooting






α -Synuclein aggregation monitored by thioflavin-T (ThT) fluorescence in a plate reader

14m

- 1 Mix reagents to a final concentration of [M] 200 micromolar (μ M) α -Synuclein, 0.05 % NaN_3 ,
[M] 10 micromolar (μ M) ThT in α -Synuclein fibril buffer. Prepare a mix for 4.5 reactions per condition (per condition, four technical replicates in each plate). Final volume is  80 μ L per well.

Note

Molecular or chemical chaperones can be included in the reaction in order to study their effect on α -Synuclein aggregation.

- 2 Dispense  80 μ L of the mix into a well of a 96 well half-area plate of black polystyrene with a clear bottom.
- 3 If possible, the outer wells should not be used and should be filled them with water.
- 4 Seal the plate with parafilm to avoid evaporation.
- 5 Set the following parameters in a SPARK multimode microplate reader (TECAN) and start the reaction:
 - Fluorescence measurement: ThT signal, excitation 440 nm, emission 480 nm, (use gain regulation) measured every  00:10:00 .
 - Temperature  37 $^{\circ}\text{C}$
 - Constant shaking ( 00:02:00 linear shaking: amplitude 1.5 mm, frequency 1080 rpm -  00:02:00 orbital shaking: amplitude 1 mm, frequency 510 rpm).

14m

Note

- SPARK multimode microplate reader (TECAN) is highly recommended due to its high sensitivity and the gain regulation mode that increases the fluorescence detection window. When using other plate readers, the sample ThT signal easily gets saturated even when reducing initial gain to the minimum gain capacity.
- Under these conditions α -Synuclein (A53T) slowly aggregates, reaching the ThT plateau after approximately 50 h aggregation.

