



Jan 09, 2026

α -Synuclein binding assay to isolated synaptic vesicles

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqp8mxvk5/v1>

Akio Mori¹, Robert Edwards¹

¹University of California San Francisco



Akio Mori

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqp8mxvk5/v1>

Protocol Citation: Akio Mori, Robert Edwards 2026. α -Synuclein binding assay to isolated synaptic vesicles. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgqp8mxvk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 19, 2025



Last Modified: January 09, 2026

Protocol Integer ID: 232990

Keywords: Synaptic vesicle, Pull-down assay, α -synuclein, assay to isolated synaptic vesicle, isolated synaptic vesicle, synaptic vesicle, synuclein, syn, binding assay, mouse brain

Abstract

The binding of α -synuclein (α -syn) to synaptic vesicles was assessed by pull-down using C-terminally His-tagged α -syn and Co Dynabeads (10103D, Invitrogen). Synaptic vesicles (SVs) were isolated from mouse brain (LP2 fraction), incubated with recombinant α -syn, and the bound fractions were analyzed by immunoblotting.

Materials

- Co Dynabeads (10103D, Invitrogen)
- Synaptic vesicles (SVs; LP2 fraction) purified from mouse brain
- Reaction buffer: 150 mM KCl, 1 mM MgCl_2 , 20 mM HEPES, pH 7.4
- BSA (A8806, Sigma)
- SHT buffer: 0.3 M sucrose, 1 mM Mg-EGTA, 10 mM HEPES-Tris, pH 7.4
- Recombinant His-tagged α -synuclein
- Wash buffer: 150 mM KCl, 1 mM MgCl_2 , 5 mM imidazole, 20 mM HEPES, pH 7.4
- Elution buffer: 150 mM KCl, 1 mM MgCl_2 , 200 mM imidazole, 20 mM HEPES, pH 7.4
- Antibody: synaptophysin (mouse monoclonal, 1:1,000; S5768, Sigma-Aldrich)

Beads Preparation

- 1 12.5 μ L per reaction of Co Dynabeads (10103D, Invitrogen) are transferred to a 1.5 mL low-adhesion microcentrifuge tube.
- 2 Wash the beads three times with 200 μ L reaction buffer (150 mM KCl, 1 mM MgCl₂, 20 mM HEPES, pH 7.4).
- 3 Pre-block Co Dynabeads with 1% BSA in reaction buffer for 30 min at room temperature with gentle rotation.
 - This step reduces nonspecific binding of SVs to the beads.
- 4 Wash the beads three times with reaction buffer.



Binding reaction

- 5 Resuspend the LP2 fraction in SHT buffer (1mg/mL)
- 6 Incubate purified recombinant His-tagged α -syn (0–25 μ g) with 25 μ g of LP2 in reaction buffer (total volume 100 μ L) for 3 h at room temperature with gentle agitation. Collect a mixture as input for Western blotting analysis.
- 7 Add the mixture to the pre-blocked beads and incubate for 15 min at room temperature with gentle rotation.



Washing and elution

- 8 Place tubes on a magnetic stand and collect the beads. Wash three times with 1 mL wash buffer.
- 9 Elute bound proteins by incubating beads in 25 μ L elution buffer for 15 min at room temperature with gentle rotation.
- 10 Spin briefly at 1,000 $\times$ g for 30 s.
- 11 Place on a magnetic stand and collect the eluted fraction.





Analysis

- 12 Analyze pull-down efficiency by immunoblotting.
- 13 Quantify band intensities using ImageJ.
- 14 Normalize synaptophysin levels in the elution fraction to the input.

