

Feb 03, 2025

α -synuclein PFF treatment

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7m821gwz/v1>

Maria Jose Perez J.¹, Federico Bertoli²

¹Institut Imagine; ²ASAP

Team Deleidi



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7m821gwz/v1>

Protocol Citation: Maria Jose Perez J., Federico Bertoli 2025. α -synuclein PFF treatment. protocols.io
<https://dx.doi.org/10.17504/protocols.io.q26g7m821gwz/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 26, 2024

Last Modified: February 03, 2025

Protocol Integer ID: 112831

Keywords: treatment

Funders Acknowledgements:

ASAP

Abstract

α -synuclein PFF treatment

α -synuclein PFF treatment

- 1 Shake Alpha-Synuclein Monomers in PBS for 7 days at 1,000 RPM. Solution should turn turbid during this period.
- 2 Spin fibrils for 10 minutes at max speed in benchtop centrifuge at Room temperature.
- 3 Discard supernatant and resuspend fibrils in PBS at a concentration of 5 ug/ul.
- 4 Validate α -synuclein PFFs via:
 - Electron microscopy
 - Sedimentation assay
 - Thioflavin-T assay (Cat. number T3516; Sigma–Aldrich)
- 5 Label α -synuclein PFFs with an Alexa Fluor 594 kit (Cat. number A10239; Thermo Fisher Scientific) following the manufacturer's instructions.

For microglial treatment:

- 6
 - Dilute α -synuclein PFFs in PBS to 100 μ g/mL.
 - Sonicate the solution (10-second pulses at 30% amplitude, six times with 2-minute intervals) using a Q700-220 sonicator (Qsonica, Newtown, CT, USA).
 - Dilute sonicated PFFs in microglial medium to a final concentration of 1 μ g/ μ L.
- 7 Add α -synuclein PFFs to iPSC-derived microglia cultures and cotreat with CDDO-Me, etoposide, or MCHA as needed.
- 8 Fix cells with 4% PFA for 10 minutes after 30, 60, and 90 minutes of treatment.
- 9 Perform immunostaining as described in the paper.
- 10 Acquire images:
 - For immunostained samples: use a Leica TCS SP8 confocal microscope with a 60 $\times$ /1.4 numerical aperture oil immersion objective.
 - For live-cell imaging: use an Evident APX100-HCU box-type microscope with a 20 $\times$ /0.8 numerical aperture objective. Capture images every 20 minutes for 3 hours.



- 11 For quantitative analysis:
 - Select three random fields containing 5–15 cells per field.
 - Image and analyze these fields across 3 independent experiments.
- 12 Perform particle analysis using Speckle Counting and Object Tracking pipelines in CellProfiler Image Analysis Software (version 4.2.6, RRID: SCR_007358) according to recommended guidelines.
- 13 Detect cell bodies using CD68 staining or brightfield imaging.