



T Cell Activation with anti-CD3 Antibodies Protocol - Human V.4

DOI

dx.doi.org/10.17504/protocols.io.793hr8n



Sam Li¹

¹BioLegend

BioLegend

Tech. support email: tech@biolegend.com



Sam Li

BioLegend

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account



DOI: <https://dx.doi.org/10.17504/protocols.io.793hr8n>

External link: <https://www.biolegend.com/protocols/t-cell-activation-with-anti-cd3-antibodies-protocol-human/4271/>

Protocol Citation: Sam Li . T Cell Activation with anti-CD3 Antibodies Protocol - Human. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.793hr8n>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: October 16, 2019



Last Modified: October 16, 2019

Protocol Integer ID: 28699

Keywords: cell

Guidelines

What if I want to stimulate my cells to a particular T-cell lineage?

- As a guide for this, you can take a look at our [Activation Bundles page](#), which will provide you with the listing of our products you will need to differentiate into different T-helper cell lineages.

Materials

MATERIALS

⊗ Anti-human CD3 Antibody: Clone UCHT1 (Ultra-LEAF™ format) **BioLegend Catalog #300437**

⊗ Anti-human CD3 Antibody: Clone HIT3a (Ultra-LEAF™ format) **BioLegend Catalog #300331**

STEP MATERIALS

⊗ 96-well flat-bottom tissue culture plates with lids **Corning Catalog #3596**

- Sterile PBS
- Cell culture medium (e.g., RPMI-1640 or IMDM supplemented with 10% FBS and 2mM L-glutamine)
- Sterile single-cell suspension of Ficoll-Hypaque-purified peripheral blood mononuclear cells, isolated T cells, or T cell subsets
- 96-well flat-bottom tissue culture plates with lids (e.g., Costar® Cat. No. 3596)
- **Note:** Soluble forms of Ultra-LEAF™ purified UCHT1 (1µg/ml) or Ultra-LEAF™ purified HIT3a (0.01 - 0.1µg/ml) may be used to activate T cells from PBMC cell populations.

Protocol materials

⊗ Anti-human CD3 Antibody: Clone UCHT1 (Ultra-LEAF™ format) **BioLegend Catalog #300437**

⊗ Anti-human CD3 Antibody: Clone HIT3a (Ultra-LEAF™ format) **BioLegend Catalog #300331**


⊗ 96-well flat-bottom tissue culture plates with lids **Corning Catalog #3596**


⊗ 96-well flat-bottom tissue culture plates with lids **Corning Catalog #3596**

Troubleshooting



- 1 Prepare a 10 µg/ml solution of anti-CD3 (clone UCHT1, OKT3, or HIT3a) in sterile PBS.
- 2 Dispense 50 µl of the antibody solution to each microwell of the 96-well assay plate. For the unstimulated control wells, add 50 µl of sterile PBS.

 96-well flat-bottom tissue culture plates with lids **Corning Catalog #3596**
- 3 Seal plate. Incubate at 37°C for 2 hours or 4°C overnight.

 02:00:00
- 4 Aseptically decant antibody solution from the microwell plate.
- 5 Wash plate microwells 3 times with sterile PBS. Discard liquid.
- 6 Prepare single cell suspension of cells of interest in supplemented cell culture medium to $1-2 \times 10^6/\text{ml}$.
- 7 Aliquot 200 µl cell suspension into plate microwells. Cover with lid. Incubate at 37°C in 5% CO₂ and 100% humidity for 3 days.