

T Cell Activation with anti-CD3 Antibodies Protocol -Human V.3

DOI

dx.doi.org/10.17504/protocols.io.tkbeksn



Kelsey Knight¹

¹BioLegend

BioLegend

Tech. support email: tech@biolegend.com



Sam Li BioLegend

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account



DOI: https://dx.doi.org/10.17504/protocols.io.tkbeksn

External link: https://www.biolegend.com/protocols/t-cell-activation-with-anti-cd3-antibodies-protocol-human/4271/

Protocol Citation: Kelsey Knight . T Cell Activation with anti-CD3 Antibodies Protocol - Human. **protocols.io** https://dx.doi.org/10.17504/protocols.io.tkbeksn

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: September 14, 2018



Last Modified: August 30, 2019

Protocol Integer ID: 15715

Keywords: cell

Materials

STEP MATERIALS

8 96-well flat-bottom tissue culture plates with lids Corning Catalog #3596

- Sterile PBS
- Anti-human CD3 Antibody
- Clone UCHT1 (LEAF™ format, Cat. No. 300413/300414/300432; Ultra-LEAF™ format, Cat. No. 300437/300438)
 - Clone OKT3 (LEAF™ format, Cat. No. 317303/317304/317315; Ultra-LEAF™ format, Cat. No. 317325/317326)
 - Clone HIT3a (LEAF™ format, Cat. No. 300313/300314; Ultra-LEAF™ format, Cat. No. 300331/300332)
- Cell culture medium (e.g., RPMI-1640 or IMDM supplemented with 10% FBS and 2mM L-glutamine)
- Sterile single-cell suspension of Ficoll-Hypaque-purified peripheral blood mononuclear cells, isolated T cells, or T cell subsets
- 96-well flat-bottom tissue culture plates with lids (e.g., Costar® Cat. No. 3596)
- * Soluble forms of LEAF™ purified UCHT1 (1 μg/ml) or LEAF™ purified HIT3a (0.01 0.1 μg/ml) may be used to activate T cells from PBMC cell populations.

Protocol materials

- Ø 96-well flat-bottom tissue culture plates with lids Corning Catalog #3596
- 8 96-well flat-bottom tissue culture plates with lids Corning Catalog #3596

Troubleshooting



- 1 Prepare a 10 μg/ml solution of anti-CD3 (clone UCHT1, OKT3, or HIT3a) in sterile PBS.
- 2 Dispense 50 µl of the antibody solution to each microwell of the 96-well assay plate. For the unstimulated control wells, add 50 µl of sterile PBS.
 - **⊠** 96-well flat-bottom tissue culture plates with lids **Corning Catalog #**3596
- 3 Seal plate. Incubate at 37°C for 2 hours or 4°C overnight.
 - **(2)** 02:00:00
- 4 Aseptically decant antibody solution from the microwell plate.
- 5 Wash plate microwells 3 times with sterile PBS. Discard liquid.
- 6 Prepare single cell suspension of cells of interest in supplemented cell culture medium to $1-2 \times 10^{6}$ ml.
- 7 Aliquot 200 µl cell suspension into plate microwells. Cover with lid. Incubate at 37°C in $5\% CO_2$ and 100% humidity for 3 days.