



Feb 29, 2024

T cell differentiation from mice spleen tissue

DOI

dx.doi.org/10.17504/protocols.io.261gedy8ov47/v1

Ningbo Zheng^{1,2}, Weiyi Peng^{1,2}

¹Department of Biology and Biochemistry, University of Houston, Houston, TX, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network



Sarfraz Ahmed

Mass General Hospital

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261gedy8ov47/v1>

Protocol Citation: Ningbo Zheng, Weiyi Peng 2024. T cell differentiation from mice spleen tissue. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261gedy8ov47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: February 29, 2024

Last Modified: February 29, 2024

Protocol Integer ID: 95978

Keywords: mice spleen tissue, spleen tissue, including differentiation, tissue

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: Grant ID: ASAP-000312

Abstract

This protocol is for T cell differentiation from mice spleen tissue to investigate T cell function (including differentiation and proliferation) in vitro.



Materials

1. Pre coating of 96 well plate with antibodies using following reagents

	A	B	C
	Reagent	Final concentration	Amount
	PBS	n/a	100mL/well
	Anti-CD3e (0.5mg/mL)	2.5mg/mL	0.5mL/well
	Anti-CD28 (1mg/mL)	1mg/mL	0.1mL/well

2. Prepare FACS buffer and T cell culture media (See Appendix 1 below)

	A	B	C
	FACS buffer		
	Reagent	Final concentration	Amount
	PBS	n/a	490mL
	FBS	2%	10mL
	Total	n/a	500mL

Note: Can be prepared in advance and stored at 4°C for at least 4 months. Filtration through a 0.2 µm vacuum filter is recommended.

	A	B	C
	T cell culture media		
	Reagent	Final concentration	Final concentration
	MEMα	n/a	429mL



	A	B	C
	FBS	10%	50mL
	GlutaMAX 100×	2×	10mL
	HEPES, 1M	0.02M	10mL
	2-mercaptoethanol (55 mM)	55μM	0.5mL
	Normocin, 1g (50mg/ml)	50μg /ml	0.5mL
	Total	n/a	500mL

Note: Can be prepared in advance and stored at 4°C for up to 1 month.

2. Prepare T-cell differentiation cocktails [1-4] (Appendixes 2 below)

	A	B	C
	Th0 differentiation cocktail (2.5×		
	Reagent	Final concentration	Amount
	T cell culture media	n/a	1mL
	IL-2 (5×10 ⁵ U/mL)	100U/mL	0.5mL

	A	B	C
	Th1 differentiation cocktail (2.5×		
	Reagent	Final concentration	Amount
	T cell culture media	n/a	1mL



	A	B	C
	IL-2 (5×105U/mL)	100U/mL	0.5mL
	IL-12 (10µg/mL)	4ng/mL	1mL
	Anti-IL-4 (10mg/mL)	10µg/mL	2.5mL

	A	B	C
	Th2 differentiation cocktail (2.5×)		
	Reagent	Final concentration	Amount
	T cell culture media	n/a	1mL
	IL-2 (5×105U/mL)	100U/mL	0.5mL
	IL-4 (10µg/mL)	10ng/mL	2.5mL
	Anti-IFNg (10mg/mL)	10µg/mL	2.5mL

	A	B	C
	Th9 differentiation cocktail (2.5×)		
	Reagent	Final concentration	Amount
	T cell culture media	n/a	1mL
	IL-2 (5×105U/mL)	100U/mL	0.5mL
	IL-4 (10µg/mL)	10ng/mL	2.5mL
	Anti-IFNg (10mg/mL)	10µg/mL	2.5mL

	A	B	C
	TGF- β (5 μ g/mL)	1ng/mL	0.5mL
	IL-1 β (10 μ g/mL)	10ng/mL	2.5mL

	A	B	C
	Treg differentiation cocktail (2.5\times)		
	Reagent	Final concentration	Amount
	T cell culture media	n/a	1mL
	IL-2 (5 \times 105U/mL)	300U/mL	1.5mL
	Anti-IFN γ (10mg/mL)	10 μ g/mL	2.5mL
	TGF- β (5 μ g/mL)	5ng/mL	2.5mL
	Anti-IL-4 (10mg/mL)	10 μ g/mL	2.5mL

Note: Cocktails can be prepared in advance and stored at 4°C for several days, seed 100mL in each well.

3. Reagents or Resources

	A	B	C
	MEM α	Fisher Scientific	Cat#32-571-101
	FBS	R&D Systems	Cat#S11150
	GlutaMAX 100 \times	Fisher Scientific	Cat#35050061
	HEPES	Fisher Scientific	Cat#MT25060CI
	2-mercaptoethanol	Fisher Scientific	Cat#21985023

	A	B	C
	Anti-CD3e	Cytek Biosciences	Cat#145-2C11
	Anti-CD28	BD Biosciences	Cat#553294
	IL-2	Prometheus Laboratories	Cat#NDC65483-116-07
	IL-12	PeproTech	Cat#210-12
	TGF- β 1	BioLegend	Cat#763104
	IL-1 β	R&D Systems	Cat#401-ML/CF
	IL-4	PeproTech	Cat#214-14
	Anti-IL-4	BioXCell	Cat#BE0045
	Anti-IFN γ	BioXCell	Cat#BE0055
	Mouse CD4+ T-cell Isolation Kit	STEMCELL Technologies	Cat#19852
	Mouse IFN-gamma DuoSet ELISA	R&D Systems	Cat#DY485
	Mouse IL-4 DuoSet ELISA	R&D Systems	Cat#DY404
	Mouse IL-9 DuoSet ELISA	R&D Systems	Cat#DY409
	Foxp3/transcription factor staining buffer set	Fisher Scientific	Cat#00-5523-00
	Anti-CD4-eFluor 450	eBioscience	Cat#48-0042-82
	Anti-CD25-PerCP	TONBO Biosciences	Cat#65-0251



	A	B	C
	Anti-Foxp3-PE	Fisher Scientific	Cat#12-5773-82
	PMA	Sigma-Aldrich	Cat#P1585

Troubleshooting



Day 0

- 1 Prepare precoated 96-well plate using antibodies described in materials section.

Day 1

- 2 Prepare single-cell suspensions of mice spleen tissues as per steps below
- 3 Euthanize mice by CO₂ inhalation or other means of euthanasia.
- 4 Collect the spleen from the mice.
- 5 Put a 40µm vacuum filter on a 50mL tube. Mesh the spleen on the strainer using the back of a syringe.
- 6 Add 5mL of RPMI medium to the mashed spleen to collect splenocytes into a 50mL tube.
- 7 Spin down cells at 1500rpm for 3min. Discard supernatant.
- 8 Add 0.5mL ACK lysis buffer to the red cell pellet, resuspend, and wait 30 seconds to lyse RBC.
- 9 Directly add 5mL T cell medium.
- 10 Spin down cells at 1500rpm for 3min. Discard supernatant.
- 11 Prepare a new 50mL tube with a new 40µm vacuum filter.
- 12 Resuspend cells in 10mL T cell medium and filter through a new 40µm vacuum filter.



13 Count the cells.

Isolation of naïve CD4⁺ T cells from single-cell suspensions from spleen tissue on Day 1

14 Naïve CD4⁺ T cells were isolated from single-cell suspensions from spleen tissue by negative selection using the EasySep Mouse CD4⁺ T-cell Isolation Kit (#19852, STEMCELL Technologies, Vancouver, Canada).

T cell differentiation on Day 1

15 Seed the cells into plate and culture with different conditions

16 Regulate the cell concentration at 5×10^5 /mL with T cell medium

	A	B	C
	Naïve CD4 ⁺ T cells	5×10^4 /well	100μL/well
	T cell medium		50μL/well
	Total	Total	150μL/well

17 Take out the pre-coated plate from 4°C.

18 Remove the supernatant and add 200μL PBS to the pre-coated plate and then remove PBS to wash the coated plate. Repeat 2 times.

19 Seed wells of a 96-well plate with 150μL cells cocktail (100μL cells and 50μL T cell media.

20 Seed 100μL differentiation cocktail.

21 Incubate for 5 days at 37°C.



Detect T cell differentiation by ELISA on Day 5 and Day 6 (Note: The chosen time point is not fixed, according to the experimental purpose)

22 Collect cells and count.

23 Regulate cell concentration at $1 \times 10^6/\text{mL}$.

24 Seed 100 μL cells into a new 96-well plate (U bottle)

25 Add 100 μL PMA cocktail and culture overnight. For non-stimulate cells, add 100 μL T cell medium.

	A	B	C
		Final concentration	Amount
	T cells	$1 \times 10^6/\text{mL}$	100 μL
	T cell Medium	n.a	100 μL
	PMA (50 $\mu\text{g}/\text{mL}$)	50ng/mL	0.2 μL

26 Collect supernatant for ELISA. (Note: The supernatant can be frozen at -20°C .)

Detect Treg cell differentiation by Flow cytometry on Day 5

27 Spin down cells at 1500rpm for 3min at RT.

28 Remove the supernatant over the sink in one motion.



29 Add 150µL FACS buffer and spin down cells at 1500rpm for 3min at RT. Remove the supernatant.

30 Surface staining, add 50µL FACS Buffer mAbs cocktail:

	A	B
	FACS buffer	50µL
	Anti-CD4-eFluor 450	0.3µL
	Anti-CD25-PerCP	0.3µL

31 Incubate at 4°C for 30min in the dark.

32 Directly add 150uL FACS buffer and spin down cells at 1500rpm for 3min at RT.

33 Remove the staining buffer and thoroughly resuspend cells.

34 Add 200µL of Foxp3 Fixation/Permeabilization working solution to each well for 1 hour at RT in the dark. (4×stock solution must be diluted prior to use with the Fixation/Permeabilization Diluent dilute 1 part concentrate with 3 parts diluents, make fresh).

35 Centrifuge samples at 400-600g for 5min at RT.

36 Add 200µL 1×Permeabilization Buffer (make fresh) to each well and centrifuge samples at 400-600g for 5min at RT. Discard the supernatant.

37 Add 50uL 1×Permeabilization Buffer mAb cocktail (Foxp3):

	A	B
	1×Permeabilization Buffer	50µL
	Anti-Foxp3-PE	0.3µL



- 38 Incubate at 4°C for 30min in the dark.
- 39 Directly add 200μL 1×Permeabilization Buffer to each well and centrifuge samples at 400-600g for 5min at RT. Discard the supernatant.
- 40 Add 200μL 1×Permeabilization Buffer for Flow Cytometry.

Protocol references

1. Liu X, Chen X, Zhong B, Wang A, Wang X, Chu F, Nurieva RI, Yan X, Chen P, van der Flier LG, Nakatsukasa H, Neelapu SS, Chen W, Clevers H, Tian Q, Qi H, Wei L, Dong C. Transcription factor achaete-scute homologue 2 initiates follicular T-helper-cell development. *Nature*. 2014 Mar 27;507(7493):513-8. doi: 10.1038/nature12910. Epub 2014 Jan 19. PMID: 24463518; PMCID: PMC4012617.
2. Xu T, Stewart KM, Wang X, Liu K, Xie M, Ryu JK, Li K, Ma T, Wang H, Ni L, Zhu S, Cao N, Zhu D, Zhang Y, Akassoglou K, Dong C, Driggers EM, Ding S. Metabolic control of TH17 and induced Treg cell balance by an epigenetic mechanism. *Nature*. 2017 Aug 10;548(7666):228-233. doi: 10.1038/nature23475. Epub 2017 Aug 2. PMID: 28783731; PMCID: PMC6701955.
3. Lu Y, Wang Q, Xue G, Bi E, Ma X, Wang A, Qian J, Dong C, Yi Q. Th9 Cells Represent a Unique Subset of CD4+ T Cells Endowed with the Ability to Eradicate Advanced Tumors. *Cancer Cell*. 2018 Jun 11;33(6):1048-1060.e7. doi: 10.1016/j.ccell.2018.05.004. PMID: 29894691; PMCID: PMC6072282.
4. Xue G, Zheng N, Fang J, Jin G, Li X, Dotti G, Yi Q, Lu Y. Adoptive cell therapy with tumor-specific Th9 cells induces viral mimicry to eliminate antigen-loss-variant tumor cells. *Cancer Cell*. 2021 Dec 13;39(12):1610-1622.e9. doi: 10.1016/j.ccell.2021.09.011. Epub 2021 Oct 21. PMID: 34678150; PMCID: PMC8678313.