

Sep 23, 2019

Taq PCR (Protocol for ZymoTaq™ PreMix)

DOI

dx.doi.org/10.17504/protocols.io.7cvhiw6

Alba Balletbó¹

¹Wageningen University

iGEM Wageningen 2019



Alba Balletbó

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.7cvhiw6

Protocol Citation: Alba Balletbó 2019. Taq PCR (Protocol for ZymoTaq™ PreMix). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.7cvhiw6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 15, 2019

Last Modified: September 23, 2019

Protocol Integer ID: 27765



Guidelines


ZymoTaq™ PreMix is supplied as a 2X concentrated “master mix” containing all the reagents needed to perform “hot start” PCR. The inclusion of a heat-activated, thermal-stable DNA polymerase reduces primer dimer and non-specific product formation that can occur when performing conventional PCR. This unique product is specifically designed for the amplification of bisulfite-treated DNA for methylation detection, realtime and quantitative PCR that are SYBR Green and probe based. The ZymoTaq™ PreMix yields specific amplicon formation with little or no byproducts. Simple and easy to use, just add water, primers, and template DNA to the ZymoTaq™ PreMix and then heat at 95°C for 10 minutes to initiate polymerization.

ZymoTaq™ DNA polymerase is a heat-activated, “hot start” polymerase that has 3'-terminal transferase activity. The addition of “A” overhangs to amplified DNA makes it ideal for use in TA-cloning.

Store ZymoTaq™ PreMix at -20°C for up to 12 months. Avoid repeated freeze/thawing of reagents. Prolonged storage is at -80°C

Materials

MATERIALS

 ZymoTaq™ PreMix - 50 rxns **Zymo Research Catalog #E2003**

 ZymoTaq™ PreMix - 200 rxns **Zymo Research Catalog #E2004**

	Reagent	Volume	Final concentration
	ZymoTaq™ PreMix	25 µL	1X
	Forward Primer (10 µM)	Variable	0.3 to 1 µM
	Reverse Primer (10 µM)	Variable	0.3 to 1 µM
	Template	Variable	< 200 ng/50 µL
	ddH ₂ O	to 50 µL	
	Total volume	50 µL	

- 1 Add all components in a 250 μ L tube making up to a 25 or 50 μ L reaction. If performing various PCR with different templates, a Master Mix is recommended to be done.

	Reag ent	Volu me	Final Conc entra tion
	Zymo Taq™ PreMi x	25 μ L	1X
	Forw ard Prime r (10 μ M)	Varia ble	0.3 to 1 μ M
	Rever se Prime r (10 μ M)	Varia ble	0.3 to 1 μ M
	Temp late	Varia ble	< 200 ng/50 μ L
	ddH2 O	to 50 μ L	
	Total Volu me	50 μ L	

- 2 Gently mix the PCR reactions and transfer the tubes to a thermocycler. Thermocycling conditions for a routine PCR:

	Step	Temp eratu re	Time	Cycle s
	Initial Denat uratio n		10 minut es	1
	Denat uratio n	94- 96°C	30 sec	30- 40
	Anne aling	Varia ble	30- 40 sec	
	Exten sion	72°C	30 sec - 1 minut e*	



Final Extension	72°C	7 minutes	1
Hold	4°C	> 4 minutes	1

* For ≤ 1 kb. Add an additional 15-30 seconds to the extension time for each kb > 1 kb. Make adjustments to the temperature and/or time if necessary.