

Aug 13, 2020

Target Guide Sequence Cloning Protocol Version 2

 Forked from [Target Guide Sequence Cloning Protocol](#)

DOI

dx.doi.org/10.17504/protocols.io.bjp3kmqn


Skye Waterland¹, Yang Li¹

¹Baylor College of Medicine



Skye Waterland

The University of Texas at Austin Dell Medical School



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bjp3kmqn

Protocol Citation: Skye Waterland, Yang Li 2020. Target Guide Sequence Cloning Protocol Version 2. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bjp3kmqn>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our workspace and it is working. This protocol is a revised version of the first Target Guide Sequence Cloning Protocol and is much easier to follow.

Created: August 13, 2020

Last Modified: August 13, 2020

Protocol Integer ID: 40411

Keywords: Lentivirus vector, cloning, vector digestion, oligo annealing, CRISPR,



Disclaimer

This protocol is a modified version of the Zhang Lab's *GeCKOv2* Target Guide Sequence Cloning Protocol attached below based off of Joung, J., Konermann, S., Gootenberg, J.*et al.* Genome-scale CRISPR-Cas9 knockout and transcriptional activation screening. *Nat Protoc* **12**, 828–863 (2017). <https://doi.org/10.1038/nprot.2017.016>

Other protocols modified/used in this protocol:

Bacterial Transformation, Addgene: <https://www.addgene.org/protocols/bacterial-transformation/>

More information about the specific lentiGuide-puro plasmid can be found here:

<https://www.addgene.org/52963/>.

Abstract

Create single gRNA vectors for targeted cloning utilizing CRISPR or CRISPR-based systems.

Attachments



[Lentivirus_Protocol....](#)

2.3MB



[Addgene_Protocol - ...](#)


164KB

Image Attribution


<https://www.addgene.org/52963/>

Materials


MATERIALS

 NEBuffer 3.1 - 5.0 ml **New England Biolabs Catalog #B7203S**

 T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**


 Agar


 lentiGuide-Puro **addgene Catalog #52963**

 double distilled water (ddH₂O)

 SOC Media

 1X TAE Buffer


 10X NEB T4 DNA ligase buffer **New England Biolabs**

 10X T4 PNK Reaction Buffer **New England Biolabs**


 ethanol

 10X PCR Buffer **Life Technologies Catalog #10966-034**

 LB-Broth Miller (= LB mix) **Formedium Catalog #LMM0104**

 One Shot[®]; TOP10 Chemically Competent E. coli **Thermo Fisher Catalog #C404010**

 BsmBI-v2 **New England Biolabs Catalog #R0739L**


 HotStarTaq Plus DNA Polymerase (1000) **Qiagen Catalog #203605**

 dNTP Set (100mM each A C G T) **Ge Healthcare Catalog #95038-256**


STEP MATERIALS

-

 lentiGuide-Puro **addgene Catalog #52963**

 double distilled water (ddH₂O)


 10X NEB T4 DNA ligase buffer **New England Biolabs**


 10X T4 PNK Reaction Buffer **New England Biolabs**


 BsmBI-v2 **New England Biolabs Catalog #R0739L**

 lentiGuide-Puro **addgene Catalog #52963**

 BsmBI-v2 **New England Biolabs Catalog #R0739L**

 Agar

 1X TAE Buffer

 double distilled water (ddH₂O)

 double distilled water (ddH₂O)



⊗ 10X NEB T4 DNA ligase buffer **New England Biolabs**

⊗ T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**

⊗ One Shot[®]; TOP10 Chemically Competent E. coli **Thermo Fisher Catalog #C404010**

⊗ ethanol

⊗ 10X PCR Buffer **Life Technologies Catalog #10966-034**

⊗ dNTP Set (100mM each A C G T) **Ge Healthcare Catalog #95038-256**

⊗ HotStarTaq Plus DNA Polymerase (1000) **Qiagen Catalog #203605**

⊗ ddH₂O

⊗ LB-Broth Miller (= LB mix) **Formedium Catalog #LMM0104**

⊗ SOC Media

lentiGuide-Puro: RRID:Addgene_52963

Sigma-Aldrich: RRID:SCR_008988



Protocol materials

☒ double distilled water (ddH₂O)

☒ lentiGuide-Puro **addgene Catalog #52963**

☒ ddH₂O

☒ LB-Broth Miller (= LB mix) **Formedium Catalog #LMM0104**

☒ 10X NEB T4 DNA ligase buffer **New England Biolabs**

☒ T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**

☒ ethanol

☒ dNTP Set (100mM each A C G T) **GE Healthcare Catalog #95038-256**

☒ 10X T4 PNK Reaction Buffer **New England Biolabs**

☒ Agar

☒ 1X TAE Buffer

☒ SOC Media

☒ BsmBI-v2 **New England Biolabs Catalog #R0739L**

☒ HotStarTaq Plus DNA Polymerase (1000) **Qiagen Catalog #203605**

☒ dNTP Set (100mM each A C G T) **GE Healthcare Catalog #95038-256**

☒ Agar

☒ 10X NEB T4 DNA ligase buffer **New England Biolabs**

☒ 10X T4 PNK Reaction Buffer **New England Biolabs**

☒ ethanol

☒ NEBuffer 3.1 - 5.0 ml **New England Biolabs Catalog #B7203S**

☒ T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**

☒ 10X PCR Buffer **Life Technologies Catalog #10966-034**

☒ double distilled water (ddH₂O)

☒ One Shot[®]; TOP10 Chemically Competent E. coli **Thermo Fisher Catalog #C404010**

☒ lentiGuide-Puro **addgene Catalog #52963**

☒ BsmBI-v2 **New England Biolabs Catalog #R0739L**

☒ BsmBI-v2 **New England Biolabs Catalog #R0739L**

☒ double distilled water (ddH₂O)

☒ One Shot[®]; TOP10 Chemically Competent E. coli **Thermo Fisher Catalog #C404010**

☒ 1X TAE Buffer



☒ SOC Media

☒ LB-Broth Miller (= LB mix) **Formedium Catalog #LMM0104**

☒ lentiGuide-Puro **addgene Catalog #52963**

☒ 10X PCR Buffer **Life Technologies Catalog #10966-034**

☒ double distilled water (ddH₂O)

☒ HotStarTaq Plus DNA Polymerase (1000) **Qiagen Catalog #203605**

☒ 10X NEB T4 DNA ligase buffer **New England Biolabs**

☒ lentiGuide-Puro **addgene Catalog #52963**

☒ ddH₂O

☒ NEBuffer 3.1 - 5.0 ml **New England Biolabs Catalog #B7203S**

☒ BsmBI-v2 **New England Biolabs Catalog #R0739L**

☒ lentiGuide-Puro **addgene Catalog #52963**

☒ Agar

☒ 1X TAE Buffer

☒ double distilled water (ddH₂O)

☒ 10X NEB T4 DNA ligase buffer **New England Biolabs**

☒ 10X T4 PNK Reaction Buffer **New England Biolabs**

☒ BsmBI-v2 **New England Biolabs Catalog #R0739L**

☒ T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**

☒ 10X NEB T4 DNA ligase buffer **New England Biolabs**

☒ ddH₂O

☒ double distilled water (ddH₂O)

☒ double distilled water (ddH₂O)

☒ ethanol

☒ One Shot™ TOP10 Chemically Competent E. coli **Thermo Fisher Catalog #C404010**

☒ 10X PCR Buffer **Life Technologies Catalog #10966-034**

☒ dNTP Set (100mM each A C G T) **GE Healthcare Catalog #95038-256**

☒ HotStarTaq Plus DNA Polymerase (1000) **Qiagen Catalog #203605**

☒ LB-Broth Miller (= LB mix) **Formedium Catalog #LMM0104**

☒ SOC Media

Before start

Design and order gRNA oligos from **Sigma-Aldrich** (RRID:SCR_008988).

Lentiviral vector digestion

1 Digest and dephosphorylate

lentiGuide-Puro addgene Catalog #52963

Created with SnapGene®




Image attribution: <https://www.addgene.org/52963/>

- 1.1 Add 40 μ L
- ddH2O
- to a 1.5 mL

Equipment

Snap Cap Microcentrifuge Tube or equivalent

NAME

Polypropylene Microcentrifuge Tube

TYPE

Corning Costar Snap Cap Microcentrifuge Tube

BRAND

07200210

SKU


<https://www.fishersci.com/shop/products/costar-microcentrifuge-tubes-6/07200210>

LINK


2 mL snap cap polypropylene micro tube

SPECIFICATIONS

1.2 Add  5 μ L of

 NEBuffer 3.1 - 5.0 ml **New England Biolabs Catalog #B7203S**
to solution.


1.3 Add  3 μ L of

 BsmBI-v2 **New England Biolabs Catalog #R0739L**
to solution.

Note

Note: we used BsmBI Cat #R0580, NEB, but this one has been discontinued. #R0739L is considered as an effective replacement.

1.4 Add  2 μ L of

 lentiGuide-Puro **addgene Catalog #52963**
to solution.

1.5 Close cap on microcentrifuge tube and place in

Equipment

Mini-centrifuge

NAME

Centrifuge

TYPE

Fisher

BRAND

S67601B

SKU


<https://www.fishersci.com/shop/products/fisherbrand-standard-mini-centrifuge-standard-mini-centrifuge/s67601b>

LINK

Any standard mini centrifuge with adapters for different tube sizes will suffice

SPECIFICATIONS



for  00:00:10 on until all of the solution is at the bottom of the tube.

1.6 Place microcentrifuge tube with vector digestion mixture in a



Equipment

Oven

NAME

Oven forced-air convection

TYPE

Fisher Isotemp

BRAND


15-103-0510

SKU

<http://www.fishersci.ca/shop/products/fisher-scientific-isotemp-general-purpose-heating-drying-ovens/151030510?keyword=true>

LINK


on  55 °C

1.7 Close lid and set a timer for  01:00:00

Gel purify the digested plasmid from Step 1

2 Prepare gel

2.1 Create gel concentration of 1.2-1.5%

 Agarfor  100 mL of 1X TAE Buffer

solution.

3 Run gel

3.1 Isolate 2kb and 8kb band. Collect 8kb (8318 bp) band for gel purification.



Preparing the gRNAs

4 Design gRNA sequence for CRISPR strategy using **CRISPR direct**.

5 Order oligos from **Sigma-Aldrich** (RRID:SCR_008988).

6 Create a dilution from stock oligos at a 1:10 ratio in

⊗ double distilled water (ddH₂O)

.

Note

Diluted oligo concentration should be 10 micromolar (μM) .

Phosphorylate and anneal each pair of oligos

7 Prepare phosphorylation/annealing reaction

7.1 Add 6.5 μL

⊗ double distilled water (ddH₂O)

to a microcentrifuge tube.

7.2 Add 1 μL each of 10 micromolar (μM) Oligo 1 (F), 10 micromolar (μM)

Oligo 2 (R), and

⊗ 10X NEB T4 DNA ligase buffer **New England Biolabs**

7.3

Add 0.5 μL

⊗ 10X T4 PNK Reaction Buffer **New England Biolabs**

7.4 Vortex and microcentrifuge

8 Place the phosphorylation/annealing reaction in a

Equipment

SimpliAmp Thermal Cycler

NAME

PCR

TYPE

Applied Biosystems

BRAND

A24811

SKU

<https://www.thermofisher.com/order/catalog/product/A24811>^{LINK}

Any standard PCR thermocycler will suffice

SPECIFICATIONS



Note

Settings should be as follows: 37°C for 00:30:00, 95°C for 00:05:00, and then ramp down to 25°C at $5^{\circ}\text{C} / 00:01:00$.

Setting up and incubating the ligation reaction

9 Making the ligation reaction

9.1 Place $4.8\ \mu\text{L}$ of

double distilled water (ddH₂O)





in a microcentrifuge tube.

9.2 Add  2.2 μL

 BsmBI-v2 **New England Biolabs Catalog #R0739L**

9.3 Add  1 μL each of



 10X NEB T4 DNA ligase buffer **New England Biolabs**

, diluted oligo duplex from  [go to step #7](#) , and

 T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**

.


9.4 Lightly vortex and microcentrifuge

10 Incubate the ligation reaction at room temperature for  02:00:00 -  03:00:00



Transformation into E. coli bacteria

11 Prepare LB Agar plates

11.1 Wipe down your bench with at least 70%


 ethanol

and light a bunsen burner.

11.2 Remove Agar Ampicillin Plates  250 μL from  4 $^{\circ}\text{C}$ and let warm up to room temperature.

12 E. coli competent cell transfection

12.1 Take competent cells

 One Shot™; TOP10 Chemically Competent E. coli **Thermo Fisher Catalog #C404010**



out of -80 °C and thaw on ice (00:20:00 - 00:30:00).

- 12.2 Add 100 µL of E.coli cells to 10 µL of DNA in a microcentrifuge tube next to the bunsen burner.

Note

When working with the E. coli, be very diligent and make sure you are working in the sterile area of your bench near the bunsen burner flame.

- 12.3 Gently flick tube a few times with your finger to mix.

- 12.4 Incubate the competent cell/DNA mixture on ice for 00:30:00 .

- 12.5 Place transformation tube(s) into water bath at 42 °C for 00:00:30 - 00:01:00 to heat shock E. coli cells.

- 12.6 Place the transformation tube(s) back on ice for 00:02:00 .

- 12.7 Add 250 µL LB-Broth Miller (= LB mix) **Formedium Catalog #LMM0104** (without antibiotic) or SOC Media to the tube(s).

- 12.8 Place tube(s) in 37 °C shaking incubator for 00:45:00 - 01:00:00 .

- 12.9 Plate all of transformation onto LB agar plate(s) with ampicillin. Incubate plates at 37 °C overnight.

Picking Colony for Suspension Growth

- 13 Before you get started



Note

Start this whole section (picking colony and suspension growth) sometime in the late afternoon (~4-5pm), so that you can run PCR for the bacterial plasmids you collect the next morning (~16 hours later). More than 24 hours of incubation can cause other non-ampicillin resistant bacteria to grow in the suspension tubes.

Note

Set up an aseptic area for handling the bacterial colonies.

14 Picking colony

14.1 Pick 10-15 colonies from each agar plate to suspend in an LB solution.

Note

We originally selected only 4 colonies from each plate, but did not have a successful PCR. To increase chances of successfully amplifying the plasmid vector, we suggest picking 10-15 colonies.

14.2 Place each colony in a tube with 3 mL of LB.

15 Culture the bacteria

15.1 Place tubes into shaking incubator at 37 °C overnight.

Run PCR to Identify Positive Clones


- 16 Prepare Master PCR Mix by adding each of the below reagents to a microcentrifuge tube.

Note

Before you start: multiply each volume below by the same amount (x 10, 15, etc.) according to how much master mix you need to run your sets.

16.1

 15.875 μL

 ddH₂O

16.2

 2.5 μL

 10X PCR Buffer **Life Technologies Catalog #10966-034**

16.3

 0.5 μL

[M] 10 millimolar (mM)

 dNTP Set (100mM each A C G T) **GE Healthcare Catalog #95038-256**

16.4

 2 μL

[M] 10 micromolar (μM) Forward Primer hU6-02

Note

Forward Primer hU6-02 sequence: TAATTAGAATTAATTTGACT
Ordered from **Sigma-Aldrich** (RRID:SCR_008988).

16.5

 2 μL


[M] 10 micromolar (μM) Reverse Primer (gRNA reverse primer oligo)

16.6



 0.125 μL

 HotStarTaq Plus DNA Polymerase (1000) **Qiagen Catalog #203605**

- 17 Loading and Prepping PCR wells

- 17.1 Aliquot equal proportions of master mix to each PCR well (we used  23 μL).



- 17.2 Add  2 μL of each bacterial suspension sample from  [go to step #15](#) to each PCR well.

Note

Total volume of Master Mix and Sample should be  25 μL in each well.

- 17.3 Seal off PCR wells tightly with a clear plastic cover or tube tops.

- 17.4 Centrifuge  1200 rpm, 00:01:00 .

- 18 Place PCR reaction in thermocycler and run PCR.

Note

Settings should be as follows:

- 1) 95C for 5mins (ramp up).
- 2) 95C for 30sec.
- 3) 55C for 45 sec.
- 4) 72C for 45 sec.
- x35 cycles steps 2-4.
- 5) 72 for 10mins.

Running Gel for PCR Product Verification

- 19 Prepare a gel of 2% concentration. Run gel and examine bands. Desired band length is about 200bp with the gRNA insertion.

Congrats!

- 20 You have successfully transformed a lentiviral vector with your gRNA sequence of interest! For confirmation, feel free to sequence your vector.

