

Feb 03, 2025

Targeted LC–MS metabolomics analyses

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.261ger4zdl47/v1>

Maria Jose Perez J.¹

¹Institut Imagine

Team Deleidi



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ger4zdl47/v1>

Protocol Citation: Maria Jose Perez J. 2025. Targeted LC–MS metabolomics analyses. protocols.io
<https://dx.doi.org/10.17504/protocols.io.261ger4zdl47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 26, 2024

Last Modified: February 03, 2025

Protocol Integer ID: 112830

Keywords: m, lc

Funders Acknowledgements:

ASAP

Abstract

Targeted LC–MS metabolomics analyses

- 1 Detach 5×10^5 iPSC-derived microglia using Accutase
- 2 Centrifuge the cells and immediately flash-freeze them in collection tubes on dry ice
- 3 Store samples at -80°C until metabolomics analysis
- 4 Prepare the extraction solution comprising 50% methanol, 30% acetonitrile (ACN), and 20% water. Use 1 mL of extraction solution per 1×10^7 cells
- 5 Add the extraction solution to the samples, vortex for 5 minutes at 4°C , and centrifuge at $16,000 \times g$ for 15 minutes at 4°C
- 6 Collect the supernatants and store them at -80°C until further analysis.
- 7 Collect the supernatants and store them at -80°C until further analysis.
- 8 Calibrate the mass spectrometer externally using a standard calibration mixture every seven days as per the manufacturer's instructions.
- 9 Inject 5 μL of the sample onto a ZIC-pHILIC column (150 mm $\times$ 2.1 mm; i.d. 5 μm) with a guard column (20 mm $\times$ 2.1 mm; i.d. 5 μm) (Millipore) for LC separation.
- 10 Use the following chromatographic gradient at a flow rate of 0.200 $\mu\text{L}/\text{minute}$:
 - 0–20 minutes: linear gradient from 80% to 20% of buffer B.
 - 20–20.5 minutes: linear gradient from 20% to 80% of buffer B.
 - 20.5–28 minutes: maintain 80% buffer B.
- 11
 - Prepare buffers:
 - Buffer A: 20 mM ammonium carbonate and 0.1% ammonium hydroxide (pH 9.2).
 - Buffer B: ACN.
- 12 Operate the mass spectrometer in full scan, polarity-switching mode with these parameters:
 - Spray voltage: 2.5 kV.
 - Heated capillary temperature: 320°C .

- Sheath gas flow: 20 units.
 - Auxiliary gas flow: 5 units.
 - Sweep gas flow: 0 units.
 - Mass range: 75–1,000 m/z at 35,000 resolution (at 200 m/z).
 - Automatic gain control (AGC) target: 10^6 .
 - Maximum injection time: 250 ms.
- 13 Use lock mass technique to ensure mass accuracy below 5 ppm.
- 14 Acquire data with Xcalibur software (Thermo Fisher Scientific).
- 15 Determine metabolite peak areas with TraceFinder software (Thermo Fisher Scientific) using exact mass and retention times.
- 16 Perform data analysis with MetaboAnalyst (version 6.0, <https://www.metaboanalyst.ca/>, RRID: SCR_016723):
- Statistical analysis module: conduct t-tests (unpaired, equal variance, FDR < 0.05), PCA, and heatmap generation.
 - Enrichment analysis module: perform MSEA using the Small Molecule Pathway Database (SMPDB; RRID: SCR_004844) with metabolite sets containing at least 2 entries.
 - Joint pathway analysis module: integrate metabolomics and bulk RNA-seq data using:
 - Hypergeometric Test for enrichment analysis
 - Degree Centrality for topology measure
 - Combine p-values (unweighted) for the integration method.